Service Repair Sheet SRS67866

Contact Name

Company/ Hospital Name

Department

Position

Direct Phone

General Phone

Opera Account

Email

Order Number

Date Received

Booked in By

Main Company

Type Return

Date Promised to Return: 01 Jan 1970 -3600

Ryan Filer

Musgrove Park Hospital

Med Electronics

Tech

0182334 2486

01823 333 444

00005110

Ryan.Filer@somersetft.nhs.uk

16/Jun/2021

Robert Connor

Viamed

Warranty

Notes 11/Jun/2021 Sophie Lines 11/Jun/2021 Sophie Lines 11/Jun/2021 Sophie Lines

Returning 1x 0110425

(SN: FK70799060) believed to be faulty.

16/Jun/2021 Robert Connor

MAX-250E s/n FK70799060.

Customer reports the sensor causes an O2 calibration error after 15 minutes of use.

Ready For quote

Repair Complete Signed

Screen 18:6:2021

SRN	Equipment		Serial Number	Warranty
SRN33936	Oxygen Sensor	0110425	FK70799060	Yes

0140013 ×1 0 60.

0110425 × 1 @ EO NEW S/N, SRS, SRN

UPS XI EO.



Taunton and Somerset

NHS Foundation Trust



Please investigate enclosed MAX-250 serial number FK70799061 reference SRS67866 as per attached email. Kind regards Gavin Reed

Medical Electronics Department Musgrove Park Hospital Taunton Somerset TA1 5DA

Tel: (01823) 342486

gavin.reed@somersetft.nhs.uk

With Compliments

I recently ordered 7x 0110425 sensors from you, PO 600025498. I fit these to Vapotherm Precision Plus devices which perform cell calibrations at increasing intervals.

After installing 1 of the cells, it caused an error related to the O2 calibration after 15 minutes. This error was repeatable on the same device, which was then tested with the cell that was removed and a different

I put the suspect cell in another device and generated the same error again after about 15 minutes. I think this cell, if not defective, is behaving differently to the others and is unusable.

I have now used 5 of the other 6 so it looks as though this was a one off issue. The serial number of this

I am not sure if cells are provided with a warranty, or whether you would like this one back for investigation. I do need another cell, so I either need a replacement or to order another one.

Please can you advise?

Many thanks