patient monitoring



Pippa Breathing Monitor

- Instantaneous or averaged breathing rate.
- Easy to use and accurate.
- Pre-set alarm limits.
- Dedicated masks with integrated sensor; universal sensor available for use with other masks
- Alternative to human observation or expensive multi-parameter patient monitors.



Utilising unique patented PEP sensors the Pippa breathing monitor reliably and accurately monitors patients' breathing rate.



Pippa Breathing Monitor

Applications for use:

- General patient monitoring
- Oxygen therapy
- Epidurals
- Continuous spinals
- PCA
- Nebulised medication

Areas for use include:

- Accident and Emergency
- Paramedic Ambulances
- Acute Medical Wards
- Respiratory Units
- Obstetric Units
- High Dependency Units
- Postoperative Wards
- Recovery Wards
- Operating Theatres

Display

Breathing rate: 2 digits

Resolution: 1 breath per minute

Accuracy

+/- 1 breath per minute

Modes of operation

Instantaneous breathing rate Averaged breathing rate

Display indications

AVERAGE BPM INSTANTANEOUS BPM

BATTERY LOW

Alarm activation

Audio alarm triggered when two consecutive breaths are monitored

- i) below 8 b.p.m., 20 seconds between breaths
- ii) above 30 b.p.m., 4.8 seconds between breaths or
- iii) when battery is low

Alarm deactivation

Alarm is deactivated when:

- rate returns within normal operational parameters, between 8-30 b.p.m.
- ii) Rate Select button is pressed during an alarm

Dimensions

14.0 cm (l) x 6.4 cm (w) x 3.2 cm (h)

Weight

130g

Batteries

2 x AA, 1.5 V



Patent Pending

The Pippa breating monitor is an hand held device, but can be easily attached to a pocket or fastened onto a belt by using the integral securing clip.

Pippa monitoring kit includes:

Qty	Part number	Description
1	1610100	Pippa breathing monitor
1	1610200	Patient lead, monitor-sensor
1	1620010	Carrying case
1pk	1655000	Battery pack, AA, 1.5V. Pack of 2
1pk	1650010	Sensor, (universal). Pack of 5 x 1650009
1pk	1650050	Mask (Flexicare) with integrated sensor. Pack of 5 x 1650049
1	1690000	Instruction manual



Mask with integrated PEP sensor.



Universal and integrated PEP sensors.

Specifications subject to change.

Please contact Viamed for further specifications and information.

Distributed by:



Viamed Limited · 15 Station Road · Cross Hills Keighley · West Yorkshire · BD20 7DT · United Kingdom Tel: +44 (0)1535 634542/636757 Fax: +44 (0)1535 635582 Email: info@viamed.co.uk Website: www.viamed.co.uk



BS EN ISO 9001 BS EN 46001 BS EN ISO 13485