

1. Label Layout: VIA Connector Sampling Lines

Dimension Pouch Label: 170 mm x 110 mm

Dimension Box Label: 150 mm x 127 mm

Ref.: 7090121001V, Type: 4420800

Single Port Nasal Prong, Adult, Short Term

VIAMED Viamed Limited, 15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, United Kingdom
www.viamed.co.uk, info@viamed.co.uk, +44(0)1535 634542

VersaStream®
Single Port Nasal Prong Sampling Line, Adult, 2m/7ft
Short Term Use ≤ 24h

REF 7090121001V LOT 21121.04.14 # 4420800

10 3% 95% 650 hPa 1060 hPa -13 °C -25 °C 140 °F 60 °C

MR

MONITOR

H₂O

CO₂ ≤ 75 ml/min

UDI (01) 1 42516796 0111 8
(11) 210414
(17) 240414
(10) 21121.04.14

MD

bluepoint MEDICAL GmbH & Co KG
An der Trave 15
23923 Selmsdorf
Germany
www.bluepoint-medical.com

≥ 10kg
(≥ 22 lb)

2021-04-14

2024-04-14

CE

A21Rev.3

Pouch Label

25 Single Port Nasal Prong Sampling Line, disposable, Adult, 2m/7ft

REF 7090121001V LOT 21121.04.14

cs Odběrová hadička s adaptérem pro dýchací cesty, jednorázová, pro dospělé/pediatrická

da Udfagningslange med luftvejsadapter, engangsbrug, voksne/børn

de Einweg-Probierabnahmeleitung mit Atemwegadapter, Erwachsene/Pädiatrie

et Hingamisestee adapter provivitoruga, ühekordsest kasutata, täiskasvanu/lapse

el Τηροσφαιστικός απορροφητής της γαστρίδας, βενζομετόλης/οξυγόνου, ενήλικων/παιδίων

es Línea de muestreo con adaptador para vías aéreas, desechable, adulto/pediatrico

fr Tubulure d'échantillonnage avec adaptateur pour voies aériennes, usage unique, adulte/enfant

ko 일회용 채취관 기도 연결자, 일회용, 성인/소아용

hr Linija za uzorkovanje s dišnim adapterom, jednokratna, odrasli/pedijatrija

it Linea di campionamento con adattatore per vie aeree, usa e getta, pazienti adulti/pediatrici

lv Elpošu adaptera paraugu ņemšanas līnija, vienureizējai lietojama, pieaugušajiem/bērniem

lt Magišų ėmimo linija su kvėpavimo takų adapteriu, vienkartinė, suaugusiems/vaikams

hu Légűlő adapter mintavételhez, egyszeri használatos, felnőtt/gyerek

nl Bemonsteringslijn met luchtwegadapter, voor eenmalig gebruik, volwassen/pediatrisch

ja 気道アダプタサンプリングライン、使い捨て、成人/小児用

no Provetakingsledning med luftveisadapter, til engangsbruk, voksne/pediatrisk

pl Adapter drog oddechowych dla linii próbkującej, jednorazowy, dorosły/dziecko

pt Linha de amostragem com adaptador de vias aéreas, descartável, adulto/pediatrico

pt-br Linha de amostragem com adaptador de vias aéreas, descartável, adulto/pediatrico

ro Linie de eșantionare cu adaptor pentru căile respiratorii, de unică folosință, adult/pediatric

ru Пробоотборная линия с адаптером дыхательных путей, одноразовая, для взрослых/детей

sk Hadička na odber vzoriek s adaptérom pre dýchacie cesty, jednorázová, pre dospelých/pediatrická

sv Hadika på uzorkování sa adapterom za dýšací puť, za jednokratnú upotrebu, za odrasle/detú

fi Ilmatte sovitimen näytelinja, kertakäyttöinen, aikuiset/lapset

sv Provtagningsledning för luftvägsadapter, för engangsbruk, vuxen/barn

tr Havayolu Adaptörü Numune Alma Hatı H, tek kullanımlık, Yetişkin/Çocuk

zh 气体采样管适配器，一次性，成人/儿童

4420800

≥ 10 kg / 22 lbs

Short Term Use ≤ 24h

2 10 3% 95% 650 hPa 1060 hPa -13 °C -25 °C 140 °F 60 °C

VIAMED

Viamed Limited, 15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, United Kingdom
www.viamed.co.uk, info@viamed.co.uk, +44(0)1535 634542

2021-04-14

2024-04-14

bluepoint medical GmbH & Co. KG
An der Trave 15
23923 Selmsdorf, Germany
www.bluepoint-medical.com

CE

UDI (01) 1 42516796 0111 5
(11) 210414
(17) 240414
(10) 21121.04.14


A21Rev.3

Box Label



Viamed Limited
15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, U.K.
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Email: info@viamed.co.uk
Website: www.viamed.co.uk

Ref.: 7090121002V, Type: 4420801
Nasal Sampling Line, Adult, Short Term




Pouch Label




Box Label

Ref.: 7090121003V, Type: 4420802
Nasal Sampling Line, Paediatric, Short Term



Pouch Label




Box Label



Viamed Limited
15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, U.K.
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Email: info@viamed.co.uk
Website: www.viamed.co.uk

Ref.: 7090121004V, Type: 4420803
Nasal Sampling Line, Infant, Short Term




Pouch Label



Box Label

Ref.: 7090121005V, Type: 4420804
Nasal Sampling Line with O2 Supply, Adult, Short Term



Pouch Label




Box Label



Viamed Limited
15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, U.K.
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Email: info@viamed.co.uk
Website: www.viamed.co.uk

Ref.: 7090121006V, Type: 4420805
Nasal Sampling Line with O₂ Supply, Pediatric, Short Term




Pouch Label




Box Label

Ref.: 7090121007V, Type: 4420806
Nasal Sampling Line with O₂ Supply, Infant, Short Term



Pouch Label



Box Label



VIAMED

Viamed Limited
15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, U.K.
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Email: info@viamed.co.uk
Website: www.viamed.co.uk

Ref.: 7090121008V, Type: 4420807
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121002V, Type: 4420801
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121009V, Type: 4420808
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121003V, Type: 4420802
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121010V, Type: 4420809
Nasal Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121004V, Type: 4420803
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121011V, Type: 4420810
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121005V, Type: 4420804
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121012V, Type: 4420811
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121006V, Type: 4420805
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h


Ref.: 7090121013V Type: 4420812
Nasal Sampling Line with O2 Supply, Infant, Long Term

Is identical with Ref.: 7090121007V, Type: 4420806
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h



Viamed Limited
15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, U.K.
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Email: info@viamed.co.uk
Website: www.viamed.co.uk

Ref.: 7090121014V, Type: 4420813
Oral/Nasal Sampling Line, Adult, Short Term




Pouch Label



Box Label

Ref.: 7090121015V, Type: 4420814
Oral/Nasal Sampling Line, Paediatric, Short Term



Pouch Label




Box Label



Viamed Limited
15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, U.K.
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Email: info@viamed.co.uk
Website: www.viamed.co.uk

Ref.: 7090121016V, Type: 4420815
Oral/Nasal Sampling Line with O2 Supply, Adult, Short Term




Pouch Label




Box Label

Ref.: 7090121017V, Type: 4420816
Oral/Nasal Sampling Line with O2 Supply, Paediatric, Short Term



Pouch Label



Box Label



VIAMED

Viamed Limited
15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, U.K.
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Email: info@viamed.co.uk
Website: www.viamed.co.uk

Ref.: 7090121018V, Type: 4420817
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121014V, Type: 4420813
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121019V, Type: 4420818
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121015, Type: 4420814
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121020, Type: 4420819
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121016V, Type: 4420815
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h


Ref.: 7090121021V, Type: 4420820
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121017V, Type: 4420816
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h



Viamed Limited
15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, U.K.
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Email: info@viamed.co.uk
Website: www.viamed.co.uk

Ref.: 7090121022V, Type: 4420821
Airway Adapter Sampling Line, Adult/Pediatric, Short Term




Pouch Label



Box Label

Ref.: 7090121023V, Type: 4420822
Airway Adapter Sampling Line, Infant, Short Term



Pouch Label



Box Label



VIAMED

Viamed Limited
15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, U.K.
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Email: info@viamed.co.uk
Website: www.viamed.co.uk

Ref.: 7090121024, Type: 4420823

Airway Adapter Sampling Line, Adult/Pediatric, Long Term

Is identical with Ref.: 7090121022V, Type: 4420821

Only Change: Short Term Use $\leq 24h$ $\rightarrow$ Long Term Use $\leq 72h$

Ref.: 7090121025V, Type: 4420824

Airway Adapter Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121023V, Type: 4420822

Only Change: Short Term Use $\leq 24h$ $\rightarrow$ Long Term Use $\leq 72h$



Viamed Limited
15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, U.K.
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Email: info@viamed.co.uk
Website: www.viamed.co.uk

2. Label Layout: LLM Connector Sampling Lines

Ref.: 7090121033V, Type: 4420857
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121002V, Type: 4420801
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121034V, Type: 4420858
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121003V, Type: 4420802
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121035V, Type: 4420859
Nasal Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121004V, Type: 4420803
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121036V, Type: 4420860
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121005V, Type: 4420804
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121037V, Type: 4420861
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121006V, Type: 4420805
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121038V, Type: 4420862
Nasal Sampling Line with O2 Supply, Infant, Long Term

Is identical with Ref.: 7090121007V, Type: 4420806
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h



VIAMED

Viamed Limited
15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, U.K.
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Email: info@viamed.co.uk
Website: www.viamed.co.uk

Ref.: 7090121043V, Type: 4420867
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121014V, Type: 4420813
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121044V, Type: 4420868
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121015V, Type: 4420814
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121045V, Type: 4420869
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121016V, Type: 4420815
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121046V, Type: 4420870
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121017V, Type: 4420816
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121049V, Type: 4420873
Airway Adapter Sampling Line, Adult/Pediatric, Long Term

Is identical with Ref.: 7090121022V, Type: 4420821
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121050V, Type: 4420874
Airway Adapter Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121023V, Type: 4420822
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h



VIAMED

Viamed Limited
15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, U.K.
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Email: info@viamed.co.uk
Website: www.viamed.co.uk

3. Label Layout: ORI Connector Sampling Lines

Ref.: 7090121102V, Type: 4420900
Single Port Nasal Prong, Adult, Short Term

Is identical with Ref.: 7090121001V, Type: 4420800

Ref.: 7090121103V, Type: 4420901
Nasal Sampling Line, Adult, Short Term

Is identical with Ref.: 7090121002V, Type: 4420801

Ref.: 7090121104V, Type: 4420902
Nasal Sampling Line, Pediatric, Short Term

Is identical with Ref.: 7090121003V, Type: 4420802

Ref.: 7090121105V, Type: 4420903
Nasal Sampling Line, Infant, Short Term

Is identical with Ref.: 7090121004V, Type: 4420803

Ref.: 7090121106V, Type: 4420904
Nasal Sampling Line with O2 Supply, Adult, Short Term

Is identical with Ref.: 7090121005V, Type: 4420804

Ref.: 7090121107V, Type: 4420905
Nasal Sampling Line with O2 Supply, Pediatric, Short Term

Is identical with Ref.: 7090121006V, Type: 4420805

Ref.: 7090121108V, Type: 4420906
Nasal Sampling Line with O2 Supply, Infant, Short Term

Is identical with Ref.: 7090121007V, Type: 4420806



Viamed Limited
15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, U.K.
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Email: info@viamed.co.uk
Website: www.viamed.co.uk

Ref.: 7090121109V, Type: 4420907
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121002V, Type: 4420807
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121110V, Type: 4420908
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121003V, Type: 4420808
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121111V, Type: 4420909
Nasal Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121004V, Type: 4420809
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121112V, Type: 4420910
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121005V, Type: 4420810
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121113V, Type: 4420911
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121006V, Type: 4420811
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121114V, Type: 4420912
Nasal Sampling Line with O2 Supply, Infant, Long Term

Is identical with Ref.: 7090121007V, Type: 4420812
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121115V, Type: 4420913
Oral/Nasal Sampling Line, Adult, Short Term



Viamed Limited
15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, U.K.
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Email: info@viamed.co.uk
Website: www.viamed.co.uk

Ref.: 7090121116V, Type: 4420914
Oral/Nasal Sampling Line, Pediatric, Short Term

Is identical with Ref.: 7090121015V, Type: 4420814

Ref.: 7090121117V, Type: 4420915
Oral/Nasal Sampling Line with O2 Supply, Adult, Short Term

Is identical with Ref.: 7090121016V, Type: 4420815

Ref.: 7090121118V, Type: 4420916
Oral/Nasal Sampling Line with O2 Supply, Pediatric, Short Term

Is identical with Ref.: 7090121017V, Type: 4420816

Ref.: 7090121119V, Type: 4420917
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121014V, Type: 4420817
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121120V, Type: 4420918
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121015V, Type: 4420818
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 70901210121V, Type: 4420919
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121016V, Type: 4420819
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121122V, Type: 4420920
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121017V, Type: 4420820
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h



Viamed Limited
15 Station Road, Cross Hills, Keighley
West Yorkshire, BD20 7DT, U.K.
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582
Email: info@viamed.co.uk
Website: www.viamed.co.uk

Ref.: 7090121123V, Type: 4420921
Airway Adapter Sampling Line, Adult/Pediatric, Short Term

Is identical with Ref.: 7090121022V, Type: 4420821

Ref.: 7090121124V, Type: 4420922
Airway Adapter Sampling Line, Infant, Short Term

Is identical with Ref.: 7090121023V, Type: 4420822

Ref.: 7090121125V, Type: 4420923
Airway Adapter Sampling Line, Adult/Pediatric, Long Term

Is identical with Ref.: 7090121022V, Type: 4420823
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h

Ref.: 7090121126V, Type: 4420924
Airway Adapter Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121023V, Type: 4420824
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h


Ref.: 7090121140V, Type: 4420925
Airway Adapter Sampling Line, Adult/Pediatric, Long Term

Is identical with Ref.: 7090121022V, Type: 4420825
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h High Humidity


Ref.: 7090121141V, Type: 4420926
Airway Adapter Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121023V, Type: 4420826
Only Change: Short Term Use $\leq$ 24h $\rightarrow$ Long Term Use $\leq$ 72h High Humidity

Tag Layout:



Tag for short term sampling lines



Tag for long term sampling lines