

Service Repair Sheet SRS67929

Contact Name

Company/ Hospital Name

Department

Position

Direct Phone

General Phone

Opera Account

Email

Order Number

Date Received

Booked in By

Main Company

Type Return

Date Promised to Return: 01 Jan 1970 -3600

David Lloyd

Southport District General Hospital

EBME Dept

Technical Practitioner of Specialist Services

01704 704188

00004760

davidlloyd1@nhs.net

19/Aug/2021

Robert Connor

Viamed

Quote

Notes 16/Aug/2021 Kate Griffiths

16/Aug/2021 Kate Griffiths

Customer to send in DB3 for us ``to test it side by side with a reference device to determine the extent of the drop-off and compare it with what we would expect to see on a new unit`` as per SH. S/N H0003122.

19/Aug/2021 Robert Connor

Received 1 x Microstim DB3 s/n H0003122

Ready For quote

Repair Complete Signed

CGreen. 3-9-2021

SRN	Equipment	Stock Ref	Serial Number	Warranty
SRN34089	Microstim DB3	2510000	H0003122	2

2540000 x 1 FOC
Returned unrepaired as requested
S/N, SRS, SRN

UPS x 1 FOC.

**Southport and
Ormskirk Hospital**
NHS Trust

Repairs Team SRS67929
Viamed Limited
15 Station Road
Cross Hills, Keighley
West Yorkshire
BD20 7DT

EBME Department
Southport & Formby District General Hospital
Town Lane
Southport
Merseyside
PR8 6PN

17TH August 2021

Returns reference number SRS67929

Following discussions with Steve Hardaker, I have returned the DB3 stimulator to verify its stimulation output.

When I have tested the TOF output of this device using an oscilloscope, the signal amplitude decreases and increases from one TOF to another. For example, TOF stimulation amplitude is measured at 800mv, then the next TOF stimulation amplitude can be 400mv, then the next 600mv. A pause of at least 10 seconds is left between each TOF activation.

If you require any further information, please contact myself.

Regards

David Lloyd

Tel: 01704 704188
Email: davidlloyd1@nhs.net

Service Repair Sheet 67929

Contact Name

David Lloyd

Company/ Hospital Name

Southport District General Hospital

Department

EBME Dept

Position

Technical Practitioner of Specialist Services

Direct Phone

01704 704188

Opera Account

00004760

Email

davidlloyd1@nhs.net

Date Received

19/Aug/2021

Booked in By

Robert Connor

Repair

Ref

S/N

Equipment Type

SRN34089

2510000

H0003122

Microstim DB3

No Fault Found

Time :0 Hour(s)

The Microstim DB3 has been tested. The fault could not be reproduced.
