



Viamed Limited - 15 Station Road - Cross Hills  
Keighley -West Yorkshire BD20 7DT - United Kingdom  
Tel: +44 (0)1535 634542 Fax: +44 (0)1535 635582  
Email: [info@viamed.co.uk](mailto:info@viamed.co.uk) Website: [www.viamed.co.uk](http://www.viamed.co.uk)

## **Calibration Accuracy – VM2500 Series Capnographs:**

**4410500 Capnograph VM2500-M, EU version**  
**4410501 Capnograph VM2500-M, Asia version**  
**4410520 Capnograph VM2500-S, EU version**  
**4410521 Capnograph VM2500-S, Asia version**

The capnography accuracy specification of the VM2500-M and VM2500-S are as per the current published specifications.

The IRMA sensor used with mainstream VM2500-M and the ISA sensor integrated into the VM2500-S are calibrated at:

**5% vol % CO<sub>2</sub>, balance N<sub>2</sub> & 15% vol % CO<sub>2</sub>, balance N<sub>2</sub>**

The units are further tested at **7.5% vol % CO<sub>2</sub>, balance N<sub>2</sub>**

The accuracy specification for the gases is **+/- 0.03%**