



Zoey Teal <viamed.zoey.teal@gmail.com>

Fwd: Update to PAQ2 re Infection Control

1 message

Main Account <office@viamed.co.uk>
To: Zoey Teal <zoey.teal@vmsecure.me.uk>

Fri, Sep 17, 2021 at 9:39 AM

----- Forwarded message -----

From: **Steve Hardaker** <office@viamed.co.uk>
Date: Fri, 17 Sept 2021 at 09:33
Subject: Fwd: Update to PAQ2 re Infection Control
To: Main Account <office@viamed.co.uk>

Good Friday morning, office! :-)

Please can someone raise a quotation for me for Kings Mill Hospital, attn Peter Lee, Medical Equipment Management Department (MEMD). Full contact details in his email below.

- 4x 2510091 - TOF 3D kit
- 4x 2520110 - Hand adapter
- 4x 2520113 - Pole mount
- 2x 2520112 - Eye adapters, pack/50

If it can be checked and sent to me, I will forward it to Peter as I want to clarify that these are just I.V. Pole adapters and not I.V. Poles as suggested by his comment.

Thanks,

Steve

----- Forwarded message -----

From: **LEE, Peter (SHERWOOD FOREST HOSPITALS NHS FOUNDATION TRUST)** <peter.lee3@nhs.net>
Date: Fri, 17 Sept 2021 at 08:25
Subject: Re: Update to PAQ2 re Infection Control
To: steve.hardaker@viamed.co.uk <steve.hardaker@viamed.co.uk>
Cc: BURGOYNE, Jamie-Rae (SHERWOOD FOREST HOSPITALS NHS FOUNDATION TRUST) <jamie-rae.burgoyne@nhs.net>

Hi Steve

Pls can we have a quote for x4 with IV Poles, all with hand sensor and 2 with eye with "starter" pack?

Kind regards

Peter D. Lee C.Eng MIET MIPEM

Consultant Head, Clinical Engineering Services & Trust Medical Device Safety Officer
Incorporating Medical Equipment Management Department (Memd)
Diagnostics and Outpatients Division
Kings Mill Hospital
Sherwood Forest NHS Foundation Trust
Mansfield Road
Sutton-in-Ashfield
NG17 4JL

01623 622515 x 3882

Please visit our websites on the following link <http://sfhnet.nnotts.nhs.uk/memd/default.aspx> or to report routine repairs <http://sfh-ram/maintenance5000/mainmenu.aspx?Guest=True>

Our missionMaximizing the availability of safe medical devices for your patients...Choosing it right- and using it right - will keep it right!

Service Registered with BSi continuously to Iso9001 Quality Management Standard since 1990

On 16 Sep 2021, at 16:04, Steve Hardaker <office@viamed.co.uk> wrote:

Hi Peter,

Following the receipt of further information from MIPM, I have been able to update the PAQ2 form with additional information concerning infection control, section 6.

Please see attached revision, dated today to supersede previous version. The updates are as follows:

The manufacturer, MIPM, has provisionally approved Gama Healthcare Clinell Universal Wipes based on assessment of the composition of the wipes and the device but has not yet issued a written declaration of compatibility. Once we receive this, we will forward it on.

Q6b: Also, please provide specific detail of the plastics used in all rigid external surfaces:

A: 1)The upper and lower casing including the battery cover is produced from amorphous thermoplastic alloy (polycarbonate + ASA)

2) The upper exposed surface of the membrane keypad is produced from textured hardcoat polyester film.

3) The transparent display protection window is produced from 2mm PMMA (Polymethyl methacrylate - acrylic glass).

The device is just going through reprocessing, which includes cleaning and disinfection. We should be able to despatch tomorrow or Monday to be with you by Tuesday at the latest. I will address it directly to you, care of MEMD.

Regards,

Steve Hardaker
Technical Support Manager
Viamed Ltd.

Please note: Viamed is enacting a coronavirus contingency plan to allow sales and admin staff to work from home, and I am now working remotely. Telephone calls to the main office will be answered remotely, but please continue to use email where possible.

Website: www.viamed.co.uk
email: steve.hardaker@viamed.co.uk
Tel: +44 (0)1535 634542
Fax: +44 (0)1535 635582

Sign up to our mailing list, go to www.viamed.co.uk/sign-up
Twitter: twitter.com/ViamedLtd
Facebook: facebook.com/ViamedLtdUK
Company registered in England, No. 1291765.

E-mail Disclaimer

The information in this email is confidential and may be legally privileged.

It is intended solely for the addressee. Access to this email by anyone else is unauthorized.

If you are not the intended recipient, any disclosure, copying, distribution or any action taken or omitted to be taken in reliance on it, is prohibited and may be unlawful.

If you feel you have received this email in error or would like to be removed from our mailing list, please reply with delete in the subject line or call +44 (0)1535 634542

GDPR / Privacy Policy available at www.viamed.co.uk/policy

This message originated from outside of NHSmail. Please do not click links or open attachments unless you recognise the sender and know the content is safe.

<paq2 PPQ ADDITIONAL INFORMATION v4 06.06.19 (SFHNHSFT) 20210916.doc>

This message may contain confidential information. If you are not the intended recipient please:

- i) inform the sender that you have received the message in error before deleting it; and
 - ii) do not disclose, copy or distribute information in this e-mail or take any action in relation to its content (to do so is strictly prohibited and may be unlawful).
- Thank you for your co-operation.

NHSmail is the secure email, collaboration and directory service available for all NHS staff in England. NHSmail is approved for exchanging patient data and other sensitive information with NHSmail and other accredited email services.

For more information and to find out how you can switch visit [Joining NHSmail – NHSmail Support](#)