

**Residual Risk Analysis Report.****Rotating Acrylic Bed**

The centre of the rotating bed is a bolt which is a source of dirt collection.

When the bed is cleaned this bolt should also be cleaned and checked that the bed rotates freely.

The mattress is round but apart from its shape should be cleaned and handled as in the general body of the manual on bed care and cleaning

The bed should never be left between these positions either in use or in storage.

Care should be taken to prevent the operators fingers becoming trapped in the mechanism

Care should be taken to ensure that the patients limbs do not protrude out of the bed perimeter when the bed is being rotated

The patient dropping off. This should never happen as the patient should never be left unattended

**Electrical**

If the heater case is removed with the unit switched on

Wiring failures if wires become disconnected

Element is touched.

Component failure

**Heat**

Patient overheat

User over heat

**Physical Injury**

.head of user coming into contact with the heater

See individual risk reports

**WARNINGS**

See label section and user manual

J.S.lamb

11/02/11