

VIAMED COMPANY OPERATING PROCEDURE

Q.A. Packing Nova Oxygen TentsTitle


VM3/COP/ 59.01

Date: 26/03/2007

Revision Date:10-Jun-11

Issue: 1


1. Unpack the Oxygen Hood from its original packing, and check that:
 - a. The product is the correct item
 - b. The product is the correct size:
 - c. Remove Nova leaflet
 - d. Remove Nova label from tent and ensure no residue is left on the hood



Add Viamed Label (0195001 Viamed clear) where the nova label was removed from.

Add Serial Number label to the top right hand corner of the front

Check the product visually for damage.
Add the insert leaflet 239001. One per hood.
Ensure the product and outer carton are clean.
Add labels to original box



Pt No. 2390020
(Cover Nova
Logo)
Barcode label
Remove label
Address label

