

	Compatibility of Components	Rev. 1
Internal Document	Multi Parameter Monitors VM-2500-M / S	Page 1 of 2

Compatibility of Components - VM-2500-M / S

Hardware main / power PCB (Modification)	Firmware	PC-Software	Product CD	Remark
V1.9 / V1.1 (M1,M2)	V4.5 – V5.0	V1.3 & V1.4	V1.2 - V1.4	Move power on/off to confirm key. Move alarm reset to alarm-mute key.
V1.9 / V1.1 (M1,M2,M3)	V5.1	V1.3 & V1.4	V1.2 - V1.4	Changed the RTC routine for the new RTC IC (PCF8563)
V1.9 / V1.1 (M1,M2,M3)	V5.2 – V5.8	V1.6	V1.5	Updated the alarms log list. Masked the alarm start and end time
V1.9 / V1.1 (M1,M2,M3)	V 5.9	V1.7	V1.6	Added second to the alarms log. If Firmware V5.9 is used with PC software V1.6, the alarms are not displayed in the “ Patient Info” Tab

