

“TEE” Adapter**Product Life.****Warranty.**

No Warranty once the product has been used

Life of Product.

No claim is made to life expectancy.

The information that could be used to be supportive is too subjective.

It is dependent upon the method of sterilization, application and general handling

The Blue Polyethylene Tee is offered as a single use product but historically it has been cleaned and used many times before discarding.

As it is normally used outside the area of patient contact usually on the output of ventilators etc. simple cleaning is all that is required.

Polycarbonate and Polyethylene are impervious to most hospital cleaners.

It can be scratched and may break or crack if dropped or misused

Supply of spares after production has ceased.

Not Applicable

Service after production has ceased.

Not Applicable