

Packaging Trial & Validation

The following should be checked:

- The unit being sent is to the specification as detailed on the official order.
- The instructions for use are included.

Packaging:

- Use one purpose made Cardboard box
- Pack the leads and the battery into the base of the box.
- Pack the instrument into the centre compartment of the box.
- Pack the “Instructions for Use” into the top of the box.
- Place in a cardboard box with airbags surrounding – applied with use of the Specific airbag machine – (Cell-O-Ez – Air Cushion Machine).
- Alternatively, for small orders they should be despatched in Jiffy bags

Postage:

- For single and multiple orders – Despatch is via an applicable courier with adequate insurance.
- Jiffy bagged product can be sent either by courier or by Post
- Alternative Delivery is via a Hospital visit by a Sales Representative.