

Master Blank**Key Rings for Tek-Ox
VM3/COP/40.46**

Date: 30-Sep-03

Revision Date: 24-Jan-
12

Issue 1

Only to be attempted by trained personel

Parts Required

Quantity	Description	Part No.	
----------	-------------	----------	--

Procedure

1. , There is a slight difference in the key ring pillar threads on the Microstim and the Tek-ox
2. The cure is to run a tap down the thread and a die through the insert.
3. Both Tools can be found in the Vandagraph Stock room 10. Row 4 Block 3 Shelf 2