

Email-shot – Export distributors – May 2107

Subject: Viamed Oxygen Hoods

Dear ***,

Viamed Oxygen Hoods - Designed for the short-term stable delivery of oxygen, for use with neonatal and paediatric patients without the complication of applying and using nasal cannulae.

Disposable Oxygen Hoods

See through environment provides full patient visibility
Drape neck opening provides easy access to the patient
Single patient use means reduced risk of cross-infection
Flat packed for easy storage and economical transportation

The standard oxygen hood has two 22mm inlet ports, a full vinyl bottom and a drape neck opening, and comes in a variety of sizes.

Other style variants of the disposable oxygen hoods are also available, including; No-Bottom variations, NICU, Low Birthweight

Reusable Oxygen Hoods

These Perspex® reusable oxygen hoods feature a sliding front aperture, and can be supplied with or without side access ports.

Standard Hood


This hood can be used in conjunction with our phototherapy lightshields for use during phototherapy.

Hood With Side Ports

Side ports enable the user to have instant access to the patient.

Please Note: *Oxygen Hoods should only to be used in conjunction with oxygen monitoring.*

Please find attached information relating to the Teledyne AX300 Oxygen Analyser and the Teledyne MX300 Oxygen Monitor, available from Viamed.


R. Swain
9/5/17