



Steve Nixon <steve.nixon.viamed@googlemail.com>

Fwd: Polymer resin oxygen sensor

1 message

Steve Nixon <steve.nixon@viamed.co.uk>

23 February 2018 at 14:59

To: Jens Schwarz <j.schwarz@sensatronic.com>

Cc: Bernd Lindner <b.lindner@bluepoint-medical.com>, John Lamb <john.lamb@viamed.co.uk>

Bcc: Derek Lamb <derek@lamb.uk.net>

Hi Jens

Further to our meeting at EC Sense, would it be possible to draw up the framework for a provisional datasheet for the sensor; in terms of dimensions, power requirements, interface, range, resolution, accuracy, performance, output..?

From discussions we have three initial target applications:

- 1) Use as an on-board oxygen alarm, i.e. to alarm if there is a leakage of oxygen. So, in principle the sensor is functional now as it works up to circa 40%. We can also use the demo box to initially prove the principle to some OEMs.
- 2) Medical applications up to 100%. Requires more testing and prototypes.
- 3) Rebreather applications, we need to see if the sensor can work in the region of 1.1 to 1.3 bar ppO2 and for testing purposes it needs to operate at say 2.0 bar ppO2 to a maximum of 2.3 bar ppO2.

For 1) above I initially see the sensor mounted on an interface board. The problem is that we don't know how long the sensor will last, it could last too long and we significantly reduce our revenue stream.

Perhaps we need to limit the life due to say the wear of the electrode or the capillaries potentially blocking up and being compromised? Could we use a Maxim Dallas chip with an accumulative actual or 'in use' work time clock? Or should we power the board and sensor with an integral battery, which in turn limits the life of the sensor and or the interface PCB. The whole PCB could be conformal coated (for protection).

Do you have any thoughts on the above?

How do we finance the projects?

1 above) Would in principle need just a little investment, plus input from Sensatronic.

2 above) Would need Peter/Axel's development time, plus input from Sensatronic.

3 above) Is open ended at the moment, until we prove in principle that it works. Would it be possible for Vamed/VST to have exclusivity for this market? We already have the long established contacts for the technical and military markets.

Form the meeting at EC sense it sounded as though they need an extra person, not to carry out R&D, but to run the production processes. This would then free up Peter an Axel's time to carry out the R&D evaluations.

Regards

Steve

----- Forwarded message -----

From: **Steve Nixon** <steve.nixon@viamed.co.uk>

Date: 19 January 2018 at 17:01

Subject: Re: Polymer resin oxygen sensor

To: J Schwarz <j.schwarz@sensatronic.com>

Thanks Jens

So, we only need say 5V supply to power the unit? In principle it could be powered by standard batteries?

Do you know what the T90 response time of the sensor is?

Regards

Steve

On 19 January 2018 at 14:46, J Schwarz <j.schwarz@sensatronic.com> wrote:

Hi Steve,

I am sending attached a short instruction for the demo unit.

Regarding the sensor specs, here are some facts:

Power requirement: Depends on the driving electronics plus sensor. The sensor itself requires a few micro amps only;

Signal output: The sensor is giving few μA @21%O₂. The output of the electronics is V or mV. The sensor requires a potentiostatic driving circuit with I/U transducer.

Drive electronics: Yes, I think it makes sense that we supply the drive electronics. They can also develop their own, but it may take much time. We have some experience in the meantime. However, first we need to understand the application and then we have to adjust the electronics for this.

Communication interface: We can go with RS232 or USB, whatever is required. The USB interface at the demo unit is for charging only. There is actually no usable protocol implemented.

Function of LEDs: see the attached IfU

Regarding sample sensors, I will try to find some and send them to you on Monday.

Have a good weekend.

Mit freundlichen Grüßen,

Kind regards,

Jens Schwarz

Sensatronic GmbH

Am Ring 9

D-23970 Wismar, Germany

Tel. +49-3841-2252 440

Fax +49-3841-2252 442

j.schwarz@sensatronic.com

Diese E-Mail enthält vertrauliche und/oder rechtlich geschützte Informationen. Wenn Sie nicht der richtige Adressat sind oder diese E-Mail irrtümlich erhalten haben, informieren Sie bitte sofort den Absender und vernichten Sie diese E-Mail. Das unerlaubte Kopieren sowie die unbefugte Weitergabe dieser E-Mail sind nicht gestattet.

This e-mail may contain confidential and/or privileged information. If you are not the intended recipient (or have received this e-mail in error) please notify the sender immediately and destroy this e-mail. Any unauthorized copying, disclosure or distribution of the material in this e-mail is strictly forbidden.

Von: steve.nixon.viamed@googlemail.com [mailto:steve.nixon.viamed@googlemail.com] **Im Auftrag von** Steve Nixon

Gesendet: Freitag, 19. Januar 2018 13:38

An: J Schwarz <j.schwarz@sensatronic.com>

Betreff: Re: Polymer resin oxygen sensor

Sorry Jens, the unit has just arrived.

What are the power requirements of the sensor?

What is the signal output?

Would we supply the drive electronics?

What would be the communications interface?

At some stage can you please let me have a couple of examples of the actual sensors. They don't need to be working.

I'll charge up the unit and see what it does.

In the meantime what are the functions of the LEDs?

Blue

Amber

Green

Red

Steve

On 19 January 2018 at 12:01, Steve Nixon <steve.nixon@viamed.co.uk> wrote:

Hi Jens, did you ship the demo unit and samples sensors?

Steve

On 11 January 2018 at 15:38, Steve Nixon <steve.nixon@viamed.co.uk> wrote:

That's great news, thank you Jens.

Steve

On 11 January 2018 at 15:13, J Schwarz <j.schwarz@sensatronic.com> wrote:

Hi Steve,

Sorry for being late with the demo unit. We had to solve a problem with an electronic component that we used for the old (CO2) demo unit (charging circuit), but which no longer available. We had to substitute this by another component which twice killed the charging circuit. This problem is solved now. We have the unit running for tests since three days and will send it out tomorrow morning.

I will send a short instruction separate by email.

Mit freundlichen Grüßen,

Kind regards,

Jens Schwarz

Sensatronic GmbH

Am Ring 9

D-23970 Wismar, Germany

Tel. +49-3841-2252 440

Fax +49-3841-2252 442

j.schwarz@sensatronic.com

Diese E-Mail enthaelt vertrauliche und/oder rechtlich geschuetzte Informationen. Wenn Sie nicht der richtige Adressat sind oder diese E-Mail irrtuemlich erhalten haben, informieren Sie bitte sofort den Absender und vernichten Sie diese E-Mail. Das unerlaubte Kopieren sowie die unbefugte Weitergabe dieser E-Mail sind nicht gestattet.

This e-mail may contain confidential and/or privileged information. If you are not the intended recipient (or have received this e-mail in error) please notify the sender immediately and destroy this e-mail. Any unauthorized copying, disclosure or distribution of the material in this e-mail is strictly forbidden.

Von: steve.nixon.viamed@googlemail.com [mailto:steve.nixon.viamed@googlemail.com] **Im Auftrag von** Steve Nixon

Gesendet: Donnerstag, 11. Januar 2018 11:26

An: J Schwarz <j.schwarz@sensatronic.com>

Cc: Bernd Lindner <b.lindner@bluepoint-medical.com>

Betreff: Polymer resin oxygen sensor

Hi Jens

Sorry to hassle you, but what is the schedule for the demo unit of the polymer resin oxygen sensors? The reason why I'm enquiring is that I am planning more customer visits, so I can delay a little if the demo unit is due soon.

--

Steve

Steve Nixon

Director - Viamed Ltd.

Office: [+44 \(0\)1535 634542](tel:+44201535634542)

Mobile: [+44 \(0\)7850 252267](tel:+442017850252267)

--

Steve

Steve Nixon

Director - Viamed Ltd.

Office: [+44 \(0\)1535 634542](tel:+44201535634542)

Mobile: [+44 \(0\)7850 252267](tel:+442017850252267)

01/05/2018

Gmail - Fwd: Polymer resin oxygen sensor

--

Steve

Steve Nixon

Director - Viamed Ltd.

Office: [+44 \(0\)1535 634542](tel:+44201535634542)

Mobile: [+44 \(0\)7850 252267](tel:+44207850252267)

--

Steve

Steve Nixon

Director - Viamed Ltd.

Office: [+44 \(0\)1535 634542](tel:+44201535634542)

Mobile: [+44 \(0\)7850 252267](tel:+44207850252267)

--

Steve

Steve Nixon

Director - Viamed Ltd.

Office: [+44 \(0\)1535 634542](tel:+44201535634542)

Mobile: [+44 \(0\)7850 252267](tel:+44207850252267)

--

Steve

Steve Nixon

Director - Viamed Ltd.

Office: [+44 \(0\)1535 634542](tel:+44201535634542)

Mobile: [+44 \(0\)7850 252267](tel:+44207850252267)