

NeoPEEP Neonatal Resuscitation Circuits – December 2017

Opera part number	Description	Notes	Status
3210010	Pink, 1.0m. Single		Discontinued.
3210011	Blue, 1.0m. Single	Ordered as singles, supplied in boxes of 20.	Retained as current variant. F&P adapter retained in the packaging.
3210012	Pink, 1.0m. Box of 20		Discontinued. No need for colour variants, eases stock control. Minimizes/clarifies options for customers.
3210013	Blue, 1.0m. Box of 20		Retained as current variant. F&P adapter retained in the packaging.
3210014	Blue, 1.0m. Box of 20	Tom Thumb adapter included	To discontinue. In accordance with CE not allowed to add our adapter.
3210020	Pink, 1.5m. Single	Ordered as singles, supplied in boxes of 20.	Discontinued. No need for colour variants, eases stock control. Minimizes/clarifies options for customers.
3210021	Blue, 1.5m. Single	Ordered as singles, supplied in boxes of 20.	Retained as current variant. F&P adapter retained in the packaging.
3210022	Pink, 1.5m. Box of 20		Discontinued. No need for colour variants, eases stock control. Minimizes/clarifies options for customers.
3210023	Blue, 1.5m. Box of 20		Retained as current variant. F&P adapter retained in the packaging.

Historically the NeoPEEPs were sold with a soft silicone adapter to be used to connect the NeoPEEP to Fisher & Paykel NeoPuff neonatal patient circuit. We removed this as not for use with our Tom Thumb units. In addition this adapter has a sampling port which is supplied open.

In accordance with current CE we are not allowed to remove this adapter, so it will be included as standard. I have asked the manufacturer Marshall Products to ensure that the adapter is supplied in the closed position (to minimize clinical risk) and they have agreed that this would be good practice.

Also note that historically we used to unlock the PEEP valve. Again we are not allowed to alter the product, but this is not an issue as the NeoPEEP circuits are now supplied set to 5cm H₂O.

NeoPEEP circuit used without adapter – connects directly to standard resuscitators.

NeoPEEP circuit used with included adapter – connects to Fisher & Paykel NeoPuff

We also have an adapter that can be used with the NeoPuff.

We also have adapters (single use & reusable) that can be used to connect the NeoPEEP circuit to Viamed Tom Thumb Infant Resuscitators.