

«TITLE» «FORENAME» «SURNAME»  
«JOB\_TITLE»  
«ORGANISATION»  
«ADDRESS\_1»  
«ADDRESS\_2»  
«ADDRESS\_3»  
«CITY»  
«POST\_CODE»

April 2018

Dear «TITLE» «SURNAME»,

**Introducing VersaStream Gas Sampling Lines  
compatible for use with a variety of Capnographs and Multigas Monitors**

Viamed are pleased to introduce VersaStream gas sampling lines as an innovative solution for sidestream gas monitoring, offering a one-piece sampling line solution with increased performance and low cost.

VersaStream gas sampling lines use the latest material technologies to allow the removal of moisture and water - even in cases of low temperatures and high humidity. They are very flexible and are suitable for both short-term and long-term applications.

**Oridion**

Nellcor, Covidien,  
Medtronics

Compatible to the Filterline range  
of sampling lines and for use with  
Microstream enabled capnographs



**Respironics**

Philips, Dräger

For use with  
LoFlo™ sidestream modules



**Luer Lock Male**

For use with Capnographs  
and Multigas Monitors that  
incorporate a Luer Lock  
sampling port



**Viamed**

Masimo, Phase-In

For use with VM-2500 Sidestream  
Capnographs and  
Masimo/Phase-In technology  
enabled capnographs



**Why VersaStream gas sampling lines?**

- One-piece sampling line solution
- **Cost-effective** alternative to Nafion® tubes
- The gas-selective non-porous material facilitates fast and removal of condensation from the sampling line
- The integrated bacterial filter **prevents contamination**
- **Reliable measured gas values**, due to negligible dead space volume of the filter - even in cases of high respiration rates
- **Integrated bacteria filter and water removal** reduces the number of patient circuit parts and speeds up patient monitoring preparation time

These lines are available in a number of variants including: single nasal sampling line, nasal cannula, nasal/oral cannula and a one piece airway adapter sampling line. They are compatible with a wide range of OEM equipment, please see our cross reference and selection guides for further information: [www.viamed.co.uk/linesxref](http://www.viamed.co.uk/linesxref). Should you be unable to find the equipment you are using please contact us so that we can search our extensive database.



Please find enclosed a leaflet outlining the range along with a price list. Should you have any further questions regarding this, or would like to obtain samples for evaluation, please do not hesitate in contacting our sales team on (01535) 634542 between 9:00am and 5:30pm Monday - Friday.

Yours sincerely



Catrin Hollings  
Sales & Marketing Co-ordinator  
sales@viamed.co.uk