

NomoLine™ Airway Adapter Sets

No-moisture sampling for extended monitoring of intubated patients



NomoLine HH Multi-Patient Use Adapter
and NomoLine Adult/Pediatric Airway
Adapter Set with Luer



NomoLine HH Infant/Neonatal
Airway Adapter Set



NomoLine LH Adult/Pediatric
Airway Adapter Set

Advanced technology enables cost-effective, hassle-free sidestream capnography and gas measurements

- Available in both low-humidity (LH) and high-humidity (HH) configurations for use in a variety of clinical scenarios
- Patented polymer allows water in the sampling line to evaporate into the surrounding air while leaving oxygen, carbon dioxide, and anesthetic gas measurements unaffected, eliminating the need for a water trap and issues related to their handling
- Enables extended monitoring time in high-humidity applications to reduce disposable costs
- Designed for low sampling flow applications, with functionality in any orientation, supporting use on patients with low tidal volumes and high breath rates in all patient populations, from neonates to adults

Using NomoLine With Masimo Devices

Root capabilities expand with the ISA™ family of modules to include real-time capnography and gas measurements



The NomoLine Airway Adapter Set connects to the ISA module

The ISA module easily mounts to the back of the Root monitor and connects via the Masimo Open Connect™ (MOC-9™) ports

Rad-97™ with NomoLine Capnography is available with an integrated port for direct connection of NomoLine sampling lines



Light Emitting Gas Inlet (LEGI) indicator changes color to indicate the NomoLine sampling line is correctly connected or occluded

NomoLine Specifications

Dead space

Adult/Pediatric	≤ 6 ml
Infant and Infant/Neonatal	≤ 0.7 ml

Dimensions

Sampling line inner diameter	1 mm
Standard length	2 m

Ordering Information

Root and ISA

	PN
Root	9515
ISA CO ₂ Kit (module, bracket)	3727
ISA AX+ Kit (module, bracket)	4376
ISA OR+ Kit (module, bracket)	9685

PN

	PN
NomoLine HH Airway Adapter Set	
Adult/Pediatric, 2 m	3827
Adult/Pediatric, 3 m	3828
Infant, 2 m	3829
Infant/Neonatal, 2 m	4367

PN

	PN
Rad-97 with NomoLine Capnography	9799

	PN
NomoLine HH Multi-Patient Use Adapter	3800

PN

	PN
NomoLine LH Airway Adapter Set	
Adult/Pediatric, 2 m	3814
Adult/Pediatric, 3 m	3815
Infant, 2 m	3816
Infant/Neonatal, 2 m	4369

	PN
NomoLine Airway Adapter Set with Luer	
Adult/Pediatric, 2 m	3801
Adult/Pediatric, 3 m	3802
Infant, 2 m	3803
Infant/Neonatal, 2 m	4368

PN

Neonatal NomoLine products, ISA modules for neonatal use, and Rad-97 are not licensed for sale in Canada.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information including indications, contraindications, warnings, and precautions.