

Sister Catrin Hollings
Senior Sister
Paediatrics Department
Vandagraph Sensor Technologies
15 Station Road
Cross Hills
Keighley, West Yorkshire BD20 7DT

September 2014

Dear Sister Hollings,

Reusable and Single Patient Use Oxygen Hoods

Designed for use with neonatal and paediatric patients - for the short-term, stable delivery of oxygen

Reusable Oxygen Hoods

An oxygen hood is probably the most common traditional short-term method of delivering extra oxygen to neonatal/paediatric patients that require it; maintaining stable FiO_2 around the patient's head, whilst leaving the body accessible for care. These Perspex® reusable oxygen hoods feature a sliding front aperture.

Standard Hoods

This hood can be used in conjunction with Viamed's phototherapy light shields for use during phototherapy.



Hood With Side Ports


Side ports enable the user to have easier access to the patient.



Viamed manufactures a range of oxygen hoods either with or without side ports, and in addition to the range detailed in the enclosed leaflet, we can manufacture oxygen hoods to your specific requirements.

Single Patient Use Oxygen Hoods

In addition to our range of reusable oxygen hoods Viamed provides a range of single patient use oxygen hoods, available in several formats and sizes. These oxygen hoods are designed for single patient use, thereby reducing the risk of cross infection and are flat packed for easy storage and transportation.



Should the above be of interest, please do not hesitate to contact our sales team on **(01535) 634542** or email **info@viamed.co.uk** if you have any queries or require further information.

Yours sincerely

Catrin

Catrin Hollings
Sales & Marketing Co-ordinator