

Hospital: *Royal Preston*

Room: 7

Date: *5/10/20*

Cabinet Body Serial number: *D008*

Check for overall signs of damage Comments ✓
Check door catch - swivel catch on door ✓
Check door catch retaining clasp on cabinet ✓
Check integrity of welds on the door hinges ✓
Check hinges for damage or deformation ✓
Check the restraining bars for damage or deformation ✓
Check the cabinet is securely fixed to wall ✓
Check medirails and tighten any loose rails ✓
Check grommets on the hose inlets are in place and intact ✓
Check the bumper stops that prevent door banging during closing ✓
Check caution sticker is in place on the platform ✓
Check bed for damage i.e. cracks, chips ✓
Check all the bed mounting screws are in place ✓
Check mattress for damage i.e. tears, split seams ✓

Tom Thumb - Serial Number: *VMEF08*

Minimum reading less 7.85 mbar: *1.9*

Adjustable valve 42.1 mbar +/-1: *42.1*

Precision Valve 44.1 mbar +/-1: *44.2*

Gauge @30cmh₂o (28.4-30.4 mbar): *30.1*

Check oxygen hose, if present, for wear, and replace if necessary (at least every 4 years)

Date of last hose replacement:

Low Suction Controller - Serial number: *031144527*

Change filter assembly ✓

Check vacuum hose for wear, and replace if necessary (at least every 4 years) ✓

Date of last hose replacement: *5/10/20*

Check rail clamp is securely attached to suction controller ✓

Replace suction tubing ✓

Replace tapered connector ✓

Check receiving canister for damage, replaced? ✓

Radiant Warmer Serial number: *WD2059RW*

Is Radiant Warmer hard-wired? ✓ *y*

PAT test performed? ✓

APGAR Timer - Serial number: *00948*

Change batteries ✓

Check storage bins present ✓ - *No, not required*

Tidy any trailing hoses, using tie-wraps if necessary ✓

Check Operator's Manuals are in place ✓

Clean cabinet with isopropyl alcohol ✓

Blender Serial number: *B1551981*

Date of last hose replacement: *5/10/20*

@20.9: @30: @60: @90: @100: (+/-3)

Check Alarm for no Oxygen

Check Alarm for no Air

Check O₂/Air Balance Alarms before output change

Check Duck Bill Valves

Check for Leaks

Blender requires overhaul

Hospital: Royal Preston

Room: 4

Date: 12/10/20

Cabinet Body Serial number: D003

Check for overall signs of damage Comments ✓
Check door catch - swivel catch on door ✓ - slightly loose, recommend replacement
Check door catch retaining clasp on cabinet ✓
Check integrity of welds on the door hinges ✓
Check hinges for damage or deformation ✓
Check the restraining bars for damage or deformation ✓
Check the cabinet is securely fixed to wall ✓
Check medirails and tighten any loose rails ✓
Check grommets on the hose inlets are in place and intact ✓
Check the bumper stops that prevent door banging during closing ✓
Check caution sticker is in place on the platform ✓
Check bed for damage i.e. cracks, chips ✓
Check all the bed mounting screws are in place ✓
Check mattress for damage i.e. tears, split seams ✓

Tom Thumb - Serial Number: VMEFO1

Minimum reading less 7.85 mbar: 20 mbar

Adjustable valve 42.1 mbar +/-1: 41.2

Precision Valve 44.1 mbar +/-1: 44.4

Gauge @30cmh2o (28.4-30.4 mbar): 29.4

Check oxygen hose, if present, for wear, and replace if necessary (at least every 4 years)

Date of last hose replacement:

Low Suction Controller - Serial number: 0311 80205

Change filter assembly ✓

Check vacuum hose for wear, and replace if necessary (at least every 4 years) ✓

Date of last hose replacement: 12/10/20

Check rail clamp is securely attached to suction controller ✓

Replace suction tubing ✓

Replace tapered connector ✓

Check receiving canister for damage, replaced? ✓ - OK

Radiant Warmer Serial number: YD2062 RW

Is Radiant Warmer hard-wired? ✓/✓

PAT test performed? ✓/Y

APGAR Timer - Serial number: N/A - Old black type timer

Change batteries ✓

Check storage bins present x - Not required

Tidy any trailing hoses, using tie-wraps if necessary ✓

Check Operator's Manuals are in place ✓

Clean cabinet with isopropyl alcohol ✓

Blender Serial number:

Date of last hose replacement:

@20.9: @30: @60: @90: @100: (+/-3)

Check Alarm for no Oxygen

Check Alarm for no Air

Check O2/Air Balance Alarms before output change

Check Duck Bill Valves

Check for Leaks

Requires overhaul.

Hospital: Royal Preston

Room: Mat B

Date: 12/10/20

Cabinet Body Serial number: E347

- Check for overall signs of damage Comments ✓
- Check door catch - swivel catch on door ✓
- Check door catch retaining clasp on cabinet ✓
- Check integrity of welds on the door hinges ✓
- Check hinges for damage or deformation ✓
- Check the restraining bars for damage or deformation ✓
- Check the cabinet is securely fixed to wall ✓
- Check medirails and tighten any loose rails ✓
- Check grommets on the hose inlets are in place and intact ✓
- Check the bumper stops that prevent door banging during closing ✓
- Check caution sticker is in place on the platform ✓
- Check bed for damage i.e. cracks, chips ✓
- Check all the bed mounting screws are in place ✓
- Check mattress for damage i.e. tears, split seams ✓

Tom Thumb - Serial Number: VMDI25

Minimum reading less 7.85 mbar: 4.2

Adjustable valve 42.1 mbar +/-1: 41.2

Precision Valve 44.1 mbar +/-1: 44.5

Gauge @30cmh₂o (28.4-30.4 mbar): 30.0

Check oxygen hose, if present, for wear, and replace if necessary (at least every 4 years)

Date of last hose replacement:

Low Suction Controller - Serial number: 031144519

- Change filter assembly ✓
- Check vacuum hose for wear, and replace if necessary (at least every 4 years) ✓
- Date of last hose replacement: 12/10/20
- Check rail clamp is securely attached to suction controller ✓
- Replace suction tubing ✓
- Replace tapered connector ✓
- Check receiving canister for damage, replaced? ✓ OK

Radiant Warmer Serial number: WD2038 RW

Is Radiant Warmer hard-wired? ✓ Y

PAT test performed? ✓ Y

APGAR Timer - Serial number: -N/A - old style timer

- Change batteries ✓
- Check storage bins present ✓ XPS
- Tidy any trailing hoses, using tie-wraps if necessary ✓
- Check Operator's Manuals are in place ✓
- Clean cabinet with isopropyl alcohol ✓

Blender Serial number: B1551980 - Brought back for overhaul

Date of last hose replacement: 12/10/20

@20.9: @30: @60: @90: @100: (+/-3)

- Check Alarm for no Oxygen
- Check Alarm for no Air
- Check O2/Air Balance Alarms before output change
- Check Duck Bill Valves
- Check for Leaks

Installed blender B1551978 (just overhauled)

Hospital: Royal Preston

Room: BIRTH CENTRE Date: 12/10/20

Cabinet Body Serial number: E357

- Check for overall signs of damage Comments ✓
- Check door catch - swivel catch on door ✓
- Check door catch retaining clasp on cabinet ✓
- Check integrity of welds on the door hinges ✓
- Check hinges for damage or deformation ✓
- Check the restraining bars for damage or deformation ✓
- Check the cabinet is securely fixed to wall ✓
- Check medirails and tighten any loose rails ✓
- Check grommets on the hose inlets are in place and intact ✓
- Check the bumper stops that prevent door banging during closing ✓
- Check caution sticker is in place on the platform ✓
- Check bed for damage i.e. cracks, chips ✓ 14"
- Check all the bed mounting screws are in place ✓
- Check mattress for damage i.e. tears, split seams ✓ 14"

Tom Thumb - Serial Number: 0401172

Minimum reading less 7.85 mbar: 1.5

Adjustable valve 42.1 mbar +/-1: 42.5

Precision Valve 44.1 mbar +/-1: 43.6

Gauge @30cmh₂o (28.4-30.4 mbar): 29.5

Check oxygen hose, if present, for wear, and replace if necessary (at least every 4 years) ✓

Date of last hose replacement: ~~25/3/19~~ N/A (TT480)

Low Suction Controller - Serial number: 25370002

Change filter assembly ✓

Check vacuum hose for wear, and replace if necessary (at least every 4 years) ✓

Date of last hose replacement: Check service logs - 25/3/19

Check rail clamp is securely attached to suction controller ✓

Replace suction tubing ✓

Replace tapered connector ✓

Check receiving canister for damage, replaced? ✓

Radiant Warmer Serial number: 475-3-0314

Is Radiant Warmer hard-wired? ✓ y

PAT test performed? ✓

APGAR Timer - Serial number: 00740

Change batteries ✓

Check storage bins present ✓ y

Tidy any trailing hoses, using tie-wraps if necessary ✓

Check Operator's Manuals are in place ✓

Clean cabinet with isopropyl alcohol ✓

Blender Serial number: B1404670

Date of last hose replacement: Check service logs - 25/3/19

@20.9: 25 @30: @60: @90: @100: 97 (+/-3)

Check Alarm for no Oxygen

Check Alarm for no Air

Check O2/Air Balance Alarms before output change

Check Duck Bill Valves

Check for Leaks

Needs flowmeter outlet. - Sorted locally by EBME Dept. ✓

Blender not calibrated, only has range from 25% - 97%

Installed Blender s/n B1551983 (just overhauled) to allow us to bring back for cal. Note: this one was only due an annual service, parts have been replaced already.

Hospital: ROYAL PRESTON

Room: 3

Date: 14-10-20

Cabinet Body Serial number: D007

Check for overall signs of damage Comments ✓
Check door catch - swivel catch on door ✓
Check door catch retaining clasp on cabinet ✓
Check integrity of welds on the door hinges ✓
Check hinges for damage or deformation ✓
Check the restraining bars for damage or deformation ✓ — FAILING WELD photo taken
Check the cabinet is securely fixed to wall ✓
Check medirails and tighten any loose rails ✓
Check grommets on the hose inlets are in place and intact ✓
Check the bumper stops that prevent door banging during closing ✓
Check caution sticker is in place on the platform ✓
Check bed for damage i.e. cracks, chips ✓
Check all the bed mounting screws are in place ✓
Check mattress for damage i.e. tears, split seams ✓

Tom Thumb - Serial Number: VMEF03 14072

Minimum reading less 7.85 mbar: ~~22.0~~ 0.4

Adjustable valve 42.1 mbar +/-1: 41.5

Precision Valve 44.1 mbar +/-1: 44.0

Gauge @30cmh2o (28.4-30.4 mbar): 29.9

Check oxygen hose, if present, for wear, and replace if necessary (at least every 4 years)

Date of last hose replacement:

Low Suction Controller - Serial number: 031144509

Change filter assembly ✓

Check vacuum hose for wear, and replace if necessary (at least every 4 years) ✓

Date of last hose replacement:

Check rail clamp is securely attached to suction controller ✓

Replace suction tubing ✓

Replace tapered connector ✓

Check receiving canister for damage, replaced? OK

Radiant Warmer Serial number: YD2060 RW 14071

Is Radiant Warmer hard-wired? NO

PAT test performed? YES

APGAR Timer - Serial number:

Change batteries

Check storage bins present x

Tidy any trailing hoses, using tie-wraps if necessary ✓

Check Operator's Manuals are in place ✓

Clean cabinet with isopropyl alcohol

Blender Serial number: B1551983

55604

Date of last hose replacement:

@20.9: @30: @60: @90: @100: (+/-3)

Check Alarm for no Oxygen

Check Alarm for no Air

Check O2/Air Balance Alarms before output change

Check Duck Bill Valves

Check for Leaks

Hospital: ROYAL PRESTON

Room: 2

Date: 1/10/20

Cabinet Body Serial number: D004

Check for overall signs of damage Comments - signs of wear, chipped paint
Check door catch - swivel catch on door ✓
Check door catch retaining clasp on cabinet ✓
Check integrity of welds on the door hinges ✓ - No cracks
Check hinges for damage or deformation - Slight bending, nothing significant
Check the restraining bars for damage or deformation - ✓
Check the cabinet is securely fixed to wall ✓
Check medirails and tighten any loose rails ✓
Check grommets on the hose inlets are in place and intact ✓
Check the bumper stops that prevent door banging during closing ✓
Check caution sticker is in place on the platform ✓ replaced.
Check bed for damage i.e. cracks, chips ✓ 12"
Check all the bed mounting screws are in place ✓
Check mattress for damage i.e. tears, split seams ✓ 12"

Tom Thumb - Serial Number: ~~D00~~ VMEF02

Minimum reading less 7.85 mbar: 2.1

Adjustable valve 42.1 mbar +/-1: 41.6

Precision Valve 44.1 mbar +/-1: 43.7

Gauge @30cmh2o (28.4-30.4 mbar): 30.4

Check oxygen hose, if present, for wear, and replace if necessary (at least every 4 years)

Date of last hose replacement:

Low Suction Controller - Serial number: 031144520 - Cracked gauge face
Change filter assembly ✓ Advised engineer.

Check vacuum hose for wear, and replace if necessary (at least every 4 years) ✓

Date of last hose replacement: 1/10/20

Check rail clamp is securely attached to suction controller ✓

Replace suction tubing ✓

Replace tapered connector ✓

Check receiving canister for damage, replaced? ✓ -ok

Radiant Warmer Serial number: YD2064RW

Is Radiant Warmer hard-wired? Y

PAT test performed? Y

- Cracked case by lamp, advised inspection by EBME.

APGAR Timer - Serial number: None - Old style black timer

Change batteries ✓

Check storage bins present - No

Tidy any trailing hoses, using tie-wraps if necessary ✓

Check Operator's Manuals are in place ✓

Clean cabinet with isopropyl alcohol ✓

Blender Serial number: B1551978

Date of last hose replacement: 1/10/20

@20.9: @30: @60: @90: @100: (+/-3)

Check Alarm for no Oxygen

Check Alarm for no Air

Check O2/Air Balance Alarms before output change

Check Duck Bill Valves

Check for Leaks

Returned for overhaul service. SRS67608
Service exchange blender installed

Hospital: Royal Preston

Room: 6

Date: 1/10/20

Cabinet Body Serial number: D006

- Check for overall signs of damage Comments ✓
- Check door catch - swivel catch on door ✓
- Check door catch retaining clasp on cabinet ✓
- Check integrity of welds on the door hinges ✓
- Check hinges for damage or deformation ✓
- Check the restraining bars for damage or deformation ✓
- Check the cabinet is securely fixed to wall ✓
- Check medirails and tighten any loose rails ✓
- Check grommets on the hose inlets are in place and intact ✓
- Check the bumper stops that prevent door banging during closing ✓
- Check caution sticker is in place on the platform ✓
- Check bed for damage i.e. cracks, chips ✓
- Check all the bed mounting screws are in place ✓
- Check mattress for damage i.e. tears, split seams ✓

Tom Thumb - Serial Number: VMEF 009

Minimum reading less 7.85 mbar: 2.2

Adjustable valve 42.1 mbar +/-1: 41.3

Precision Valve 44.1 mbar +/-1: 43.4

Gauge @30cmh2o (28.4-30.4 mbar): 28.9

Check oxygen hose, if present, for wear, and replace if necessary (at least every 4 years)

Date of last hose replacement:

Low Suction Controller - Serial number: 031144526

- Change filter assembly ✓
- Check vacuum hose for wear, and replace if necessary (at least every 4 years) ✓
- Date of last hose replacement: 1/10/20
- Check rail clamp is securely attached to suction controller ✓
- Replace suction tubing ✓
- Replace tapered connector ✓
- Check receiving canister for damage, replaced? ✓

Radiant Warmer Serial number: WD2042RW

Is Radiant Warmer hard-wired? ✓ Y

PAT test performed? ✓ Y

APGAR Timer - Serial number: None. Old style black timer

- Change batteries ✓
- Check storage bins present - None
- Tidy any trailing hoses, using tie-wraps if necessary ✓
- Check Operator's Manuals are in place ✓
- Clean cabinet with isopropyl alcohol ✓

Blender Serial number: B1551985

Date of last hose replacement: 1/10/20

@20.9: @30: @60: @90: @100: (+/-3)

- Check Alarm for no Oxygen
- Check Alarm for no Air
- Check O2/Air Balance Alarms before output change
- Check Duck Bill Valves
- Check for Leaks

Requires Overhaul to be returned to Vicamed.

Hospital: Royal Preston

Room: 5

Date: 1/10/20

Cabinet Body Serial number: E346

Check for overall signs of damage Comments ✓ Shuttleworths cabinet (replacement)
Check door catch - swivel catch on door ✓
Check door catch retaining clasp on cabinet ✓
Check integrity of welds on the door hinges ✓
Check hinges for damage or deformation ✓
Check the restraining bars for damage or deformation ✓
Check the cabinet is securely fixed to wall ✓
Check medirails and tighten any loose rails ✓
Check grommets on the hose inlets are in place and intact ✓
Check the bumper stops that prevent door banging during closing - No. Gap wouldn't allow th
Check caution sticker is in place on the platform ✓
Check bed for damage i.e. cracks, chips ✓ 12" bed
Check all the bed mounting screws are in place ✓
Check mattress for damage i.e. tears, split seams 12" mattress

Tom Thumb - Serial Number: VMDI23

Minimum reading less 7.85 mbar: 2.1

Adjustable valve 42.1 mbar +/-1: 42.5

Precision Valve 44.1 mbar +/-1: 44.9

Gauge @30cmh2o (28.4-30.4 mbar): 28.7

Check oxygen hose, if present, for wear, and replace if necessary (at least every 4 years) - Replace

Date of last hose replacement:

Low Suction Controller - Serial number: 031144522

Change filter assembly ✓

Check vacuum hose for wear, and replace if necessary (at least every 4 years) ✓

Date of last hose replacement: 1/10/20

Check rail clamp is securely attached to suction controller ✓

Replace suction tubing ✓

Replace tapered connector ✓

Check receiving canister for damage, replaced? ✓ - ok

Radiant Warmer Serial number: YD2068RW

Is Radiant Warmer hard-wired? ✓

PAT test performed? ✓

APGAR Timer - Serial number: 00947

Change batteries ✓

Check storage bins present ✗ - No

Tidy any trailing hoses, using tie-wraps if necessary ✓

Check Operator's Manuals are in place ✓

Clean cabinet with isopropyl alcohol ✓

Blender Serial number: B1551983

Date of last hose replacement: 1/10/20

@20.9: @30: @60: @90: @100: (+/-3)

Check Alarm for no Oxygen

Check Alarm for no Air

Check O2/Air Balance Alarms before output change

Check Duck Bill Valves

Check for Leaks

Service exchange Blender installed

Returned for overhaul service - SRS67608

Hospital: *Royal Preston*

Room: *69*

Date: *1/10/20*

Cabinet Body Serial number: *DOOC*

Check for overall signs of damage Comments ✓
Check door catch - swivel catch on door ✓
Check door catch retaining clasp on cabinet ✓
Check integrity of welds on the door hinges ✓
Check hinges for damage or deformation ✓
Check the restraining bars for damage or deformation ✓
Check the cabinet is securely fixed to wall ✓
Check medirails and tighten any loose rails ✓
Check grommets on the hose inlets are in place and intact ✓
Check the bumper stops that prevent door banging during closing ✓
Check caution sticker is in place on the platform ✓
Check bed for damage i.e. cracks, chips ✓ *12"*
Check all the bed mounting screws are in place ✓
Check mattress for damage i.e. tears, split seams ✓ *12"*

Tom Thumb - Serial Number: *VMDI22*

Minimum reading less 7.85 mbar: *2.6*

Adjustable valve 42.1 mbar +/-1: *41.7*

Precision Valve 44.1 mbar +/-1: *43.7*

Gauge @30cmh2o (28.4-30.4 mbar): *30.2*

Check oxygen hose, if present, for wear, and replace if necessary (at least every 4 years)

Date of last hose replacement:

Low Suction Controller - Serial number: *031163907*

Change filter assembly ✓
Check vacuum hose for wear, and replace if necessary (at least every 4 years) ✓
Date of last hose replacement: *1/10/20*
Check rail clamp is securely attached to suction controller ✓
Replace suction tubing ✓
Replace tapered connector ✓
Check receiving canister for damage, replaced? ✓ *OK*

Radiant Warmer Serial number: *YD2066 RW*

Is Radiant Warmer hard-wired? *Y*

PAT test performed? *Y*

APGAR Timer - Serial number: *-N/a old black timer type*

Change batteries ✓
Check storage bins present *-No*
Tidy any trailing hoses, using tie-wraps if necessary ✓
Check Operator's Manuals are in place ✓
Clean cabinet with isopropyl alcohol ✓

Blender Serial number: *B1551977*

Date of last hose replacement: *1/10/20*

@20.9: @30: @60: @90: @100: (+/-3)

Check Alarm for no Oxygen
Check Alarm for no Air
Check O2/Air Balance Alarms before output change
Check Duck Bill Valves
Check for Leaks

Blender to be sent back for overhaul.

Hospital: ROYAL PRESTON

Room: 1

Date: 5/10/20

Cabinet Body Serial number: E337

- Check for overall signs of damage Comments ✓
- Check door catch - swivel catch on door ✓
- Check door catch retaining clasp on cabinet ✓
- Check integrity of welds on the door hinges ✓
- Check hinges for damage or deformation ✓
- Check the restraining bars for damage or deformation ✓
- Check the cabinet is securely fixed to wall ✓
- Check medirails and tighten any loose rails ✓
- Check grommets on the hose inlets are in place and intact ✓
- Check the bumper stops that prevent door banging during closing ✓
- Check caution sticker is in place on the platform ✓
- Check bed for damage i.e. cracks, chips ✓
- Check all the bed mounting screws are in place ✓
- Check mattress for damage i.e. tears, split seams ✓

Tom Thumb - Serial Number: VMDI24

Minimum reading less 7.85 mbar: 4.4

Adjustable valve 42.1 mbar +/-1: 42.3

Precision Valve 44.1 mbar +/-1: 44.6

Gauge @30cmh2o (28.4-30.4 mbar): 28.7

Check oxygen hose, if present, for wear, and replace if necessary (at least every 4 years)

Date of last hose replacement: 4

Low Suction Controller - Serial number: 16431002

- Change filter assembly ✓
- Check vacuum hose for wear, and replace if necessary (at least every 4 years) ✓
- Date of last hose replacement: 5/10/20
- Check rail clamp is securely attached to suction controller ✓
- Replace suction tubing ✓
- Replace tapered connector ✓
- Check receiving canister for damage, replaced? ✓ OK

Radiant Warmer Serial number: - YD 2063 RW

Is Radiant Warmer hard-wired? ✓

PAT test performed? ✓

APGAR Timer - Serial number: - N/A, old black timer

- Change batteries ✓
- Check storage bins present - No
- Tidy any trailing hoses, using tie-wraps if necessary ✓
- Check Operator's Manuals are in place ✓
- Clean cabinet with isopropyl alcohol ✓

Blender Serial number: B1551986

Date of last hose replacement: 5/10/20

@20.9: @30: @60: @90: @100: (+/-3)

- Check Alarm for no Oxygen
- Check Alarm for no Air
- Check O2/Air Balance Alarms before output change
- Check Duck Bill Valves
- Check for Leaks

Blender Requires Overhaul

Hospital: Royal Preston

Room: 10

Date: 5/10/20

Cabinet Body Serial number: D002

- Check for overall signs of damage Comments ✓
- Check door catch - swivel catch on door ✓
- Check door catch retaining clasp on cabinet ✓
- Check integrity of welds on the door hinges ✓
- Check hinges for damage or deformation ✓
- Check the restraining bars for damage or deformation ✓
- Check the cabinet is securely fixed to wall ✓
- Check medirails and tighten any loose rails ✓
- Check grommets on the hose inlets are in place and intact ✓
- Check the bumper stops that prevent door banging during closing ✓
- Check caution sticker is in place on the platform ✓
- Check bed for damage i.e. cracks, chips ✓
- Check all the bed mounting screws are in place ✓
- Check mattress for damage i.e. tears, split seams ✓

Tom Thumb - Serial Number: VMEF05

Minimum reading less 7.85 mbar: 1.7

Adjustable valve 42.1 mbar +/-1: 42.0

Precision Valve 44.1 mbar +/-1: 43.3

Gauge @30cmh2o (28.4-30.4 mbar): 29.5

Check oxygen hose, if present, for wear, and replace if necessary (at least every 4 years)

Date of last hose replacement:

Low Suction Controller - Serial number: 031144530

- Change filter assembly ✓
- Check vacuum hose for wear, and replace if necessary (at least every 4 years) ✓
- Date of last hose replacement: 5/10/20
- Check rail clamp is securely attached to suction controller ✓
- Replace suction tubing ✓
- Replace tapered connector ✓
- Check receiving canister for damage, replaced? ✓

Radiant Warmer Serial number: YD2061RW

Is Radiant Warmer hard-wired? ✓✓

PAT test performed? ✓

APGAR Timer - Serial number: - N/A old black style time

- Change batteries ✓
- Check storage bins present ✗ - No
- Tidy any trailing hoses, using tie-wraps if necessary ✓
- Check Operator's Manuals are in place ✓
- Clean cabinet with isopropyl alcohol ✓

Blender Serial number: B1551982

Date of last hose replacement: 5/10/20

@20.9: @30: @60: @90: @100: (+/-3)

- Check Alarm for no Oxygen
- Check Alarm for no Air
- Check O2/Air Balance Alarms before output change
- Check Duck Bill Valves
- Check for Leaks

Requires overhaul

Hospital: Royal Preston

Room: 11

Date: 5/10/20

Cabinet Body Serial number: DO10

Check for overall signs of damage ✓
Check door catch - swivel catch on door ✓
Check door catch retaining clasp on cabinet ✓
Check integrity of welds on the door hinges ✓
Check hinges for damage or deformation ✓
Check the restraining bars for damage or deformation ✓
Check the cabinet is securely fixed to wall ✓
Check medirails and tighten any loose rails ✓
Check grommets on the hose inlets are in place and intact ✓
Check the bumper stops that prevent door banging during closing ✓
Check caution sticker is in place on the platform ✓
Check bed for damage i.e. cracks, chips ✓
Check all the bed mounting screws are in place ✓
Check mattress for damage i.e. tears, split seams ✓

Tom Thumb - Serial Number: VMDI21

Minimum reading less 7.85 mbar: 1.3

Adjustable valve 42.1 mbar +/-1: 41.5

Precision Valve 44.1 mbar +/-1: 43.6

Gauge @30cmh2o (28.4-30.4 mbar): 28.4

Check oxygen hose, if present, for wear, and replace if necessary (at least every 4 years)

Date of last hose replacement:

Low Suction Controller - Serial number: 031144510

Change filter assembly ✓

Check vacuum hose for wear, and replace if necessary (at least every 4 years) ✓

Date of last hose replacement: 5/10/20

Check rail clamp is securely attached to suction controller ✓

Replace suction tubing ✓

Replace tapered connector ✓

Check receiving canister for damage, replaced?

Radiant Warmer Serial number: YD206SRW

Is Radiant Warmer hard-wired? ✓

PAT test performed? ✓

APGAR Timer - Serial number: 0091644510

Change batteries ✓

Check storage bins present ✓ N

Tidy any trailing hoses, using tie-wraps if necessary ✓

Check Operator's Manuals are in place ✓

Clean cabinet with isopropyl alcohol ✓

Blender Serial number: B1551984

Date of last hose replacement: 5/10/20

@20.9: @30: @60: @90: @100: (+/-3)

Check Alarm for no Oxygen

Check Alarm for no Air

Check O2/Air Balance Alarms before output change

Check Duck Bill Valves

Check for Leaks

Requires overhaul

Hospital: *Royal Preston*

Room: 12

Date: *5/10/20*

Cabinet Body Serial number: *D100*

- Check for overall signs of damage Comments ✓
- Check door catch - swivel catch on door ✓
- Check door catch retaining clasp on cabinet ✓
- Check integrity of welds on the door hinges ✓
- Check hinges for damage or deformation ✓
- Check the restraining bars for damage or deformation ✓
- Check the cabinet is securely fixed to wall ✓
- Check medirails and tighten any loose rails ✓
- Check grommets on the hose inlets are in place and intact ✓
- Check the bumper stops that prevent door banging during closing ✓
- Check caution sticker is in place on the platform ✓
- Check bed for damage i.e. cracks, chips ✓
- Check all the bed mounting screws are in place ✓
- Check mattress for damage i.e. tears, split seams ✓

Tom Thumb - Serial Number: *VMEF04*

Minimum reading less 7.85 mbar: *3.1*

Adjustable valve 42.1 mbar +/-1: *42.3*

Precision Valve 44.1 mbar +/-1: *43.8*

Gauge @30cmh2o (28.4-30.4 mbar): *29.9*

Check oxygen hose, if present, for wear, and replace if necessary (at least every 4 years)

Date of last hose replacement:

Low Suction Controller - Serial number: *031144525*

- Change filter assembly ✓
- Check vacuum hose for wear, and replace if necessary (at least every 4 years) ✓
- Date of last hose replacement: *5/10/20*
- Check rail clamp is securely attached to suction controller ✓
- Replace suction tubing ✓
- Replace tapered connector ✓
- Check receiving canister for damage, replaced? ✓ *No, ok.*

Radiant Warmer Serial number: *ZW2198ZW*

- Is Radiant Warmer hard-wired? ✓
- PAT test performed? ✓

APGAR Timer - Serial number: *-N/A - old style timer*

- Change batteries ✓
- Check storage bins present *X No - Not required*
- Tidy any trailing hoses, using tie-wraps if necessary ✓
- Check Operator's Manuals are in place ✓
- Clean cabinet with isopropyl alcohol ✓

Blender Serial number: *B1551971*

Date of last hose replacement: *05/10/20*

@20.9: @30: @60: @90: @100: (+/-3)

- Check Alarm for no Oxygen
- Check Alarm for no Air
- Check O2/Air Balance Alarms before output change
- Check Duck Bill Valves
- Check for Leaks

Requires over-haul.