

Supplies Department  
Nuffield Hospital Derby  
Rykneld Road  
Littleover  
Derby  
DE23 4SN

Ref: 20191113/0191235/544

13 November 2019

Dear Sir or Madam,

**BlenderBuddy 2**  
**A compact and clean addition to your blender**

Viamed is pleased to introduce the **BlenderBuddy 2** to their Blender and accessories range. The **BlenderBuddy 2** is a flowmeter intended for use with an air/oxygen blender and is to be used by physicians, respiratory therapists and other authorised hospital personnel to administer selected doses of medical gases to a patient.

The **BlenderBuddy 2** utilizes a Maxtec Designed for Blenders (DFB) flowmeter which is designed and calibrated specifically for air/oxygen blenders to provide increased accuracy compared to standard flowmeters. By utilizing dual-scale graduations, the flowmeter provides two flowmeters in one for increased accuracy at lower flow rates. The flowmeter also features a high-quality acrylic block body and precision valve. It provides a sensor analysis port which may be used with a suitable oxygen analyzer to measure the gas concentration from the blender. A small continuous sensor bleed is provided which enables gas analysis regardless of which outlet port is used on the blender.



The **BlenderBuddy 2** can be supplied on its own or with a MAXO2+A Oxygen Analyser.

**Part Numbers**

BlenderBuddy 2, DFB, Low Flow 0 - 15 lpm

0310211

BlenderBuddy 2, DFB with MaxO2+A Oxygen Analyser, Low Flow 0 - 15 lpm

0310216

Please find enclosed leaflet on the **Maxtec BlenderBuddy 2** for further information, which I believe will be of interest to you.

Should you have any further questions regarding this, please do not hesitate in contacting our sales team on 01535 634542 between 9:00am and 5:30pm Monday – Friday or email [blenderbuddy@viamed.co.uk](mailto:blenderbuddy@viamed.co.uk)

Yours sincerely

*Catrin*

Catrin Hollings  
Marketing

Should you wish to unsubscribe from our mailings  
please email [unsubscribe@viamed.co.uk](mailto:unsubscribe@viamed.co.uk)