





Single Patient Use Oxygen Hoods


Designed for use with neonatal and paediatric patients, for the short-term, stable delivery of oxygen




**Standard
Oxygen Hood**



**No Bottom Panel
Oxygen Hood**



**Partial Bottom Panel
Oxygen Hood**



**NICU
Oxygen Hood**

These Oxygen Hoods consist of flat sheets of transparent vinyl with built-in support frames and a drape neck opening.

- See-through environment provides full patient visibility.
- Drape neck opening provides easy access to the patient.
- Single patient use means reduced risk of cross-infection.
- Flat packed for easy storage and transportation.

Standard Hood

The standard Oxygen Hood has two 22 mm inlet ports, a full vinyl bottom, a drape neck opening and is available in a variety of sizes.


Other variants include:

- Partial bottom - remains stable whilst reducing heat build-up under the hood
- No-bottom - provides easier access to the patient


NICU

The NICU single patient use hoods are bottomless with a 75 mm x 125 mm draped access port on the rear panel, plus a drape to the front. Offers maximum patient access for free flow (blow-by) procedures or adjustments.


Single Patient Use Oxygen Hoods




Standard



Partial bottom



No bottom



NICU

Dimensions (mm) L x W x H	Standard 22 mm inlet on top and rear panels, drape to front, fixed bottom panel	Partial bottom 22 mm inlet on top and rear panels, drape to the front, partial bottom panel	No bottom 22 mm inlet on top and rear panels, drape to the front, no bottom panel	NICU 22 mm inlet on top and rear panels, drape to the front, no bottom panel. 75 mm x 125 mm draped access port on rear panel
178 x 178 x 127			2310025	
178 x 178 x 178	2310000		2310026	2310075
216 x 216 x 216	2310001			2310076
305 x 305 x 254			2310031	2310080
305 x 305 x 305		2310050		2310081
457 x 457 x 457		2310052		
610 x 610 x 610	2310013	2310055		

Specifications Subject to Change.

Distributed by:



VIAMED

Viamed Limited · 15 Station Road · Cross Hills
Keighley · West Yorkshire · BD20 7DT · United Kingdom
Tel: +44 (0)1535 634 542 Fax: +44 (0)1535 635 582
Email: info@viamed.co.uk Website: www.viamed.co.uk