

Dispatch date 14.04.20  
 Our Reference NMDM  
 Your Reference  
 Delivery mode UPS  
 Mandate Rep. 2 x Radiant warmer 600-3  
 Cust. number 00158  
 Special receiver Steve Hardaker  
 Invoice Addr. Viamed Limited

Viamed Limited  
 15 Station Road  
 Steve Hardaker  
 Cross Hills  
 BD20 7DT Keighley, W. Yorks - GB

## DELIVERY NOTE 2000649

Page: 001

For order 2000237

| Article Number | Article description | Quant. | Unit |
|----------------|---------------------|--------|------|
|----------------|---------------------|--------|------|

RMA:28022020 / Viamed RN: RTS1100719

Radiant warmer Ceratherm 600-3

S/No.: 155330218

Inventory Number: 294070

Failure: Continuous alarm tone. Alarm cannot be silenced. Alarm does not reset upon restoring power.

Radiant warmer Ceratherm 600-3

S/No.: 182231018

Inventory Number: 294066

Failure: not restoring previous settings even after 2 seconds with power off.

### Materials:

|        |               |   |   |
|--------|---------------|---|---|
| 099-1  | Working hours | 2 | 1 |
| 099-1K | Goodwill      | 2 | 1 |

### Remarks:

Continuous alarm / power failure alarm for both devices.

Device with S/No.: 15533218 replaced the CPU-print and terminal X2 on power-print.

Device with S/No.: 182231018 error could not be reproduced, CPU-Print was replaced for safety reasons.

Function check, safety check, test run, cleaning, administrative work.

All work and material FOC.

Performed by: 5726

Kind Regards  
 Nufer Medical AG

Weight

Charges:

Package: