

Department Manager
«DEPARTMENT»
«COMPANY»
«ADDRESS1»
«ADDRESS2»
«ADDRESS3»
«CITY» «STATE»
«ZIP»

Dear Sir or Madam,

AlcoTrue M Breathalyser provides rapid and accurate determination of breath alcohol concentration



The AlcoTrue M is a professional handheld medical breathalyser for accurate determination of exhaled breath alcohol concentration which rivals its competitors and closely matches evidential devices.

The AlcoTrue M provides midwives and health professionals with quick and accurate determination of breath alcohol of patients, when there are concerns regarding blood alcohol level or alcohol consumption.

The main benefits include:

- **Fast warm-up time** < 6 seconds; allowing for quick and easy diagnosis.
- **Rapid, precise and reliable clinically accurate measurements**
- **Easy to read display.**
- This high spec product provides up to five measurement modes, ideal for any situation;
 - **Active mode** – To determine the precise alcohol level in the end-expiratory breath of the person to be tested after achieving a minimal breath volume from 1.2 L. It is ensured that only deep lung breath (alveolar breath) is analysed.
 - **Fast mode** – Intended for persons suffering from COPD, Asthma or reduced lung capacity, whom might not be able to provide the standard minimum sample volume.
 - **Passive mode*** - To obtain a precise measurement, the sample uptake must be triggered when the person is breathing out. The measurement mode is to determine the quantitative value of alcohol concentration.
 - **Scan mode*** – Qualitative testing of the ambient air alcohol content, qualitative determination of alcohol content in exhaled air of a person who is not able to provide a breath sample (e.g. unconscious person), qualitative test of “non-alcoholic” drinks for alcohol content.
 - **Mouth alcohol mode*** - Consists of two tests with a determined waiting period between first and second measurement.
- Also available with a portable printing system.

* Measurement modes only available if activated.



Please find enclosed further information regarding the AlcoTrue M. Should you wish to discuss this in further detail please do not hesitate to contact a member of our sales team on 01535 634542 or by emailing alcotrue@viamed.co.uk.

Yours faithfully

Catrin

Catrin Hollings
Marketing

Should you no longer wish to receive marketing material via post, sign-up to our email mailing list at www.viamed.co.uk/sign-up