







# 4000 Series Disposable SpO<sub>2</sub> Sensors

Minimising the risk of cross infection



- **Patient types available**  
Adult, Paediatric, Infant and Neonatal
- Plaster tape
- **Cable lengths available**  
Adult & Paediatric - 0.45m (18")  
Infant & Neonatal - 0.9m (3')
- **Compatibles available for**  
Nellcor, BCI, Nonin and Viamed equipment
- Bio-compatible and latex free to help prevent allergic reactions
- Durable and reliable design



For accurate monitoring of oxygen saturation and pulse rate. Compatible for use with a range of original equipment manufacturers' pulse oximeters and patient monitors.

# 4000 Series Disposable SpO<sub>2</sub> Sensors

For use with original equipment manufacturers pulse oximeters and patient monitors including **Nellcor\***, **BCI**, **Nonin** and **Viamed**

\*Non-Oximax technology. All registered trademarks recognised.



Adult



Paediatric



Infant



Neonatal

## Model and Part Numbers

Patient Type	Nellcor*	BCI	Nonin	Viamed**
<b>Adult</b>				
Model Number	1-AP	5-AP	7-AP	10-AP
Part Number	0014909	0014957	0014981	0015010
<b>Paediatric</b>				
Model Number	1-PP	5-PP	7-PP	10-PP
Part Number	0014910	0014958	0014982	0015011
<b>Infant</b>				
Model Number	1-IP	5-IP	7-IP	10-IP
Part Number	0014911	0014959	0014983	0015012
<b>Neonatal</b>				
Model Number	1-NP	5-NP	7-NP	10-NP
Part Number	0014912***	0014960***	0014984***	0015013

\*Non-Oximax technology. Compatible to Oxisensor II range of adhesive sensors.

\*\* Now incorporating SMARTsat SpO<sub>2</sub> technology

\*\*\* Available whilst stocks last only

All disposable SpO<sub>2</sub> sensors are supplied in boxes of 24.

## Environmental Requirements

- **Operating temperature range**  
10°C - 40°C (50°F - 104°F)
- **Storage and transportation temperature range**  
0°C - 50°C (32°F - 122°F)
- **Relative humidity**  
0 - 75% non-condensing

### SpO<sub>2</sub> Accuracy

- 70 - 100% A<sub>RMS</sub> <3.0
- 0 - 69%: Unspecified

Specifications Subject to Change

Distributed by:



# VIAMED

**Viamed Limited** · 15 Station Road · Cross Hills  
Keighley · West Yorkshire · BD20 7DT · United Kingdom  
Tel: +44 (0)1535 634 542 Fax: +44 (0)1535 635 582  
Email: [info@viamed.co.uk](mailto:info@viamed.co.uk) Website: [www.viamed.co.uk](http://www.viamed.co.uk)



Part Number: 0092840  
Date: 03/21