

bluepoint MEDICAL	PHA VersaStream Sampling Lines	Revision 3
Internal Document	Viamed Labeling	Seite 1 von 18

Issued	Signature	Date
Rayk Jürgensen		09.06.2021
Reviewed	Signature	Date

Document History Log

Rev.	Date / By	Modifications
0	13.04.2015 / RJ	Initial version
1	15.08.2018 / AM	Bluepoint as manufacturer
2	14.04.2021 / RJ	Updated label layout
3	09.06.2021 /RJ	Added Oridion versions

Issued: R.Jürgensen Date: 09.06.2021	Approval: R. Jürgensen Date: 09.06.2021
Datei: 4.9_VersaStream_Labeling_Viamed_rev3.docx	




	PHA VersaStream Sampling Lines	Revision 3
Internal Document	Viamed Labeling	Seite 2 von 18

1. Label Layout: VIA Connector Sampling Lines

Dimension Pouch Label: 170 mm x 110 mm
Dimension Box Label: 150 mm x 127 mm

Ref.: 7090121001V, Type: 4420800
Single Port Nasal Prong, Adult, Short Term




Pouch Label




Box Label

Ref.: 7090121002V, Type: 4420801
Nasal Sampling Line, Adult, Short Term




Pouch Label



Box Label

Ref.: 7090121003V, Type: 4420802
Nasal Sampling Line, Paediatric, Short Term




Pouch Label



Box Label

Ref.: 7090121004V, Type: 4420803
Nasal Sampling Line, Infant, Short Term




Pouch Label



Box Label


Ref.: 7090121005V, Type: 4420804
Nasal Sampling Line with O2 Supply, Adult, Short Term



Pouch Label



Box Label

10 Nasal Sampling Line with O2 Supply, disposable, Paediatric, 2m/ft		LOT		22221.04_14	
REF	7090121006V				
cs	Odobrova nosní hadička s O ₂ , jednorázová, pediatrická	no	Nosální pravetakéslingslínie s O ₂ , til engjarnsbúrk, pediatrisk		
da	Nasaludtagingslinje med O ₂ , engangsbrug, barn	pl	Nosowa linia próbkująca z O ₂ , jednorazowa, dziecko		
de	Einweg-Probenableitung mit O ₂ , nasal, Pädiatrie	pt	Linha de amostragem nasal com O ₂ , descartável, pediátrico		
es	Nosa provintur O ₂ -ga, uhokorvelat kasutatav, laps	pt-BR	Linha de amostragem nasal com O ₂ , descartável, pediátrico		
et	Grasuti pihvits, õhukorraltootav, laste	ro	Linie de eşantionare nazală cu O ₂ , de unică folosință, pediatrică		
es	Linea de muestreo nasal con O ₂ descartable, pediatríco	ru	Насальная проботоотводящая линия с O ₂ , одноразовая, для детей		
fr	Tubulure d'échantillonnage nasale avec O ₂ , usage unique, enfant	sk	Názdňňa hadička na odber vzorky s O ₂ , jednorázová, pediatrická		
ko	산소 배관 검정 재질, 일회용, 소아용	sr	Vod za uzorkovanje za nazalnu primenu sa O ₂ , za jednokratnu upotrebu, za decu		
hr	Linija za uzorkovanje nosnog sekreta s O ₂ , jednokratna, pedijatrija	si	Nenānāryēlinja O ₂ nā karssa, kertakāyitōnā, lapset		
it	Linea di campionamento nasale con O ₂ , usa e getta, pazienti pediatrici	sv	Nasal provtagningssliding med O ₂ , för engångsbruk, barn		
iv	Deguna parakuq gēnasnas tīnja ar O ₂ , viereizlietojama, bērniem	tr	Burun Numune Alma Hattı, O ₂ li, tek kullanımlık, çocuk		
it	Nosinė mēginēji ēmimo linija s O ₂ , vienkartinē, vaikams	zh	带 O ₂ 的鼻黏膜采样管，一次性，儿童		
hu	Nasális mintavételi vezeték oxigénnel, egyszer használatos, gyermek	MD			
nl	Nasale bemonsteringslijn met O ₂ , voor eenmalig gebruik, pediatrisch	UDI			
ja	O ₂ 鼻サンプリングライン、使い捨て、小児用				
#	4420805	 ≥ 10 kg / 22 lbs		Short Term Use ≤ 24h	
					
					
					
Viamed Limited, 15 Station Road, Croyde Hills, Kington West Yorkshire, BD20 7DT, United Kingdom Email: sales@viamed.co.uk Tel: +44(0)1534 634542					
 blueprint medical GmbH & Co. KG Am der Trave 15 23923 Selmsdorf, Germany www.blueprintmedical.com					

Pouch Label

Box Label

Ref.: 7090121007V, Type: 4420806
Nasal Sampling Line with O2 Supply, Infant, Short Term

21 Nasal Sampling Line with O2 Supply, disposable, Infant, 2m/ftt

 REF 7090121007V	 LOT 22321.04.14
CS Odběrová nosní hadička s O ₂ , jednorázová, kojenecká	no Nasal pravetakingstangente med O ₂ , til engangsbruk, spedtarn
da Nasaludtagingslange med O ₂ , engangsbrug, spædbørn	pl Nosowa linia próbkująca z O ₂ , jednorazowa, niemowlę
de Einweg-Probenabnahmeleitung mit O ₂ , nasal, Kleinkinder	pt Linha de amostragem nasal com O ₂ , descartável, infante
es Nosoproviador O ₂ , úhoredokel kasutaval, väikepõlv	pe Línea de muestreo nasal con O ₂ , descartable, bebé
et Graafik pinnkõõ õhkuvõetavõttevõrk ju O ₂ , üksiktootmine, väikelaps	ro Linie de explorare nazală cu O ₂ , de unică folosință, bebeluși
es Linie de muestreo nasal con O ₂ , desechable, infante	ru Установочный проводотборный линии для O ₂ , одноразового, для младенцев
fr Tubulure d'échantillonnage nasale avec O ₂ , usage unique, nourisson	sk Násálna hadička na orálnu vzorku s O ₂ , jednorázová, pre dojčatá
KO 산소 배관 검체 채취관, 일회용, 소아용	sr Voda za uzorkovanje za nazalnu primenu sa O ₂ , za jednokratnu upotrebu, za novorođenčad
hr Linija za uzorkovanje nosnog sekreta s O ₂ , jednokratna, djeca	fi Nenäinlyetäjä O ₂ kanssa, kertakäyttöinen, pikkupotilas
it Linea di campionamento nasale con O ₂ , usa e getta, neonatale	sv Nasal provtagningssystem med O ₂ för engångsbruk, spädbarn
IV Deguna paruogu pēdasnas linija ar O ₂ , vireizotiekšmā, zīdaiņam	tr Burun Numune Alma Hattı, O ₂ li, tek kullanımlık, bebek
it Nasale mágicim éimio linia ar O ₂ , vlenkártná, kúdkákims	zh 市 O ₂ 的鼻黏膜采样管，一次性，婴幼儿
hu Nasális mintavétel vételező oxigénnel, egyszeres használatú, csecsemő	MD
nl Nasale bemonsteringslijn met O ₂ , voor eenmalig gebruik, zuigeling	UDI
ja O ₂ 鼻サンプリングライン、使い捨て、幼児用	 (01) 2 42516706 0117 7 (1) 210414 (17) 240414 (10) 22321.04.14

 # 4420806	 < 10 kg / 22 lbs	 Short Term Use ≤ 24h
--	---	---

 VIAMED Limited, 15 Station Road, Cross Halls, Kighley West Yorkshire, BD20 7PT, United Kingdom www.viamed.co.uk info@viamed.co.uk +44(0)1535 634542	 bluepoint medical GmbH & Co. KG Am der Trave 15 23923 Selmsdorf, Germany www.bluepointmedical.com
---	---

Pouch Label

Box Label

	PHA VersaStream Sampling Lines	Revision 3
Internal Document	Viamed Labeling	Seite 6 von 18

Ref.: 7090121008V, Type: 4420807
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121002V, Type: 4420801
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121009V, Type: 4420808
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121003V, Type: 4420802
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121010V, Type: 4420809
Nasal Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121004V, Type: 4420803
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121011V, Type: 4420810
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121005V, Type: 4420804
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h


Ref.: 7090121012V, Type: 4420811
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121006V, Type: 4420805
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h


Ref.: 7090121013V Type: 4420812
Nasal Sampling Line with O2 Supply, Infant, Long Term

Is identical with Ref.: 7090121007V, Type: 4420806
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121014V, Type: 4420813
Oral/Nasal Sampling Line, Adult, Short Term




Pouch Label



Box Label

Ref.: 7090121015V, Type: 4420814
Oral/Nasal Sampling Line, Paediatric, Short Term




Pouch Label



Box Label

	PHA VersaStream Sampling Lines	Revision 3
Internal Document	Viamed Labeling	Seite 8 von 18

Ref.: 7090121016V, Type: 4420815
Oral/Nasal Sampling Line with O2 Supply, Adult, Short Term




Pouch Label



Box Label

Ref.: 7090121017V, Type: 4420816
Oral/Nasal Sampling Line with O2 Supply, Paediatric, Short Term



Pouch Label



Box Label

	PHA VersaStream Sampling Lines	Revision 3
Internal Document	Viamed Labeling	Seite 9 von 18

Ref.: 7090121018V, Type: 4420817
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121014V, Type: 4420813
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121019V, Type: 4420818
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121015, Type: 4420814
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h


Ref.: 7090121020, Type: 4420819
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121016V, Type: 4420815
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h


Ref.: 7090121021V, Type: 4420820
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121017V, Type: 4420816
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121022V, Type: 4420821
Airway Adapter Sampling Line, Adult/Pediatric, Short Term




Pouch Label



Box Label

Ref.: 7090121023V, Type: 4420822
Airway Adapter Sampling Line, Infant, Short Term



Pouch Label



Box Label

	PHA VersaStream Sampling Lines	Revision 3
Internal Document	Viamed Labeling	Seite 11 von 18

Ref.: 7090121024, Type: 4420823
Airway Adapter Sampling Line, Adult/Pediatric, Long Term

Is identical with Ref.: 7090121022V, Type: 4420821
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121025V, Type: 4420824
Airway Adapter Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121023V, Type: 4420822
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

	PHA VersaStream Sampling Lines	Revision 3
Internal Document	Viamed Labeling	Seite 12 von 18

2. Label Layout: LLM Connector Sampling Lines

Ref.: 7090121033V, Type: 4420857
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121002V, Type: 4420801
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121034V, Type: 4420858
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121003V, Type: 4420802
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121035V, Type: 4420859
Nasal Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121004V, Type: 4420803
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121036V, Type: 4420860
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121005V, Type: 4420804
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121037V, Type: 4420861
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121006V, Type: 4420805
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121038V, Type: 4420862
Nasal Sampling Line with O2 Supply, Infant, Long Term

Is identical with Ref.: 7090121007V, Type: 4420806
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

	PHA VersaStream Sampling Lines	Revision 3
Internal Document	Viamed Labeling	Seite 13 von 18

Ref.: 7090121043V, Type: 4420867
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121014V, Type: 4420813
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121044V, Type: 4420868
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121015V, Type: 4420814
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121045V, Type: 4420869
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121016V, Type: 4420815
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121046V, Type: 4420870
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121017V, Type: 4420816
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121049V, Type: 4420873
Airway Adapter Sampling Line, Adult/Pediatric, Long Term

Is identical with Ref.: 7090121022V, Type: 4420821
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121050V, Type: 4420874
Airway Adapter Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121023V, Type: 4420822
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

	PHA VersaStream Sampling Lines	Revision 3
Internal Document	Viamed Labeling	Seite 14 von 18

3. Label Layout: ORI Connector Sampling Lines

Ref.: 7090121102V, Type: 4420900
Single Port Nasal Prong, Adult, Short Term

Is identical with Ref.: 7090121001V, Type: 4420800

Ref.: 7090121103V, Type: 4420901
Nasal Sampling Line, Adult, Short Term

Is identical with Ref.: 7090121002V, Type: 4420801

Ref.: 7090121104V, Type: 4420902
Nasal Sampling Line, Pediatric, Short Term

Is identical with Ref.: 7090121003V, Type: 4420802

Ref.: 7090121105V, Type: 4420903
Nasal Sampling Line, Infant, Short Term

Is identical with Ref.: 7090121004V, Type: 4420803

Ref.: 7090121106V, Type: 4420904
Nasal Sampling Line with O2 Supply, Adult, Short Term

Is identical with Ref.: 7090121005V, Type: 4420804

Ref.: 7090121107V, Type: 4420905
Nasal Sampling Line with O2 Supply, Pediatric, Short Term

Is identical with Ref.: 7090121006V, Type: 4420805

Ref.: 7090121108V, Type: 4420906
Nasal Sampling Line with O2 Supply, Infant, Short Term

Is identical with Ref.: 7090121007V, Type: 4420806

	PHA VersaStream Sampling Lines	Revision 3
Internal Document	Viamed Labeling	Seite 15 von 18

Ref.: 7090121109V, Type: 4420907
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121002V, Type: 4420807
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121110V, Type: 4420908
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121003V, Type: 4420808
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121111V, Type: 4420909
Nasal Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121004V, Type: 4420809
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121112V, Type: 4420910
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121005V, Type: 4420810
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121113V, Type: 4420911
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121006V, Type: 4420811
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121114V, Type: 4420912
Nasal Sampling Line with O2 Supply, Infant, Long Term

Is identical with Ref.: 7090121007V, Type: 4420812
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121115V, Type: 4420913
Oral/Nasal Sampling Line, Adult, Short Term

Is identical with Ref.: 7090121014V, Type: 4420813

	PHA VersaStream Sampling Lines	Revision 3
Internal Document	Viamed Labeling	Seite 16 von 18

Ref.: 7090121116V, Type: 4420914
Oral/Nasal Sampling Line, Pediatric, Short Term

Is identical with Ref.: 7090121015V, Type: 4420814

Ref.: 7090121117V, Type: 4420915
Oral/Nasal Sampling Line with O2 Supply, Adult, Short Term

Is identical with Ref.: 7090121016V, Type: 4420815

Ref.: 7090121118V, Type: 4420916
Oral/Nasal Sampling Line with O2 Supply, Pediatric, Short Term

Is identical with Ref.: 7090121017V, Type: 4420816

Ref.: 7090121119V, Type: 4420917
Nasal Sampling Line, Adult, Long Term

Is identical with Ref.: 7090121014V, Type: 4420817
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121120V, Type: 4420918
Nasal Sampling Line, Pediatric, Long Term

Is identical with Ref.: 7090121015V, Type: 4420818
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 70901210121V, Type: 4420919
Nasal Sampling Line with O2 Supply, Adult, Long Term

Is identical with Ref.: 7090121016V, Type: 4420819
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121122V, Type: 4420920
Nasal Sampling Line with O2 Supply, Pediatric, Long Term

Is identical with Ref.: 7090121017V, Type: 4420820
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

	PHA VersaStream Sampling Lines	Revision 3
Internal Document	Viamed Labeling	Seite 17 von 18

Ref.: 7090121123V, Type: 4420921
Airway Adapter Sampling Line, Adult/Pediatric, Short Term

Is identical with Ref.: 7090121022V, Type: 4420821

Ref.: 7090121124V, Type: 4420922
Airway Adapter Sampling Line, Infant, Short Term

Is identical with Ref.: 7090121023V, Type: 4420822

Ref.: 7090121125V, Type: 4420923
Airway Adapter Sampling Line, Adult/Pediatric, Long Term

Is identical with Ref.: 7090121022V, Type: 4420823
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121126V, Type: 4420924
Airway Adapter Sampling Line, Infant, Long Term

Is identical with Ref.: 7090121023V, Type: 4420824
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h

Ref.: 7090121140V, Type: 4420925
Airway Adapter Sampling Line, Adult/Pediatric, Long Term


Is identical with Ref.: 7090121022V, Type: 4420825
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h High Humidity

Ref.: 7090121141V, Type: 4420926
Airway Adapter Sampling Line, Infant, Long Term


Is identical with Ref.: 7090121023V, Type: 4420826
Only Change: Short Term Use $\leq$ 24h → Long Term Use $\leq$ 72h High Humidity

bluepoint MEDICAL	PHA VersaStream Sampling Lines	Revision 3
Internal Document	Viamed Labeling	Seite 18 von 18

Tag Layout:



Tag for short term sampling lines



Tag for long term sampling lines