



Catrin Hollings &lt;catrin.hollings@gmail.com&gt;

**[Test] VersaStream: Oridion Compatible Gas Sampling Lines**

1 message

**Viamed Ltd. - Sales** <oridion@viamed.co.uk>

Tue, Sep 13, 2022 at 9:24 AM


Reply-To: us12-9a7ee8662e-5f6d7567af@inbound.mailchimp.com

To: catrin.hollings@gmail.com




**VersaStream Gas Sampling Lines**  
**compatible for use with**  
**Capnographs and patient monitors with Oridion**  
**(Covidien, Nellcor, Medtronics, Philips, Masimo,**  
**Mindray, Zoll, Nihon Kohden...)**

Viamed are pleased to continually support Intensive Care Units with the supply of VersaStream gas sampling lines for high humidity applications. An innovative solution for sidestream capnography, offering a **one-piece sampling solution** **with increased performance**. These lines are compatible to the FilterLine<sup>®</sup> range of sampling lines and for use with Microstream<sup>®</sup> enabled capnograph and patient monitors.

Part Number: **4420925****VersaStream Adult/Paediatric Airway Adapter****High Humidity Sampling Line**

Long Term &lt; 72 hours Use

Box of 25

**£12.80 per patient**Part Number: **4420926****VersaStream Infant Airway Adapter****High Humidity Sampling Line**

Long Term &lt; 72 hours Use

Box of 25

**£13.60 per patient**

VersaStream gas sampling lines use the latest material technologies to allow the removal of moisture and water - even in cases of low temperatures and high humidity. They are very flexible and are suitable for both short-term and long-term applications.

## Why VersaStream gas sampling lines?


- One-piece sampling line solution
- **Cost-effective** alternative to Nafion<sup>®</sup> tubes - **from £6.80 per patient\***
- The gas-selective non-porous material facilitates fast and removal of condensation from the sampling line
- The integrated bacterial filter **prevents contamination**
- **Reliable measured gas values**, due to negligible dead space volume of the filter - even in cases of high respiration rates
- **Integrated bacteria filter and water removal** reduces the number of patient circuit parts and speeds up patient monitoring preparation time

For use with Oridion (Medtronic/Coviden/Nellcor) capnographs, as well as capnographs and patient monitors from other manufacturers that incorporate Microstream<sup>®</sup> technology, such as Philips, Dräger, Physio Control, Zoll...

These lines are available in a number of variants including: single nasal sampling line, nasal cannula, nasal/oral cannula and a one piece airway adapter sampling line.

They are compatible with a wide range of OEM equipment, please see our cross reference and selection guides for further information:

[www.viamed.co.uk/linesxref](http://www.viamed.co.uk/linesxref). Should you be unable to find the equipment you are using, please contact us so that we can search our extensive database.



Should you have any further questions regarding this, or would like to obtain samples for evaluation, please do not hesitate in contacting our sales team on **01535 634542** between 9:00am and 5:00pm Monday – Friday or email [oridion@viamed.co.uk](mailto:oridion@viamed.co.uk)

See our website for further information...

## Viamed's VersaStream range includes the following compatible gas sampling lines:

### Oridion Nelcor, Covidien, Medtronics

Compatible to the Filterline range of sampling lines and for use with Microstream enabled capnographs



### Luer Lock Male

For use with Capnographs and Multigas Monitors that incorporate a Luer Lock sampling port



### Viamed Masimo, Phase-In

For use with VM-2500 Sidestream Capnographs and Masimo/Phase-In technology enabled capnographs



Go to our Website for further information...

For information on compatability for the full range see our Cross Reference...

---

\* Type and size dependant

*Copyright © 2022 Viamed Ltd, All rights reserved.*

You are receiving this email because you or one of your colleagues has purchased a previous version of this product.

**Our mailing address is:**

Viamed Ltd

[15 Station Road](#)

[Cross Hills](#)

Keighley, West Yorkshire BD20 7DT

United Kingdom

[Add us to your address book](#)

Want to change how you receive these emails?

You can update your preferences or unsubscribe from this list.

