

Calibration Verification

The OxiQuant-S Oxygen Analyser has been tested for accuracy using medical grade Oxygen and air supplied by BOC. The output values were measured in accordance with procedure: VM3COP44.02.

Medical test gases used:

Medical Oxygen (BOC)	PL 0735/5000
Purity	99.5% (min)
Batch No.	
Medical Air (BOC)	PL 0735/5002R
Purity	21.0%
Batch No.	

Serial number of the UltraMaxO2
REF

Oxygen input:	UltraMaxO2 reading:
99.5%	%
Air input	UltraMaxO2 reading:
21.0%	%
Air Flow Test:	
Pressure Test:	

Tested by:

Signed: