



CID3988 Opera 00009147 Type:Company Supplier Customer

Teledyne Analytical Instruments 16830 Chestnut Street, , City Of Industry, California, 91748-1020, USA

Tel 0016269341589 [Log Call](#)

- [Account Home](#) [Main Details](#) [Account Memos](#) [Sales History](#) [Purchase Orders](#) [Supplier Returns](#) [Stock References](#)
- [Customer Returns](#) [Service Visits](#) [Distributor Agreements](#) [Documents Files](#) [Linked Issues](#) [Regulatory Agreements](#)

Teledyne Analytical Instruments - Returns Status Page

Summary Goods Waiting to Return to Teledyne Analytical Instruments

Stock Ref Description	Total To Return Supplier
`0110017 Oxygen sensor R-17MED	42
`0110051 Oxygen sensor T-1	3
`0110022 Oxygen sensor R-22MED	7
`0110024 Oxygen sensor R-24MED	2
`0110122 Oxygen sensor - Automotive R-22A	4
`0110117 Oxygen sensor - Automotive R-17A	2
`0110127 Oxygen sensor - Automotive R-17AV	3
`0110137 Oxygen sensor - Automotive R-17AH	15
`0110361 Oxygen sensor R-22AHJR Jikco OEM sensor	12
`0110121 Oxygen sensor - Automotive R-21A	13
`0110132 Oxygen sensor - Automotive R-22AV	10

Summary Goods Back at Teledyne Analytical Instruments

Stock Ref Description	Total With Supplier
0110081 Oxygen sensor UFO 130-2	2
0110024 Oxygen sensor R-24MED	66
0110017 Oxygen sensor R-17MED	16
0110022 Oxygen sensor R-22MED	26

0110051	Oxygen sensor T-1	1
0110057	Oxygen sensor T-7	1
0110023	Oxygen sensor R-23V	1
0110122	Oxygen sensor - Automotive R-22A	5
0110137	Oxygen sensor - Automotive R-17AH	3
0110361	Oxygen sensor R-22AHJR Jikco OEM sensor 2	

All Outstanding Items for Teledyne Analytical Instruments Supplier return File

Count	Return Box	RMA	Date Shipped to Supplier	Stock Reference	Description	Barcode	Serialnumber	Fault	Customer SRS
1	BOX778 Receipts		29 Dec 2021	0110081	Oxygen sensor UFO 130-2	 1578917	479035	Untested	 SRS67687
2	BOX778 Receipts		29 Dec 2021	0110081	Oxygen sensor UFO 130-2	 1573307	500812	Untested	 SRS67687
3	BOX781 Receipts		29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738169	609181	Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user	
4	BOX781 Receipts		29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738179	609191	Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user	
5	BOX781 Receipts		29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738175	609187	Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user	
6	BOX781 Receipts		29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738134	608460	Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user	
7	BOX781 Receipts		29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738150	608476	Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Output (10mV - 17mV) ERRATIC n Move Item to	

[illegible]

16	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738154 608480	Did not reach end user Output (10mV - 17mV) erratic n Move Item to location 126561 n Failed in Q.A. n Did not reach end user
17	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738149 608475	Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user
18	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738140 608466	Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user
19	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738474 609250	Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user
20	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738468 609244	Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user
21	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738481 609257	Pass Yes n Move Item to location 126561 n Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user
22	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738467 609243	Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user
23	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738473 609249	Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n

24	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738463 609239	Did not reach end user Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user
25	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738488 609264	Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user
26	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738462 609238	
27	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738469 609245	Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user
28	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738472 609248	Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user
29	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738487 609263	Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user
30	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738450 609226	Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user
31	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738458 609234	Output (10mV - 17mV) 1738458 n Pass Yes n Move Item to location 126561 n Output (10mV - 17mV)

						ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
32	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738484 609260	ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
33	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738506 609282	ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
34	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738433 609209	ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
35	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738511 609287	ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
36	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738440 609216	ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
37	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738489 609265	ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
38	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738452 609228	ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
39	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738445 609221	ERRATIC n

						Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV) erratic n
40	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738421 609197	Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
41	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738449 609225	ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
42	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738492 609268	ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
43	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738427 609203	ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
44	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738460 609236	ERRATIC n Move Item to location 126561 n Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
45	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738457 609233	ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
46	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738444 609220	ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user
47	BOX781	29 Dec	0110024	Oxygen	609261	Output (10mV -

	Reciepts	2021		sensor R-24MED	 1738485		17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
48	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738432	609208	17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
49	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738430	609206	17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
50	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738493	609269	17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
51	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738461	609237	17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
52	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738477	609253	17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
53	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738431	609207	17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
54	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R-24MED	 1738455	609231	17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)
55	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R-	 1738494	609270	17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV)

24MED

						ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Output (10mV - 17mV) ERRATIC n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV) erratic n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (10mV - 17mV) erratic n Move Item to location 126561 n Failed in Q.A. n Did not reach end user Output (7mV - 13mV) 13.6 Failed in Q.A. Did not reach end user Output (7mV - 13mV) erratic Failed in Q.A. Did not reach end user		
56	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738498 609274			
57	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738422 609198			
58	BOX781 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1738416 609192			
59	BOX801 Reciepts	29 Dec 2021	0110017	Oxygen sensor R- 17MED	 1743328 617363			
60	BOX801 Reciepts	29 Dec 2021	0110017	Oxygen sensor R- 17MED	 1766421 625571			
61	BOX801 Reciepts	29 Dec 2021	0110017	Oxygen sensor R- 17MED	 1685715 593247	Non-linear		SRS67865
62	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1604861 531997	Non-linear		SRS67851
63	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1750030 618744	Output (7mV - 13mV) 14.2 Failed in Q.A. Did not reach end user		
64	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1750105 618819	Output (7mV - 13mV) 13.5 Failed in Q.A. Did not reach end user		
65	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1750040 618754	Output (7mV - 13mV) 13.5 Failed in Q.A. Did not reach end user		

66	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1750315 609160	Output (7mV - 13mV) 13.8 Failed in Q.A. Did not reach end user
67	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1750082 618796	Output (7mV - 13mV) 13.5 Failed in Q.A. Did not reach end user
68	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1750161 618875	Output (7mV - 13mV) 13.7 n Move Item to location 126560 n Failed in Q.A. n Did not reach end user
69	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1750317 609162	Output (7mV - 13mV) 13.8 Failed in Q.A. Did not reach end user
70	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1750073 618787	Output (7mV - 13mV) 13.7 Failed in Q.A. Did not reach end user
71	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1750046 618760	Output (7mV - 13mV) 14.0 Failed in Q.A. Did not reach end user
72	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1750391 622729	Output (7mV - 13mV) 13.4 Failed in Q.A. Did not reach end user
73	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1750093 618807	Output (7mV - 13mV) 13.4 Failed in Q.A. Did not reach end user
74	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1750042 618756	Output (7mV - 13mV) 13.6 Failed in Q.A. Did not reach end user
75	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1750137 618851	Output (7mV - 13mV) 13.5 Failed in Q.A. Did not reach end user
76	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1750367 622705	Output (7mV - 13mV) 13.4 Failed in Q.A. Did not reach end user
77	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1750028 618742	Output (7mV - 13mV) 13.7 Failed in Q.A. Did not reach end user
78	BOX801 Reciepts	29 Dec 2021	0110022	Oxygen sensor R-	 1766320 625113	Output (7mV - 13mV) erratic

22MED

						Failed in Q.A. Did not reach end user
						Output (10mV - 17mV) erratic
79	BOX801 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1753021 608501	Failed in Q.A. Did not reach end user
						Output (10mV - 17mV) erratic
80	BOX801 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1753025 608505	Failed in Q.A. Did not reach end user
						Output (10mV - 17mV) erratic
81	BOX801 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1753017 608497	Failed in Q.A. Did not reach end user
						Output (10mV - 17mV) erratic
82	BOX801 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1753026 609292	Failed in Q.A. Did not reach end user
						Output (10mV - 17mV) erratic
83	BOX801 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1753018 608498	Failed in Q.A. Did not reach end user
						Output (10mV - 17mV) erratic
84	BOX801 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1753022 608502	Failed in Q.A. Did not reach end user
						Output (10mV - 17mV) erratic
85	BOX801 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1753014 608494	Failed in Q.A. Did not reach end user
						Output (10mV - 17mV) erratic
86	BOX801 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1753027 609293	Failed in Q.A. Did not reach end user
						Output (10mV - 17mV) erratic
87	BOX801 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1770420 628298	2.6 Erratic Failed in Q.A. Did not reach end user
						Output (10mV - 17mV) 2.5
88	BOX801 Reciepts	29 Dec 2021	0110024	Oxygen sensor R- 24MED	 1770416 628294	Erratic Failed in Q.A. Did not reach end user
						Leaking Failed in Q.A. Did not reach end user
89	BOX801 Reciepts	29 Dec 2021	0110051	Oxygen sensor T-1	 1642496 553531	
90	BOX801 Reciepts	29 Dec 2021	0110057	Oxygen sensor T-7	 1684862 590756	Calibration Failed
91	BOX812 Reciepts	29 Dec 2021	0110017	Oxygen sensor R- 17MED	 1823645 646457	Output (7mV - 13mV) 13.5 n Move Item to location 126559 n Failed in Q.A. n

 SRS67817

92	BOX812 Reciepts	29 Dec 2021	0110017	Oxygen sensor R- 17MED	 1823659 646471	Did not reach end user Output (7mV - 13mV) 00.1 n Move Item to location 126559 n Failed in Q.A. n Did not reach end user
93	BOX812 Reciepts	29 Dec 2021	0110017	Oxygen sensor R- 17MED	 1823649 646461	Output (7mV - 13mV) 13.2 sensor doesnt fit on tester very loose. n Move Item to location 126559 n Failed in Q.A. n Did not reach end user
94	BOX812 Reciepts	29 Dec 2021	0110017	Oxygen sensor R- 17MED	 1823626 646438	Output (7mV - 13mV) 13.4 n Move Item to location 126559 n Output (7mV - 13mV) 13.4 n Move Item to location 126559 n Failed in Q.A. n Did not reach end user
95	BOX812 Reciepts	29 Dec 2021	0110017	Oxygen sensor R- 17MED	 1823655 646467	Output (7mV - 13mV) 1823656 n Pass Yes n Move Item to location 126559 n Output (7mV - 13mV) 13.6 n Move Item to location 126559 n Failed in Q.A. n Did not reach end user
96	BOX812 Reciepts	29 Dec 2021	0110017	Oxygen sensor R- 17MED	 1823616 646428	Output (7mV - 13mV) 12.4 n Pass Yes n Move Item to location 126559 n Output (7mV - 13mV) 13.3 n Move Item to location 126559 n Failed in Q.A. n Did not reach end user
97	BOX812 Reciepts	29 Dec 2021	0110017	Oxygen sensor R- 17MED	 1823628 646440	Output (7mV - 13mV) 13.4 n Move Item to location 126559 n Failed in Q.A. n Did not reach end user
98	BOX812 Reciepts	29 Dec 2021	0110017	Oxygen sensor R- 17MED	 1598954 524654	Fluctuating O/p - During Use  SRS67918
99	BOX812 Reciepts	29 Dec 2021	0110017	Oxygen sensor R- 17MED	 1598955 524655	Fluctuating O/p - During Use  SRS67918

Output (7mV - 13mV) 12.9 n
Pass Yes n Move Item to location
126560 n Output (7mV - 13mV)
14.0 n Move Item to location
126560 n Failed in Q.A. n Did not reach end user
Output (7mV - 13mV) 12.7 n
Pass Yes n Move Item to location
126560 n Output (7mV - 13mV)
13.6 n Move Item to location
126560 n Failed in Q.A. n Did not reach end user
Output (7mV - 13mV) 9.6 n
Pass Yes n Move Item to location
126560 n Output (7mV - 13mV)
00.3 n Move Item to location
126560 n Failed in Q.A. n Did not reach end user
Output (7mV - 13mV) 11.0 n
Pass Yes n Move Item to location
126560 n Output (7mV - 13mV)
00.3 n Move Item to location
126560 n Failed in Q.A. n Did not reach end user
Output (7mV - 13mV) 10.5 n
Pass Yes n Move Item to location
126560 n Output (7mV - 13mV)
00.7 n Move Item to location
126560 n Failed

109	BOX812 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1678082 589419	in Q.A. n Did not reach end user Output (7mV - 13mV) 10.0 n Pass Yes n Move Item to location 126560 n Output (7mV - 13mV) 01.4 n Move Item to location 126560 n Failed in Q.A. n Did not reach end user Output (7mV - 13mV) 9.9 n Pass Yes n Move Item to location 126560 n Output (7mV - 13mV) 5.6 n Move Item to location 126560 n Output (7mV - 13mV) 5.5 n Move Item to location 126560 n Failed in Q.A. n Did not reach end user	 SRS67829
110	BOX812 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1678048 589385		 SRS67829
111	BOX812 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1695378 597205	Unstable	 SRS67923
112	BOX812 Reciepts	29 Dec 2021	0110022	Oxygen sensor R- 22MED	 1537243 473358	Unstable	 SRS67756
113	BOX812 Reciepts	29 Dec 2021	0110023	Oxygen sensor R- 23V	 1549281 V121059	Unstable	 SRS67756
114	BOX839 Reciepts	10 Feb 2022	0110122	Oxygen sensor - Automotive R-22A	 1849731 665124	Low Output 5.3mV Failed in Q.A. n Did not reach end user	
115	BOX839 Reciepts	10 Feb 2022	0110122	Oxygen sensor - Automotive R-22A	 1856048 668990	High Output 23.5 Failed in Q.A. n Did not reach end user	
116	BOX839 Reciepts	10 Feb 2022	0110122	Oxygen sensor - Automotive R-22A	 1849692 665085	No Output 00.1mV Failed in Q.A. Did not reach end user	
117	BOX839 Reciepts	10 Feb 2022	0110122	Oxygen sensor - Automotive R-22A	 1849599 664992	Low Output 2.6mV Failed in Q.A. Did not reach end user	
118	BOX839 Reciepts	10 Feb 2022	0110122	Oxygen sensor - Automotive R-22A	 1905108 693147	Low Output 6.4mV Failed in Q.A. Did not reach end user	
119	BOX839 Reciepts	10 Feb 2022	0110137	Oxygen sensor - Automotive R-17AH	 1910687 697102	High Output 14.5mV Failed in Q.A. Did not reach end user	
120	BOX839 Reciepts	10 Feb 2022	0110137	Oxygen sensor -	 1910677 697092	No Output 00.1mV Failed in	

121	BOX839 Reciepts	10 Feb 2022	0110137	Automotive R-17AH Oxygen sensor - Automotive R-17AH	 1910817 697227	Q.A. Did not reach end user Low Output 10.3mV Failed in Q.A. Did not reach end user
122	BOX839 Reciepts	10 Feb 2022	0110361	Oxygen sensor R- 22AHJR Jikco OEM sensor	 1906118 A695163	Low Output 9.8mV failed in Q/A did not reach end user
123	BOX839 Reciepts	10 Feb 2022	0110361	Oxygen sensor R- 22AHJR Jikco OEM sensor	 1906122 A695167	Low Output 8.6 mV Failed in Q.A. Did not reach end user
124	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1898068 688393	No Output 00.1mV erratic Failed in Q.A. Did not reach end user
125	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1898102 688427	No Output 00.1mV Failed in Q.A. Did not reach end user
126	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1902900 688483	High Output 13.3mV Failed in Q.A. Did not reach end user
127	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1889073 673681	No Output 00.1mV Failed in Q.A. Did not reach end user
128	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1897648 687973	No Output 00.1mV erratic Failed in Q.A. Did not reach end user
129	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1897675 688000	No Output 00.1mv Failed in Q.A. Did not reach end user
130	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1889168 673776	No Output 00.1mV Failed in Q.A. Did not reach end user
131	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1897802 688127	No Output 00.1mV Failed in Q.A. n Did not reach end user
132	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1897889 688214	High Output13.3mV Failed in Q.A. n Did not reach end user
133	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1897693 688018	No Output 00.1mV Failed in Q.A. Did not reach end user
134	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1897959 688284	No Output 00.1mV Failed in Q.A. Did not reach end user

135	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1897878 688203	No Output 00.1mV Failed in Q.A. Did not reach end user
136	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1889173 673781	No Output 00.1mV erratic Failed in Q.A. Did not reach end user
137	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1889157 673765	No Output 00.1mV Failed in Q.A. Did not reach end user
138	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1889385 673885	No Output 00.1mV Failed in Q.A Did not reach end user
139	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1889169 673777	No Output 00.1mVFailed in Q.A. Did not reach end user
140	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1889159 673767	No Output 00.1mVFailed in Q.A. Did not reach end user
141	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1889388 673888	No Output 00.1mVFailed in Q.A. Did not reach end user
142	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1889708 674208	No Output 00.1mV erratic Failed in Q.A. Did not reach end user
143	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1889306 673806	No Output 00.1mV Failed in Q.A. Did not reach end user
144	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1903214 688797	No Output 00.1 mV Failed in Q.A. Did not reach end user
145	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1910930 696958	Output erratic 0.1 -9.3 Failed in Q.A. Did not reach end user
146	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1910963 696991	Low Output 3.0mV Failed in Q.A. Did not reach end user
147	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1903239 688822	High Output 13.6 Failed in Q.A. Did not reach end user
148	BOX840 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1903199 688782	No Output 00.1mV Failed in Q.A. Did not reach end user
149	BOX840 Reciepts	Not Yet Shipped	0110051	Oxygen sensor T-1	 1849489 664734	High Output Did not reach end user Failed in QA
150	BOX840 Reciepts	Not Yet Shipped	0110051	Oxygen sensor T-1	 1849474 664719	No Output Failed in Q.A. Did not

							reach end user
151	BOX840 Reciepts	Not Yet Shipped	0110051	Oxygen sensor T-1	 1849395 664640		High Output Failed in QA Did not reach end user
152	BOX861 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1723755 606017	Non Linear	 SRS67996
153	BOX861 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1743254 617289	Non Linear	 SRS67971
154	BOX861 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1685751 593283	Non Linear	 SRS68040
155	BOX861 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1685753 593285	Non Linear	 SRS68040
156	BOX861 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1823523 646335	Manufacturing Error	 SRS68052
157	BOX861 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1823524 646336	Unstable	 SRS68052
158	BOX861 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1743181 617216	No Output	 SRS67999
159	BOX861 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1723668 605930	Non Linear	 SRS67960
160	BOX861 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1743173 617208	High Output	 SRS67960
161	BOX861 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1652350 566049	High Output	 SRS67984
162	BOX861 Reciepts	Not Yet Shipped	0110022	Oxygen sensor R- 22MED	 1695377 597204	Unstable	 SRS67954
163	BOX861 Reciepts	Not Yet Shipped	0110024	Oxygen sensor R- 24MED	 1753015 608495	No Output	 SRS68035
164	BOX861 Reciepts	Not Yet Shipped	0110024	Oxygen sensor R- 24MED	 1738435 609211	Non Linear	 SRS67936
165	BOX862 Reciepts	Not Yet Shipped	0110122	Oxygen sensor - Automotive R-22A	 1664051 573056	High Output	 SRS67945
166	BOX862 Reciepts	Not Yet Shipped	0110122	Oxygen sensor - Automotive R-22A	 1664213 573449	Unstable	 SRS67944
167	BOX863 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1959402 782035	Output (7mV - 13mV) 00.1 Failed in Q.A. n Did not reach end user	
168	BOX863 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1953107 774192	Output (7mV - 13mV) 00.1 Failed in Q.A. n Did not reach end user	

169	BOX863 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R-17MED	 1952884 773969	Output (7mV - 13mV) 00.1 Failed in Q.A. n Did not reach end user
170	BOX863 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R-17MED	 1952745 773830	Output (7mV - 13mV) 00.1 Failed in Q.A. n Did not reach end user
171	BOX863 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R-17MED	 1959354 781987	Output (7mV - 13mV) 5.8 Failed in Q.A. n Did not reach end user
172	BOX863 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R-17MED	 1914245 701320	Output (7mV - 13mV) 00.1 Failed in Q.A. n Did not reach end user
173	BOX863 Reciepts	Not Yet Shipped	0110022	Oxygen sensor R-22MED	 1960784 786565	Output (7mV - 13mV) 5.0 Failed in Q.A. n Did not reach end user
174	BOX863 Reciepts	Not Yet Shipped	0110022	Oxygen sensor R-22MED	 1960825 786606	Output (7mV - 13mV) 1.6 Failed in Q.A. n Did not reach end user
175	BOX863 Reciepts	Not Yet Shipped	0110022	Oxygen sensor R-22MED	 1604460 531596	Output (7mV - 13mV) 0.9 Failed in Q.A. n Did not reach end user
176	BOX864 Reciepts	Not Yet Shipped	0110117	Oxygen sensor - Automotive R-17A	 1775083 601125	Output (7mV - 13mV) 13.6 Failed in Q.A. n Did not reach end user
177	BOX864 Reciepts	Not Yet Shipped	0110117	Oxygen sensor - Automotive R-17A	 1721730 601054	Output (7mV - 13mV) 15.4 Failed in Q.A. n Did not reach end user
178	BOX864 Reciepts	Not Yet Shipped	0110127	Oxygen sensor - Automotive R-17AV	 1951898 778111	Output (9mV - 13mV) 13.8 Failed in Q.A. n Did not reach end user
179	BOX864 Reciepts	Not Yet Shipped	0110137	Oxygen sensor - Automotive R-17AH	 1910808 697218	Output (10.5mV - 13.5mV) 10.4 Failed in Q.A. n Did not reach end user
180	BOX864 Reciepts	Not Yet Shipped	0110361	Oxygen sensor R-22AHJR Jikco OEM sensor	 1944813 A771170	Output (10.5 - 13.5 mV) 10.4 Failed in Q.A. n Did not reach end user
181	BOX864 Reciepts	Not Yet Shipped	0110361	Oxygen sensor R-22AHJR	 1953574 A771192	Output (10.5 - 13.5 mV) 10.4 Failed in Q.A. n

182	BOX873 Reciepts	Not Yet Shipped	0110121	Jikco OEM sensor Oxygen sensor - Automotive R-21A	 2005686	807384	Did not reach end user High Output 15.6 mv Failed in Q.A. n Did not reach end user
183	BOX873 Reciepts	Not Yet Shipped	0110121	Oxygen sensor - Automotive R-21A	 2005782	807413	High Output 15.4mv Failed in Q.A. n Did not reach end user
184	BOX873 Reciepts	Not Yet Shipped	0110121	Oxygen sensor - Automotive R-21A	 2005777	807495	High Output 15.1 mv Failed in Q.A. n Did not reach end user
185	BOX873 Reciepts	Not Yet Shipped	0110121	Oxygen sensor - Automotive R-21A	 2005725	807443	00.5-11.1mv erratic Failed in Q.A. n Did not reach end user
186	BOX873 Reciepts	Not Yet Shipped	0110122	Oxygen sensor - Automotive R-22A	 2008347	815709	fail leaking Failed in Q.A. n Did not reach end user
187	BOX873 Reciepts	Not Yet Shipped	0110122	Oxygen sensor - Automotive R-22A	 2003907	808230	Low Output 6.3mV Failed in Q.A. n Did not reach end user
188	BOX873 Reciepts	Not Yet Shipped	0110127	Oxygen sensor - Automotive R-17AV	 2003574	807323	00.1mV Failed in Q.A. Did not reach end user
189	BOX873 Reciepts	Not Yet Shipped	0110127	Oxygen sensor - Automotive R-17AV	 2003784	807367	High Output 13.5mV Failed in Q.A. Did not reach end user
190	BOX873 Reciepts	Not Yet Shipped	0110132	Oxygen sensor - Automotive R-22AV	 2004902	810431	Low Output 01.5mV Failed in Q.A. Did not reach end user
191	BOX873 Reciepts	Not Yet Shipped	0110132	Oxygen sensor - Automotive R-22AV	 2004920	810449	High Output 14.1 mV Failed in Q.A. Did not reach end user
192	BOX873 Reciepts	Not Yet Shipped	0110132	Oxygen sensor - Automotive R-22AV	 2004921	810450	High Output 13.6mV Failed in Q.A. Did not reach end user
193	BOX873 Reciepts	Not Yet Shipped	0110132	Oxygen sensor - Automotive R-22AV	 2008101	817252	High Output 13.6mV Failed in Q.A. Did not reach end user
194	BOX873 Reciepts	Not Yet Shipped	0110132	Oxygen sensor - Automotive R-22AV	 2004933	810462	High Output 13.7 Failed in Q.A. Did not reach end user
195	BOX873 Reciepts	Not Yet Shipped	0110132	Oxygen sensor - Automotive R-22AV	 2008039	817190	High Output 13.5mV Failed in Q.A. Did not reach end user
196	BOX873 Reciepts	Not Yet Shipped	0110132	Oxygen sensor - Automotive R-22AV	 1994253	808627	High Output 13.5mV Failed in Q.A. Did not reach end user
197	BOX873	Not Yet	0110132	Oxygen		802967	Low Output (

	Reciepts	Shipped		sensor - Automotive R-22AV		1976835		9mV - 13mV) 8.6mV Failed in Q.A. Did not reach end user
198	BOX873 Reciepts	Not Yet Shipped	0110132	Oxygen sensor - Automotive R-22AV		1976726	802858	Low output 8.8mV Failed in QA did not reach end user
199	BOX873 Reciepts	Not Yet Shipped	0110132	Oxygen sensor - Automotive R-22AV		1976707	802839	Low Output 8.7mV Failed in Q.A. Did not reach end user
200	BOX873 Reciepts	Not Yet Shipped	0110137	Oxygen sensor - Automotive R-17AH		2003809	807242	High Output (10.5mV - 13.5mV) 13.7mV Failed in Q.A. Did not reach end user
201	BOX873 Reciepts	Not Yet Shipped	0110137	Oxygen sensor - Automotive R-17AH		2003832	807265	High Output (10.5mV - 13.5mV) 14.7mV Failed in Q.A. Did not reach end user
202	BOX873 Reciepts	Not Yet Shipped	0110137	Oxygen sensor - Automotive R-17AH		2003830	807263	High Output (10.5mV - 13.5mV) 14.2mV Failed in Q.A. Did not reach end user
203	BOX873 Reciepts	Not Yet Shipped	0110137	Oxygen sensor - Automotive R-17AH		2003520	807221	High Output (10.5mV - 13.5mV) 13.6mV Failed in Q.A. Did not reach end user
204	BOX873 Reciepts	Not Yet Shipped	0110137	Oxygen sensor - Automotive R-17AH		2003838	807271	High Output (10.5mV - 13.5mV) 13..9mV Failed in Q.A. Did not reach end user
205	BOX873 Reciepts	Not Yet Shipped	0110137	Oxygen sensor - Automotive R-17AH		2003807	807240	High Output (10.5mV - 13.5mV) 14.1mV Failed in Q.A. Did not reach end user
206	BOX873 Reciepts	Not Yet Shipped	0110137	Oxygen sensor - Automotive R-17AH		2003810	807243	High Output (10.5mV - 13.5mV) 14.9mV Failed in Q.A. Did not reach end user
207	BOX873 Reciepts	Not Yet Shipped	0110137	Oxygen sensor - Automotive R-17AH		2003445	807146	High Output 14.7 mV Failed in Q.A. Did not reach end user
208	BOX873 Reciepts	Not Yet Shipped	0110137	Oxygen sensor - Automotive R-17AH		2003492	807193	High Output (10.5mV - 13.5mV) 14.8mV Failed in Q.A. Did not reach end user

209	BOX873 Reciepts	Not Yet Shipped	0110137	Oxygen sensor - Automotive R-17AH	 2003490	807191	High Output (10.5mV - 13.5mV) 14.0mV Failed in Q.A.Did not reach end user
210	BOX873 Reciepts	Not Yet Shipped	0110137	Oxygen sensor - Automotive R-17AH	 2003457	807158	High Output 14.5mV Move Failed in Q.A. Did not reach end user
211	BOX873 Reciepts	Not Yet Shipped	0110137	Oxygen sensor - Automotive R-17AH	 2003448	807149	00.1mV Q.A. Did not reach end user
212	BOX873 Reciepts	Not Yet Shipped	0110137	Oxygen sensor - Automotive R-17AH	 2003484	807185	High Output (10.5mV - 13.5mV) 14.3mV Failed in Q.A. Did not reach end user
213	BOX873 Reciepts	Not Yet Shipped	0110137	Oxygen sensor - Automotive R-17AH	 2003453	807154	High Output 14.0 Failed in Q.A. Did not reach end user
214	BOX873 Reciepts	Not Yet Shipped	0110361	Oxygen sensor R- 22AHJR Jikco OEM sensor	 2003424	A809903	High Output (10.5 - 13.5 mV) 13.8mV Failed in Q.A. Did not reach end user
215	BOX873 Reciepts	Not Yet Shipped	0110361	Oxygen sensor R- 22AHJR Jikco OEM sensor	 2003293	A809924	High Output (10.5 - 13.5 mV) 13.8mV Failed in Q.A. Did not reach end user
216	BOX873 Reciepts	Not Yet Shipped	0110361	Oxygen sensor R- 22AHJR Jikco OEM sensor	 2003258	A809799	High Output (10.5 - 13.5 mV) 14.0mV Failed in Q.A. Did not reach end user
217	BOX873 Reciepts	Not Yet Shipped	0110361	Oxygen sensor R- 22AHJR Jikco OEM sensor	 2003220	A809756	High Output (10.5 - 13.5 mV) 13.9mV Failed in Q.A. Did not reach end user
218	BOX873 Reciepts	Not Yet Shipped	0110361	Oxygen sensor R- 22AHJR Jikco OEM sensor	 2003235	A809771	High Output (10.5 - 13.5 mV) 14.9mV Failed in Q.A. Did not reach end user
219	BOX873 Reciepts	Not Yet Shipped	0110361	Oxygen sensor R- 22AHJR Jikco OEM sensor	 2003400	A809879	High Output (10.5 - 13.5 mV) 14.3 Failed in Q.A. Did not reach end user
220	BOX873 Reciepts	Not Yet Shipped	0110361	Oxygen sensor R- 22AHJR Jikco OEM sensor	 2003409	A809888	Low Output (10.5 - 13.5 mV) 10.4mV Failed in Q.A. Did not reach end user
221	BOX873 Reciepts	Not Yet Shipped	0110361	Oxygen sensor R- 22AHJR	 2003352	A809831	Low Output (10.5 - 13.5 mV) 10.4mV Failed in

222	BOX873 Reciepts	Not Yet Shipped	0110361	Jikco OEM sensor Oxygen sensor R- 22AHJR Jikco OEM sensor	 2003433 A809912	Q.A. Did not reach end user Low Output (10.5 - 13.5 mV) 10.1mV Failed in Q.A. Did not reach end user	
223	BOX873 Reciepts	Not Yet Shipped	0110361	Oxygen sensor R- 22AHJR Jikco OEM sensor	 2003415 A809894	High Output (10.5 - 13.5 mV) 13.9mV Failed in Q.A. Did not reach end user	
224	BOX874 Reciepts	Not Yet Shipped	0110017	Oxygen sensor R- 17MED	 1743285 617320	Unstable	 SRS68145
225	BOX874 Reciepts	Not Yet Shipped	0110022	Oxygen sensor R- 22MED	 1750420 622758	Non Linear	 SRS68143
226	BOX874 Reciepts	Not Yet Shipped	0110022	Oxygen sensor R- 22MED	 1750421 622759	Unstable	 SRS68143
227	BOX874 Reciepts	Not Yet Shipped	0110022	Oxygen sensor R- 22MED	 1764008 624826	Unstable	 SRS68143
228	BOX874 Reciepts	Not Yet Shipped	0110121	Oxygen sensor - Automotive R-21A	 1736109 608267	High Output	 SRS68108
229	BOX874 Reciepts	Not Yet Shipped	0110121	Oxygen sensor - Automotive R-21A	 1736114 608272	High Output	 SRS68108
230	BOX874 Reciepts	Not Yet Shipped	0110121	Oxygen sensor - Automotive R-21A	 1736111 608269	No Output	 SRS68108
231	BOX874 Reciepts	Not Yet Shipped	0110121	Oxygen sensor - Automotive R-21A	 1736104 608262	High Output	 SRS68108
232	BOX874 Reciepts	Not Yet Shipped	0110121	Oxygen sensor - Automotive R-21A	 1736123 608281	Unstable	 SRS68108
233	BOX874 Reciepts	Not Yet Shipped	0110121	Oxygen sensor - Automotive R-21A	 1736203 608361	Non Linear	 SRS68108
234	BOX874 Reciepts	Not Yet Shipped	0110121	Oxygen sensor - Automotive R-21A	 1736213 608371	Unstable	 SRS68108
235	BOX874 Reciepts	Not Yet Shipped	0110121	Oxygen sensor - Automotive R-21A	 1736215 608373	Non Linear	 SRS68108
236	BOX874 Reciepts	Not Yet Shipped	0110121	Oxygen sensor - Automotive R-21A	 1856124 669463	Non Linear	 SRS68108

Current Customer SRS Stock References in Building

Stock Ref	Description	Customer Stock In House	Customer Stock At Supplier
0111235	Oxygen monitor with alarms - MX300-I.	1	
0110017	Oxygen sensor R-17MED	3	
0110057	Oxygen sensor T-7	2	
0110237	Oxygen sensor R-17S1 (high output)	1	
0110122	Oxygen sensor - Automotive R-22A	3	
0111230	Oxygen analyser w/o alarms - AX300-I.	2	

ALL Outstanding SRS/SRNS open

Stock Ref	Description	Barcode	SRS	SRN
0111235	Oxygen monitor with alarms - MX300-I.	1418471	 SRS68156	SRN35037
0110017	Oxygen sensor R-17MED	1450310	 SRS68156	SRN35039
0110057	Oxygen sensor T-7	1684907	 SRS68177	SRN35076
0110057	Oxygen sensor T-7	1684908	 SRS68177	SRN35077
0110237	Oxygen sensor R-17S1 (high output)	2025331	 SRS68230	SRN35122
0110122	Oxygen sensor - Automotive R-22A	1949159	 SRS68208	SRN35136
0110122	Oxygen sensor - Automotive R-22A	1849532	 SRS68208	SRN35135
0110122	Oxygen sensor - Automotive R-22A	1948248	 SRS68208	SRN35137
0111230	Oxygen analyser w/o alarms - AX300-I.	1234711	 SRS68219	SRN35138
0111230	Oxygen analyser w/o alarms - AX300-I.	1234713	 SRS68219	SRN35139
0110017	Oxygen sensor R-17MED	1564300	 SRS68219	SRN35142
0110017	Oxygen sensor R-17MED		 SRS68219	SRN35143