



# NeoPEEP Resuscitation Circuit

Lorem ipsum dolor sit amet, consectetur adipiscing elit. Phasellus ligula???



**For use  
in Infant  
Resuscitation  
Cabinets**

The NeoPEEP is a single patient use neonatal resuscitation circuit incorporating a variable PEEP valve.

There are three steps involved in selecting your resuscitation circuit items:

- Step One:** Select the quantity of resuscitation circuits required.
- Step Two:** Select the relevant adapter for your infant resuscitator, if required.
- Step Three:** Select the patient mask size you require.



# NeoPEEP Resuscitation Circuit

## Step One: Resuscitation Circuit

| Part Number | Model Ref.  | Description                                                                      | Qty         |
|-------------|-------------|----------------------------------------------------------------------------------|-------------|
| 3210011     | 800100/BLUE | <b>NeoPEEP Neonatal Resuscitation Circuit</b><br>Single Patient Use, 1.3m length | Single Unit |
| 3210013     | 800100/BLUE | <b>NeoPEEP Neonatal Resuscitation Circuit</b><br>Single Patient Use, 1.3m length | Box of 20   |



## Step Two: Resuscitation Circuit Connectors

The NeoPEEP will connect to any resuscitator with a 15 mm I.D. patient outlet.

The NeoPEEP will also connect directly to standard endotracheal tubing.

For connecting the NeoPEEP to a resuscitator without a standard 15 mm I.D. patient outlet various connectors are available:

| Part Number | Model Ref. | Description                                                                 | For Connecting to:      | Qty         |
|-------------|------------|-----------------------------------------------------------------------------|-------------------------|-------------|
| 0120139     | PF 1510    | <b>Adapter: 15 mm I.D. - 10 mm I.D.</b><br>Silicone, Reusable               | Fisher & Paykel NeoPuff | Single Unit |
| 0120140     | PF 1515    | <b>Adapter: 15 mm I.D. - 15 mm I.D.</b><br>Silicone, Reusable               | Viamed Tom Thumb        | Single Unit |
| 0120141     | KC 2124    | <b>Adapter: 15 mm I.D./22 mm O.D. - 15 mm I.D./22 mm O.D.</b><br>Disposable | Viamed Tom Thumb        | Single Unit |
| 0120151     | KC 2124    | <b>Adapter: 15 mm I.D./22 mm O.D. - 15 mm I.D./22 mm O.D.</b><br>Disposable | Viamed Tom Thumb        | Pack of 20  |



## Step Three: Patient Masks

| Part Number | Model Ref. | Description                                         | Qty         |
|-------------|------------|-----------------------------------------------------|-------------|
| 3210070     | 200600D    | <b>Single use silicone round facemask - Size 00</b> | Single Unit |
| 3210073     | 200600D    | <b>Single use silicone round facemask - Size 00</b> | Pack of 20  |
| 3210076     | 200600D    | <b>Single use silicone round facemask - Size 00</b> | Box of 100  |
| 3210071     | 200500D    | <b>Single use silicone round facemask - Size 01</b> | Single Unit |
| 3210074     | 200500D    | <b>Single use silicone round facemask - Size 01</b> | Pack of 20  |
| 3210077     | 200500D    | <b>Single use silicone round facemask - Size 01</b> | Box of 100  |
| 3210072     | 200400D    | <b>Single use silicone round facemask - Size 02</b> | Single Unit |
| 3210075     | 200400D    | <b>Single use silicone round facemask - Size 02</b> | Pack of 20  |
| 3210078     | 200400D    | <b>Single use silicone round facemask - Size 02</b> | Box of 100  |



Distributed by:



**Viamed Limited** · 15 Station Road · Cross Hills  
Keighley · West Yorkshire · BD20 7DT · United Kingdom  
Tel: +44 (0)1535 634 542 Fax: +44 (0)1535 635 582  
Email: [info@viamed.co.uk](mailto:info@viamed.co.uk) Website: [www.viamed.co.uk](http://www.viamed.co.uk)



Part Number: 3210011  
Date: 03/11



BS EN ISO 9001:2008  
ISO 13485:2003