

Viamed Returns Shipping Invoice

Delivery Address:
Teledyne Analytical Instruments
Luz Acosta
16830 Chestnut Street

Date: 30 Sep 2022
RMA:
Viamed RN: RTS1100801

Shipper Address:
Viamed
15 Station Road
Cross Hills
Keighley
West Yorkshire
BD20 7DT
+44 (0) 1535 634542
+44 (0) 1535 635582
info@viamed.co.uk

91748-1020
0016269341589

Returns Under Warranty Value for Customs : \$100.00
See Attached Listing for Serial numbers and Faults

Part Number	P/N / Value	Description	Quantity
C43690-R17MED	0110017 (\$22.00 Ea.) Tariff Code: 9019209000 Country of Origin: USA	Teledyne Oxygen sensor R-17MED	3
C44611-R22MED	0110022 (\$25.00 Ea.) Tariff Code: 9019200000 Country of Origin: USA	Teledyne Oxygen sensor R-22MED	17
C61340-R24MED	0110024 (\$25.00) Tariff Code: 9019200000 Country of Origin: USA	Teledyne Oxygen sensor R-24MED	10
A37016-T1	0110051 (\$27.00) Tariff Code: 9019200000 Country of Origin: USA	Teledyne Oxygen sensor T-1	1
A51327-T7	0110057 (\$46.69 Ea.) Tariff Code: 9019200000 Country of Origin: USA	Teledyne Oxygen sensor T-7	1

Part Number	Serial Number	RMA	Fault Description
C43690-R17MED	593247		Non-linear
C43690-R17MED	617363		Output (7mV - 13mV) 13.6 Failed in Q.A. Did not reach end user
C43690-R17MED	625571		Output (7mV - 13mV) erratic Failed in Q.A. Did not reach end user
C44611-R22MED	531997		Non-linear
C44611-R22MED	609160		Output (7mV - 13mV) 13.8 Failed in Q.A. Did not reach end user
C44611-R22MED	609162		Output (7mV - 13mV) 13.8 Failed in Q.A. Did not reach end user
C44611-R22MED	618742		Output (7mV - 13mV) 13.7 Failed in Q.A. Did not reach end user
C44611-R22MED	618744		Output (7mV - 13mV) 14.2 Failed in Q.A. Did not reach end user
C44611-R22MED	618754		Output (7mV - 13mV) 13.5 Failed in Q.A. Did not reach end user
C44611-R22MED	618756		Output (7mV - 13mV) 13.6 Failed in Q.A. Did not reach end user
C44611-R22MED	618760		Output (7mV - 13mV) 14.0 Failed in Q.A. Did not reach end user
C44611-R22MED	618787		Output (7mV - 13mV) 13.7 Failed in Q.A. Did not reach end user
C44611-R22MED	618796		Output (7mV - 13mV) 13.5 Failed in Q.A. Did not reach end user
C44611-R22MED	618807		Output (7mV - 13mV) 13.4 Failed in Q.A. Did not reach end user
C44611-R22MED	618819		Output (7mV - 13mV) 13.5 Failed in Q.A. Did not reach end user
C44611-R22MED	618851		Output (7mV - 13mV) 13.5 Failed in Q.A. Did not reach end user
C44611-R22MED	618875		Output (7mV - 13mV) 13.7 n Move Item to location 126560 n Failed in Q.A. n Did not reach end user
C44611-R22MED	622705		Output (7mV - 13mV) 13.4 Failed in Q.A. Did not reach end user
C44611-R22MED	622729		Output (7mV - 13mV) 13.4 Failed in Q.A. Did not reach end user
C44611-R22MED	625113		Output (7mV - 13mV) erratic Failed in Q.A. Did not reach end user
C61340-R24MED	608494		Output (10mV - 17mV) erratic Failed in Q.A. Did not reach end user
C61340-R24MED	608497		Output (10mV - 17mV) erratic Failed in Q.A. Did not reach end user
C61340-R24MED	608498		Output (10mV - 17mV) erratic Failed in Q.A. Did not reach end user
C61340-R24MED	608501		Output (10mV - 17mV) erratic Failed in Q.A. Did not reach end user
C61340-R24MED	608502		Output (10mV - 17mV) erratic Failed in Q.A. Did not reach end user
C61340-R24MED	608505		Output (10mV - 17mV) erratic Failed in Q.A. Did not reach end user
C61340-R24MED	609292		Output (10mV - 17mV) erratic Failed in Q.A. Did not reach end user
C61340-R24MED	609293		Output (10mV - 17mV) erratic Failed in Q.A. Did not reach end user
C61340-R24MED	628294		Output (10mV - 17mV) 2.5 Erratic Failed in Q.A. Did not reach end user

Part Number	Serial Number	RMA	Fault Description
C61340-R24MED	628298		Output (10mV - 17mV) 2.6 Erratic Failed in Q.A. Did not reach end user
A37016-T1	553531		Leaking Failed in Q.A. Did not reach end user
A51327-T7	590756		Calibration Failed