


VIAMED

Suppliers of Electromedical Equipment

2023



Hawthorn Garden, Ribblesdale

January

M	30	2	9	16	23
T	31	3	10	17	24
W		4	11	18	25
T		5	12	19	26
F		6	13	20	27
S		7	14	21	28
S	1	8	15	22	29

February

M		6	13	20	27
T		7	14	21	28
W	1	8	15	22	
T	2	9	16	23	
F	3	10	17	24	
S	4	11	18	25	
S	5	12	19	26	

March

M		6	13	20	27
T		7	14	21	28
W	1	8	15	22	29
T	2	9	16	23	30
F	3	10	17	24	31
S	4	11	18	25	
S	5	12	19	26	

April

M		3	10	17	24
T		4	11	18	25
W		5	12	19	26
T		6	13	20	27
F		7	14	21	28
S	1	8	15	22	29
S	2	9	16	23	30

May

M	1	8	15	22	29
T	2	9	16	23	30
W	3	10	17	24	31
T	4	11	18	25	
F	5	12	19	26	
S	6	13	20	27	
S	7	14	21	28	

June

M		5	12	19	26
T		6	13	20	27
W		7	14	21	28
T	1	8	15	22	29
F	2	9	16	23	30
S	3	10	17	24	
S	4	11	18	25	

July

M	31	3	10	17	24
T		4	11	18	25
W		5	12	19	26
T		6	13	20	27
F		7	14	21	28
S	1	8	15	22	29
S	2	9	16	23	30

August

M		7	14	21	28
T	1	8	15	22	29
W	2	9	16	23	30
T	3	10	17	24	31
F	4	11	18	25	
S	5	12	19	26	
S	6	13	20	27	

September

M		4	11	18	25
T		5	12	19	26
W		6	13	20	27
T		7	14	21	28
F	1	8	15	22	29
S	2	9	16	23	30
S	3	10	17	24	

October

M	30	2	9	16	23
T	31	3	10	17	24
W		4	11	18	25
T		5	12	19	26
F		6	13	20	27
S		7	14	21	28
S	1	8	15	22	29

November

M		6	13	20	27
T		7	14	21	28
W	1	8	15	22	29
T	2	9	16	23	30
F	3	10	17	24	
S	4	11	18	25	
S	5	12	19	26	

December

M		4	11	18	25
T		5	12	19	26
W		6	13	20	27
T		7	14	21	28
F	1	8	15	22	29
S	2	9	16	23	30
S	3	10	17	24	31



VIAMED



15 Station Road · Cross Hills · Keighley · West Yorkshire · BD20 7DT · UK
Tel: +44 (0)1535 634 542 · Fax: +44 (0)1535 635 582
info@viamed.co.uk · www.viamed.co.uk

M Oxygen Sensors

Key:
Main OEM
OEM Sensor Part Number
or Application
Output Connector

Sensor Output
at 1013 mb, 25°C dry ambient air
Response Time
to 90% of final value

0110051 T-1 Teledyne A37016 Slip rings, 2 200 ± 40 µA ≤ 30 s	0110014 R-15V Datex-Ohmeda 0237-2034-700 Slip rings, 3 14.0 - 20.0 mV with 300 Ω load ≤ 12 s	0110017 R-17MED Teledyne C43690 Phone jack, 0.141" 7.0 - 13.0 mV ≤ 6 s	0110018 R-17V Envitec OOM103-1 Phone jack, 3.5 mm 9.0 - 13.0 mV ≤ 5 s	0110019 R-17Vi Viamed Incubators Phone jack, 3.5mm 9.0 - 14.0 mV ≤ 12 s	0110022 R-22MED Teledyne C-44611 Molex, 3 pin plug 7.0 - 13.0 mV ≤ 6 s	0110021 R-22MEDV Dräger 83 620 30 Molex, 3 pin plug 9.0 - 13.0 mV ≤ 5 s
0110039 R-22VA Atom Medical C81828- R22MATOM Molex, 3 pin plug 9.0 - 13.0 mV ≤ 12 s	0110073 R-22VA (Matched Pair) Atom Medical 60859 Molex, 3 pin plug 9.0 - 13.0 mV ≤ 12 s	0110020 R-22Vi Viamed Incubators Molex, 3 pin plug 9.0 - 14.0 mV ≤ 12 s	0110072 R-22Vi (Matched Pair) Viamed Incubators Molex, 3 pin plug 9.0 - 14.0 mV ≤ 12 s	0110023 R-23V Dräger 6850645 Slip rings, 3 24 - 32 µA ≤ 12 s	0110024 R-24MED Siemens 64 19 332 E380E Molex, 3 pin plug 10.0 - 17.0 mV ≤ 10 s	0110026 R-26MED Respironics 8-100412-00 Phone jack, 0.141" 7.0 - 13.0 mV ≤ 10 s
0110040 R-30V GE Healthcare 6050-0004-110 Phone socket, 6P4C 10.0 - 12.0 mV ≤ 12 s	0110041 R-41V Airshields 6735142 Phone socket, 3.5 mm 13.0 - 16.0 mV ≤ 12 s	0110042 R-42V Stephan 130060001 AMP plug, 2 pin 9.0 - 11.5 mV ≤ 12 s	0110043 R-43V Pack of 2 Datex-Ohmeda 6600-1278-600 Slip rings, 2 28.0 - 50.0 mV ≤ 12 s	0110044 R-44V MSA 806572 Molex socket, 3 way 9.0 - 13.0 mV ≤ 12 s	0110045 R-45V Datascope 0600-00-0002 Molex, 3 pin plug 9.0 - 13.0 mV ≤ 12 s	0110047 R-47V Acutronic OOM204 Molex, 3 pin plug 9.0 - 13.5 mV (2 channels) ≤ 12 s
0110048 R-48V City Technology MOX-3 Phone socket, 4 way 10.0 - 12.0 mV ≤ 12 s	0110049 R-49V Siemens 64 19 332 E380E Molex, 3 pin plug 13.0 - 16.0 mV ≤ 12 s	0110057 T-7 Teledyne A51327 Phone socket, 4 way 0.4 ± 0.1 µA ≤ 6 s	0110058 T-7V Envitec OOM105 Phone socket, 4 way 0.4 ± 0.1 µA ≤ 5 s	0110081 UFO 130-2 Spacelabs 010-1120-01 DIL, 10 pin 3.4 ± 0.4 V @ 100% oxygen ≤ 130 ms	0110401 MAX-1 Analytical Industries PSR-11-915 Slip rings, 3 17.0 ± 3.0 mV ≤ 13 s	0110402 MAX-2 Dräger 6850645 Slip rings, 3 12.2 - 18.1 mV ≤ 15 s
0110403 MAX-3 City Technology MOX-3 Phone socket, 4 ways 9.0 - 13.0 mV ≤ 15 s	0110460 MAX-7 Maxtec R109P09 Phone jack, 3.5 mm 8.0 - 12.0 mV ≤ 12 s	0110461 MAX-7L Maxtec R109P13 Phone jack, 3.5 mm 3.5 - 7.5 mV ≤ 12 s	0110409 MAX-9 Teledyne B39636-R13 1 mm pin, x 4 26.4 - 44.0 mV ≤ 20 s	0110411 MAX-11 Dräger Fabius, Evita, Babylog, Cato Slip rings, 3 14.0 - 20.0 mV ≤ 12 s	0110441 MAX-11i Dräger 6850645 Slip rings, 3 14.0 - 20.0 mV ≤ 12 s	0110412 MAX-12 Siemens 9004979-E347E Molex, 3 pin plug 12.5 - 16.5 mV ≤ 12 s
0110462 MAX-12A Philips (Respironics) NeoPAP Molex, 3 pin plug 8.0 - 13.0 mV ≤ 15 s	0110424 MAX-12C Viasys Avea Ventilator Molex, 3 pin plug 8.0 - 13.0 mV ≤ 12 s	0110413 MAX-13 Philips (Respironics) Criterion Vent Phone jack, 3.5 mm 13.0 - 16.5 mV ≤ 12 s	0110463 MAX-13C BioMed Devices 4433 Phone jack, 2.5 mm 13.0 - 17.0 mV ≤ 13 s	0110469 MAX-13CS BioMed Devices Crossvent* Phone jack, 2.5 mm 10.5 - 15.0 mV ≤ 15 s	0110444 MAX-13-250 Maxtec R125P07 Phone jack, 3.5 mm 12.0 - 17.0 mV ≤ 15 s	0110414 MAX-14 MSA 806572 Molex, 3 pin plug 9.0 - 13.0 mV ≤ 13 s
0110420 MAX-14 (Fluke) Fluke R116P82-001 Molex SL, 2 pin 9.0 - 13.0 mV ≤ 12 s	0110416 MAX-16 Nelcor Puritan Bennett 4-072214-00 Molex, 3 pin plug 8.5 - 15.0 mV ≤ 15 s	0110421 MAX-16HT Versamed Vent Molex, 3 pin plug 9.5 - 16.0 mV ≤ 12 s	0110417 MAX-17 Teledyne T-7 Phone socket, 4 way 7.5 - 12.5 mV ≤ 8 s	0110418 MAX-18 Teledyne C-1R Slip rings, 2 200 ± 40 µA ≤ 13 s	0110419 MAX-19 Teledyne T-1 Slip rings, 2 170 - 230 µA ≤ 20 s	0110464 MAX-23 Teledyne R-22MED Molex, 3 pin plug 8.0 - 13.0 mV 1.5 - 3.5 s
0110465 MAX-24 Maxtec R109P15 AMP plug, 2 pin 8.5 - 11.5 mV ≤ 12 s	0110466 MAX-25 Figaro KE-25 Unterminated 10.0 - 15.5 mV 14 ± 2 s	0110449 MAX-43 Pack of 2 GE (Ohmeda) Giraffe Slip rings, 2 25.0 - 35.0 mV ≤ 17 s	0110438 MAX-47 Teleflex (Hudson) 5803 Phone socket, 4 way 4.0 - 11.0 mV ≤ 12 s	0110439 MAX-48 Analytical Industries All-11-60 Phone socket, 4 way 3.8 - 4.8 mV ≤ 12 s	0110450 MAX-50 Figaro KE-50 2 pin socket 47.0 - 65.0 mV 60 ± 5 s	0110470 MAX-108A R109P21 8.0 - 15.0 mV
0110468 MAX-125M Covidien 4-072214-00 Molex, 3 pin plug 10.0 - 15.5 mV ≤ 15 s	0110473 MAX-125WF Fluke R116P82-001 Molex SL, 2 pin 10.5 - 15.0 mV ≤ 15 s	0110428 MAX-250 Maxtec R125P01-002 2 pin socket 10.0 - 15.5 mV ≤ 15 s without 'O' Ring	0110425 MAX-250 Maxtec R125P01-007 2 pin socket 10.0 - 15.5 mV ≤ 15 s with 'O' Ring	0110427 MAX-250+ Maxtec R125P02-011 Molex SL, 2 pin 10.0 - 15.5 mV ≤ 15 s	0110431 MAX-250A Blease 70300001 Phone socket, 4 way 11.0 ± 2.5 mV ≤ 15 s	0110432 MAX-250(B) Maxtec R125P02-003 Unterminated 10.0 - 15.5 mV ≤ 15 s
0110429 MAX-250E Maxtec R125P03-002 Mini power socket 10.0 - 15.5 mV ≤ 15 s	0110430 MAX-250ESF Maxtec R115P85 Mini power socket 12.0 ± 2.0 mV ≤ 5 s	0110433 MAX-250K with DC Viasys 16101 Mini power plug 10.0 - 15.5 mV ≤ 15 s	0110451 MAX-250K with Panduit Viasys Vela T2 Panduit 10.0 - 15.5 mV ≤ 15 s	0110435 MAX-250TM Maxtec R125P18-012 Molex, 3 pin plug, 3P2C 11.0 ± 2.5 mV ≤ 15 s	0110453 MAX-250MS MEK R125P97 Molex 3 pin plug (female) 10.0 - 15.5 mV ≤ 15 s	0110452 MAX-550E Maxtec R140P02, MaxO2 ME Monitor Stereo phone jack, 3.5 mm 10.5 - 15.0 mV ≤ 15 s
0110454 MAX-550E Maxtec MaxBlend Devices Stereo phone jack, 3.5 mm 10.5 - 15.0 mV ≤ 15 s	0110584 OOM103-1M Maxtec MAX-250ESF Mini power socket 9.0 - 13.0 mV ≤ 5 s	0110565 OOM109 Cortex Molex, 3 pin plug 9.0 - 13.0 mV ≤ 360 ms	0110566 OOM109-LF2 Molex, 3 pin plug 9.0 - 13.0 mV ≤ 360 ms	0110560 OOM111 Envitec 01-00-0114 Stereo phone jack, 3.5 mm 11.0 - 13.0 mV ≤ 12 s	0110660 OOI103 Envitec 01-00-0095 Slip rings, 4 9.0 - 15.0 mV ≤ 13 s	0110524 KE-25 Unterminated 10.0 - 15.5 mV ≤ 14 s ± 2 s Other KE sensors also available



Viamed supply a range of Oxygen Monitors and Analysers

Free oxygen sensor disposal service
for the safe disposal of any of your depleted oxygen sensors



Viamed now supplies compatible lead-free sensors. See our online cross reference
Cross reference: www.viamed.co.uk/xref