

**Test Report OOM108 Oxygen Sensor**

|                                                           |                          |                                |                |
|-----------------------------------------------------------|--------------------------|--------------------------------|----------------|
| Test Place:                                               | Sri Jayawardana Hospital | Tested by:                     | H.G. Priyankha |
| <b>Ventilator Master Data</b>                             |                          |                                |                |
| Ventilator type:                                          | maquet servo Air         | Ventilator Serial Number:      | 22407          |
| Ventilator Software Version:                              |                          | Ventilator Manufacturing Date: |                |
| <b>Sensor Master Data (currently used / original one)</b> |                          |                                |                |
| Sensor type:                                              | maquet                   | Sensor Part Number:            | 6640044        |
| Sensor Serial Number:                                     | AA450598                 | Others:                        |                |
| <b>Sensor Master Data (tested sensor)</b>                 |                          |                                |                |
| Sensor type:                                              | OOM108                   | Sensor Part Number:            | E01-00-0075    |
| Sensor Serial Number:                                     | AA 600010                | Others:                        |                |

**Test Assignment**

Check of oxygen sensor OOM108 in MAQUET ventilator regarding to compatibility, usability, installation, calibration and measurement accuracy.

**Checklist**

| Pos. | Check Item                                                                        | Yes | No | Comments |
|------|-----------------------------------------------------------------------------------|-----|----|----------|
| 1    | <b>Sensor usability</b><br>Sensor usability is same as the original one?          | ✓   |    |          |
| 2    | <b>Sensor installation</b><br>Sensor installation is same as the original one?    | ✓   |    |          |
| 3    | <b>Sensor calibration</b><br>Sensor calibration was successful in the ventilator? | ✓   |    |          |
| 4    | <b>Oxygen Measurement</b><br>Sensor works without issues in ventilator?           | ✓   |    |          |


**Final Conclusion:**

The oxygen sensor type OOM108 tested in MAQUET ventilator is compatible to the currently used oxygen sensor made by MAQUET.

☒ Yes ☐ No

Date: 28/08/2023

Signature:



|                                                                                                                      |           |                     |               |
|----------------------------------------------------------------------------------------------------------------------|-----------|---------------------|---------------|
|                                   | Revision: | Created - Erstellt: | Page - Seite: |
| C:\Users\E339512\AppData\Local\Microsoft\Windows\NetCache\Content.Outlook\YEPULFMN\OOM108_compatibility_check-0.docx |           | MS                  | 1 / 1         |

Template: Prüfbericht.dot / Rev. 2 gültig: 05.01.2023