

«COMPANY»
«TITLE»
«DEPARTMENT»
«COMPANY»
«ADDRESS1»
«ADDRESS2»
«ADDRESS3»
«CITY»
«STATE»
«ZIP»



Ref: 4490103-4490109/20180724/500
Date

Dear «UPOSITION»,

Introducing VersaStream Gas Sampling Lines for use with Capnographs and patient monitors with Oridion (Covidien, Nellcor, Medtronic, Philips, Masimo, Mindray, Zoll, Nihon Kohden...) connectors from Viamed

Viamed are pleased to introduce VersaStream gas sampling lines as an innovative solution for sidestream gas monitoring with increased performance and low cost. These lines are compatible to the FilterLine® range of sampling lines and for use with Microstream® enabled capnograph and patient monitors.



VersaStream gas sampling lines use the latest material technologies to allow the removal of moisture and water - even in cases of low temperatures and high humidity. They are very flexible and are suitable for both short-term and long-term applications.

Why VersaStream gas sampling lines?

- One-piece sampling line solution
- **Cost-effective** alternative to Nafion® tubes
- The gas-selective non-porous material facilitates fast and removal of condensation from the sampling line
- The integrated bacterial filter **prevents contamination**
- **Reliable measured gas values**, due to negligible dead space volume of the filter - even in cases of high respiration rates
- **Integrated bacteria filter and water removal** reduces the number of patient circuit parts and speeds up patient monitoring preparation time

For use with Oridion (Medtronic/Covidien/Nellcor) capnographs, as well as capnographs and patient monitors from other manufacturers that incorporate Microstream® technology, such as Philips, Dräger, Physio Control, Zoll...

These lines are available in a number of variants including: single nasal sampling line, nasal cannula, nasal/oral cannula and a one piece airway adapter sampling line. They are compatible with a wide range of OEM equipment, please see our cross reference and selection guides for further information: www.viamed.co.uk/linesxref. Should you be unable to find the equipment you are using, please contact us so that we can search our extensive database.



Should you have any further questions regarding this, or would like to obtain samples for evaluation, please do not hesitate in contacting our sales team on **01535 634542** between 9:00am and 5:30pm Monday – Friday or email sales@viamed.co.uk

Yours sincerely

Catrin

Catrin Hollings
Marketing
Email: sales@viamed.co.uk

Our VersaStream range includes compatibles for the following:

Oridion Nellcor, Covidien, Medtronics

Compatible to the Filterline range of sampling lines and for use with Microstream enabled capnographs



Luer Lock Male

For use with Capnographs and Multigas Monitors that incorporate a Luer Lock sampling port



Viamed Masimo, Phase-In

For use with VM-2500 Sidestream Capnographs and Masimo/Phase-In technology enabled capnographs

