



15 Station Road  
Cross Hills  
Keighley  
West Yorkshire  
BD20 7DT

+44 (0)1535 634542  
[www.viamed.co.uk](http://www.viamed.co.uk)



*OxiQuant S & Handi+*  
Oxygen Analysers

2024



Used to determine oxygen concentrations in gas mixtures, delivered by oxygen concentrators, blenders...

- Flexible manual calibration
- Auto-off (after 3 minutes)
- Response time: <13 seconds to 90% of value
- Sensor (with integrated battery) can be replaced

This analyser is great for spot-checking oxygen percentage in hospitals, clinics and in homecare environments. Includes integrated sensor and battery.



- One-touch calibration
- Auto-off (after 80 seconds)
- Response time: <15 seconds to 90% of value
- Disposable at end of sensor life
- Overmould for extra protection