

# Training Feedback Form

<b>Training Course Completed:</b> Maxtec Oxygen Monitoring Technical Training			
<b>Date:</b> 19/6/23	<b>Time/Length:</b> 45 mins	<b>Trainer:</b> Steve Hardaker	
<b>Content</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Was the course content presented in a logical manner?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Was the course content and material complete and comprehensive?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Will this information be useful to you in your job role?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Relevance</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Do you feel you now have a better understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did the course challenge your thinking and understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Do you feel the training is beneficial to your team?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Trainer</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Did the trainer communicate and explain the material clearly?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did you feel the instructor was knowledgeable in the area covered?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did the trainer encourage discussions and questions?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Comments</b>			
Do you require any further training in this area?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
If so, what would you like this training to cover? <div style="text-align: center;">N/A</div>			
Further comments: <div style="text-align: center;">N/A</div>			
<b>Name:</b> Sophie Lines			
<b>Date:</b> 19-6-2023			

\*Please delete as applicable

Maxtec Oxygen Monitoring Introduction Competency Questions

1. Why would you choose a monitor over an analyser?

✓ A monitor has alarms and can be used in patient monitoring.

2. What is the difference between the A and AE models?

✓ The A model has an internal sensor  
The AE model has an external cabled sensor

3. What is a Ultramax's application?

✓ The Ultramax is used to measure the concentration, flow and pressure of Oxygen concentrators

4. How do Nitrogen analysers work?

✓ By reading the Oxygen in a sample and displaying the remaining number after the Oxygen has been ~~removed~~ subtracted.

5. What is the USD distributor price for MaxO2+AE and part number?

✓ Part number: 0111261 Price: \$337.40

6. When was the UK GBP End User price last increased for the Handi+?

✓ 10 Jan 2023

7. What is the replacement sensor for the Handi+?

✓ The sensor cannot be replaced.

8. What is the warranty of the MaxO2 ME monitor?

✓ Warranty is 24 months

9. What is the difference between the two MAX-550E oxygen sensor models and what products are they each for use with?

✓ 0110452 - Used with MaxO2 ME  
0110454 - Used with MaxBlend services  
The sensors use different flow diverters.

✓ 10. When does the Maxtec CE certificate expire? Date: 26/may/2024

✓ Is this document available to be provided to customers? Yes

# X Did not attend training due to holiday

## Training Feedback Form

<b>Training Course Completed: Maxtec Oxygen Monitoring Introduction</b>			
<b>Date:</b>	<b>Time/Length:</b>	<b>Trainer: Catrin Hollings</b>	
<b>Content</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Was the course content presented in a logical manner?			
Was the course content and material complete and comprehensive?			
Will this information be useful to you in your job role?			
<b>Relevance</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Do you feel you now have a better understanding of the product/procedure/training area*?			
Did the course challenge your thinking and understanding of the product/procedure/training area*?			
Do you feel the training is beneficial to your team?			
<b>Trainer</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Did the trainer communicate and explain the material clearly?			
Did you feel the instructor was knowledgeable in the area covered?			
Did the trainer encourage discussions and questions?			
<b>Comments</b>			
Do you require any further training in this area?			
If so, what would you like this training to cover?			
Further comments:			
Name: Michael Green			
Date:			

\*Please delete as applicable

## Maxtec Oxygen Monitoring Introduction Competency Questions

1. Why would you choose a monitor over an analyser?

✓ You can set Alarms with Monitor, (Not battery powered - DC Power Port Supplied).

2. What is the difference between the A and AE models?

✓ One has a Internal and the other has a external Sensor also Come with a Short cable or standard 3m.

3. What is a Ultramax's application?

is to measure Oxygen Concentration Flow and outlet Pressure of Oxygen Sensors.

4. How do Nitrogen analysers work?

✓ This works as the oxygen diffuses through the membrane a electrical current is generated which is proportional pressure of oxygen in gas sample. The oxygen percentage is then subtracted from 100 - Remainder being displayed as nitrogen.

5. What is the USD distributor price for MaxO2+AE and part number?

Part number: 0111261 Price: \$337.40 ✓

6. When was the UK GBP End User price last increased for the Handi+?

• 24<sup>th</sup> November 2022 ✓

7. What is the replacement sensor for the Handi+?

There is none as it has an Integrated Sensor.

8. What is the warranty of the MaxO2 ME monitor?

24 Months ✓

9. What is the difference between the two MAX-550E oxygen sensor models and what products are they each for use with?

✓ The Max-550E (0110454) is to be used with Maxtec Max blend Range of blenders whereas the 0110452 is used with Max O2 ME monitor. The difference is one of them are equipped with a bleed control which allows users to get more accuracy and waste less gas.

✓ 10. When does the Maxtec CE certificate expire? Date: 26-05-24

✓ Is this document available to be provided to customers? Yes.

# Training Feedback Form

<b>Training Course Completed:</b> Maxtec Oxygen Monitoring Introduction			
<b>Date:</b> 12/6/23	<b>Time/Length:</b> 30 mins	<b>Trainer:</b> Catrin Hollings	
<b>Content</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Was the course content presented in a logical manner?	<input checked="" type="checkbox"/>		
Was the course content and material complete and comprehensive?	<input checked="" type="checkbox"/>		
Will this information be useful to you in your job role?	<input checked="" type="checkbox"/>		
<b>Relevance</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Do you feel you now have a better understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>		
Did the course challenge your thinking and understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>		
Do you feel the training is beneficial to your team?	<input checked="" type="checkbox"/>		
<b>Trainer</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Did the trainer communicate and explain the material clearly?	<input checked="" type="checkbox"/>		
Did you feel the instructor was knowledgeable in the area covered?	<input checked="" type="checkbox"/>		
Did the trainer encourage discussions and questions?	<input checked="" type="checkbox"/>		
<b>Comments</b>			
Do you require any further training in this area?		<input checked="" type="checkbox"/>	
If so, what would you like this training to cover?			
Further comments:			
Name: Agib Majeed			
Date: 12.06.23.			

\*Please delete as applicable

## Maxtec Oxygen Monitoring Introduction Competency Questions

1. Why would you choose a monitor over an analyser?

used for measurements of gases for fractional oxygen

2. What is the difference between the A and AE models?

A has a internal sensor, adapts for tubing applications.  
AE used where an inline analyser is required ventilators / oscillator

3. What is a Ultramax's application?

Helps oxygen concentrators, ~~maxi oxygen~~

4. How do Nitrogen analysers work?

✓ Oxygen percentage is subtracted from 100 result is <sup>nitrogen</sup> ~~percentage~~  
concentrator flow and pressure

5. What is the USD distributor price for MaxO2+AE and part number?

Part number: 33740 ✓ Price: 0111262 ~~maxi~~

6. When was the UK GBP End User price last increased for the Handi+?

• NOV 22 ✓

7. What is the replacement sensor for the Handi+?

NONE ✓ ~~maxi~~

8. What is the warranty of the MaxO2 ME monitor?

36 month x ~~maxi~~

9. What is the difference between the two MAX-550E oxygen sensor models and what products are they each for use with?

one is for MAXO2ME monitor and other is  
maxibend, difference is design and size of flow  
diverters

~~maxi 550E from maxtec~~ ~~maxi 550E from maxtec~~ ~~maxi 550E from maxtec~~  
~~maxi 550E from maxtec~~ ~~maxi 550E from maxtec~~ ~~maxi 550E from maxtec~~

10. When does the Maxtec CE certificate expire? Date: 26/5/2024 ✓

Is this document available to be provided to customers? yes ✓

# Training Feedback Form

<b>Training Course Completed: Maxtec Oxygen Monitoring Introduction</b>			
<b>Date:</b> 12/6/23	<b>Time/Length:</b> 30 min	<b>Trainer:</b> Catrin Hollings	
<b>Content</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Was the course content presented in a logical manner?	<input checked="" type="checkbox"/>		
Was the course content and material complete and comprehensive?	<input checked="" type="checkbox"/>		
Will this information be useful to you in your job role?	<input checked="" type="checkbox"/>		
<b>Relevance</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Do you feel you now have a better understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>		
Did the course challenge your thinking and understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>		
Do you feel the training is beneficial to your team?	<input checked="" type="checkbox"/>		
<b>Trainer</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Did the trainer communicate and explain the material clearly?	<input checked="" type="checkbox"/>		
Did you feel the instructor was knowledgeable in the area covered?	<input checked="" type="checkbox"/>		
Did the trainer encourage discussions and questions?	<input checked="" type="checkbox"/>		
<b>Comments</b>			
Do you require any further training in this area?		<input checked="" type="checkbox"/>	
If so, what would you like this training to cover?			
Further comments:			
Name: <i>EMH</i>			
Date:			

\*Please delete as applicable

skyped

## Maxtec Oxygen Monitoring Introduction Competency Questions

1. Why would you choose a monitor over an analyser?

Monitors have adjustable alarms and may be used for monitoring  $\text{SpO}_2$  gases being delivered to a patient. Analyser is for checking  $\text{SpO}_2$  ✓

2. What is the difference between the A and AE models?

A has internal sensor included, AE has a separate sensor supplied with it ✓

3. What is a Ultramax's application?

Not for use on patients. For use in testing oxygen concentrators; does not need a sensor. Tests oxygen,  $\text{N}_2$  and argon ✓

4. How do Nitrogen analysers work?

It measures the oxygen concentration in the same way as an oxygen analyser, but then subtracts the measured oxygen from 100% and displays the remainder as  $\text{N}_2$  ✓

5. What is the USD distributor price for MaxO2+AE and part number?

Part number: 0111261 Price: 337.40 USD ✓

6. When was the UK GBP End User price last increased for the Handi+?

• 10/01/23 ~~24th/2022 price increased~~  
price not ✓

7. What is the replacement sensor for the Handi+?

No replacement sensor, replace the whole unit when sensor <sup>on battery</sup> depleted ✓

8. What is the warranty of the MaxO2 ME monitor?

24 months ✓

9. What is the difference between the two MAX-550E oxygen sensor models and what products are they each for use with?

P/N 0110452 MAX550E with standard flow diverter for use with MAXO2 ME.  
0110454 MAX550E with blender flow diverter for use with MaxBlend type of blenders. ✓

10. When does the Maxtec CE certificate expire? Date: 26/05/2024 ✓

Is this document available to be provided to customers? Yes ✓



# Training Feedback Form

<b>Training Course Completed:</b> Maxtec Oxygen Monitoring Introduction			
<b>Date:</b> 12/06/23	<b>Time/Length:</b> 40 min	<b>Trainer:</b> Catrin Hollings	
<b>Content</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Was the course content presented in a logical manner?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Was the course content and material complete and comprehensive?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Will this information be useful to you in your job role?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Relevance</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Do you feel you now have a better understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did the course challenge your thinking and understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Do you feel the training is beneficial to your team?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Trainer</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Did the trainer communicate and explain the material clearly?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did you feel the instructor was knowledgeable in the area covered?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did the trainer encourage discussions and questions?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Comments</b>			
Do you require any further training in this area?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
If so, what would you like this training to cover?			
Further comments:			
Name: ROBERT CONNER			
Date: 12/06/23			

\*Please delete as applicable

STYLED

## Maxtec Oxygen Monitoring Introduction Competency Questions

1. Why would you choose a monitor over an analyser?

The monitor has alarms ✓

2. What is the difference between the A and AE models?

The A has an internal sensor while the AE has an external sensor. ✓

3. What is a Ultramax's application?

As test oxygen concentrators. not for patient use. ✓  
oxygen, argon and nitrogen

4. How do Nitrogen analysers work?

The monitor measures oxygen percentage and then remove it from 100% to get the remaining percentage of nitrogen. ✓

5. What is the USD distributor price for MaxO2+AE and part number?

Part number: 0111261 Price: \$337.40 ✓

6. When was the UK GBP End User price last increased for the Handi+?

• Jun 2023 ✓  
~~24th Feb 2024~~ price increased  
parcel list

7. What is the replacement sensor for the Handi+?

you replace once the sensor is depleted ✓

8. What is the warranty of the MaxO2 ME monitor?

24 months ✓

9. What is the difference between the two MAX-550E oxygen sensor models and what products are they each for use with?

max blend, max 550E with flow - 0110454 ✓

max O2 ME, max 550E with standard - 0110452 ✓

10. When does the Maxtec CE certificate expire? Date:

26/05/2024

Is this document available to be provided to customers?

Yes ✓

# Training Feedback Form

<b>Training Course Completed: Maxtec Oxygen Monitoring Introduction</b>			
<b>Date:</b> 12/06/23	<b>Time/Length:</b>	<b>Trainer: Catrin Hollings</b>	
<b>Content</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Was the course content presented in a logical manner?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Was the course content and material complete and comprehensive?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Will this information be useful to you in your job role?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
<b>Relevance</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Do you feel you now have a better understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did the course challenge your thinking and understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Do you feel the training is beneficial to your team?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Trainer</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Did the trainer communicate and explain the material clearly?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did you feel the instructor was knowledgeable in the area covered?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did the trainer encourage discussions and questions?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Comments</b>			
Do you require any further training in this area?			
If so, what would you like this training to cover?			
Further comments:			
Name: <i>Michael Lamb</i>			
Date: 12/06/23			

\*Please delete as applicable

### Maxtec Oxygen Monitoring Introduction Competency Questions

1. Why would you choose a monitor over an analyser?

you would need a monitor if you require alarms ✓

2. What is the difference between the A and AE models?

A - internal sensor AE - external sensor ✓

3. What is a Ultramax's application?

used by service technicians to check ox concentrators ✓

4. How do Nitrogen analysers work?

They measure the oxygen + calculate the amount of gas that is not oxygen. eg nitrogen + maybe other gases ✓

5. What is the USD distributor price for MaxO2+AE and part number?

Part number: 0111261 Price: \$337.40 ✓

6. When was the UK GBP End User price last increased for the Handi+?

• 24.11.22 ✗

7. What is the replacement sensor for the Handi+?

none - it is disposable at the end of sensor life ✓

8. What is the warranty of the MaxO2 ME monitor?

24 months from shipping/invoice ✓

9. What is the difference between the two MAX-550E oxygen sensor models and what products are they each for use with?

0110452 (R140P02) is for MAXO2 ME ox monitors - supplied with standard flow diverter (0120095)

0110454 (R140P02-001) is for MaxBlend2/Lite/Blender Buddy 2 the flow div is narrower with a small bleed hole (0120086) ✓

10. When does the Maxtec CE certificate expire? Date: 26-5-2024

Is this document available to be provided to customers? YES ✓

# Training Feedback Form

<b>Training Course Completed:</b> Maxtec Oxygen Monitoring Introduction			
<b>Date:</b> 12-6-23	<b>Time/Length:</b> 30 mins	<b>Trainer:</b> Catrin Hollings	
<b>Content</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Was the course content presented in a logical manner?	✓		
Was the course content and material complete and comprehensive?	✓		
Will this information be useful to you in your job role?	✓		
<b>Relevance</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Do you feel you now have a better understanding of the product/procedure/training area*?	✓		
Did the course challenge your thinking and understanding of the product/procedure/training area*?	✓		
Do you feel the training is beneficial to your team?	✓		
<b>Trainer</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Did the trainer communicate and explain the material clearly?	✓		
Did you feel the instructor was knowledgeable in the area covered?	✓		
Did the trainer encourage discussions and questions?	✓		
<b>Comments</b>			
Do you require any further training in this area?	✓	no	more to follow ...
If so, what would you like this training to cover?			
Further comments:			
Name: KATE GRIFATHS			
Date: 12-6-23			

\*Please delete as applicable

**Maxtec Oxygen Monitoring Introduction Competency Questions**

1. Why would you choose a monitor over an analyser?

It has adjustable alarms ✓

2. What is the difference between the A and AE models?

A model has an internal sensor ✓  
AE model has an external sensor ✓

3. What is a Ultramax's application?

for checking Oxygen concentration, flow, pressure of Oxygen concentrators - technicians use them ✓

4. How do Nitrogen analysers work?

Not for patient use ✓

The display shows whatever remaining percentage of gas (presumably nitrogen) remains after the oxygen percentage has been subtracted from 100. ✓

5. What is the USD distributor price for MaxO2+AE and part number?

Part number: 0111261 Price: \$337.40 ✓

6. When was the UK GBP End User price last increased for the Handi+?

• 24/11/2022 X

7. What is the replacement sensor for the Handi+?

whole unit is

There isn't a replacement - disposable at end of sensor life. ✓

8. What is the warranty of the MaxO2 ME monitor?

24 months from invoice ✓

9. What is the difference between the two MAX-550E oxygen sensor models and what products are they each for use with?

The difference is the design and size of the flow diverter supplied with each. One is for the MaxO2 ME monitor, the other is for the MaxBlend. ✓

10. When does the Maxtec CE certificate expire? Date: 26/5/2024 ✓

Is this document available to be provided to customers? YES ✓

# Training Feedback Form

<b>Training Course Completed:</b> Maxtec Oxygen Monitoring Introduction			
<b>Date:</b> 12/6/2023	<b>Time/Length:</b> 30 mins	<b>Trainer:</b> Catrin Hollings	
<b>Content</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Was the course content presented in a logical manner?	✓		
Was the course content and material complete and comprehensive?	✓		
Will this information be useful to you in your job role?	✓		
<b>Relevance</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Do you feel you now have a better understanding of the product/procedure/training area*?	✓		
Did the course challenge your thinking and understanding of the product/procedure/training area*?	✓		
Do you feel the training is beneficial to your team?	✓		
<b>Trainer</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Did the trainer communicate and explain the material clearly?	✓		
Did you feel the instructor was knowledgeable in the area covered?	✓		
Did the trainer encourage discussions and questions?	✓		
<b>Comments</b>			
Do you require any further training in this area?	✓		
If so, what would you like this training to cover? SH to provide technical training			
Further comments:			
Name: Sophie Lines			
Date: 12/6/2023			

\*Please delete as applicable

Maxtec Oxygen Monitoring Introduction Competency Questions

1. Why would you choose a monitor over an analyser?

Alarms ✓

2. What is the difference between the A and AE models?

AE have external sensors. ✓

3. What is a Ultramax's application?

measures different gases. oxygen concentration

4. How do Nitrogen analysers work?

The analysers measure oxygen. The oxygen percentage is subtracted from 100 with the remainder being displayed as percent ~~any~~ nitrogen. ✓

5. What is the USD distributor price for MaxO2+AE and part number?

Part number: 0111261 Price: \$337.40 ✓

6. When was the UK GBP End User price last increased for the Handi+?

• 0111250 - history 24.11.22 X

7. What is the replacement sensor for the Handi+?

No replacement ✓

8. What is the warranty of the MaxO2 ME monitor?

24 months ✓

9. What is the difference between the two MAX-550E oxygen sensor models and what products are they each for use with?

• External sensor / non external sensor

• Size of flow diverter

0110452 standard flow diverter for MaxO2 ME  
0110454 flow diverter with hde for MaxRend

10. When does the Maxtec CE certificate expire? Date: May 2024 ✓

Is this document available to be provided to customers? yes. ✓



# Training Feedback Form

<b>Training Course Completed:</b> Maxtec Oxygen Monitoring Introduction			
<b>Date:</b> 12.6.23.	<b>Time/Length:</b> 30 mins	<b>Trainer:</b> Catrin Hollings	
<b>Content</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Was the course content presented in a logical manner?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Was the course content and material complete and comprehensive?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Will this information be useful to you in your job role?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Relevance</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Do you feel you now have a better understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did the course challenge your thinking and understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Do you feel the training is beneficial to your team?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Trainer</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Did the trainer communicate and explain the material clearly?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did you feel the instructor was knowledgeable in the area covered?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did the trainer encourage discussions and questions?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Comments</b>			
Do you require any further training in this area?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
If so, what would you like this training to cover?			
Further comments:			
Name: Jamie Gill			
Date: 12.6.23.			

\*Please delete as applicable

## Maxtec Oxygen Monitoring Introduction Competency Questions

1. Why would you choose a monitor over an analyser?

MONITOR WAS ACCURATE ✓

2. What is the difference between the A and AE models?

A IS INTERNAL AE IS EXTERNAL SENSOR ✓

3. What is a Ultramax's application?

ITS FOR TECHNICIANS TO CHECK OXYGEN CONCENTRATION ✓

4. How do Nitrogen analysers work?

THEY TEST THE OXYGEN CONCENTRATION AND CALCULATE  
FROM 100% THEN DISPLAY THE REMAINDER ✓

5. What is the USD distributor price for MaxO2+AE and part number?

Part number: EN11261 Price: \$337.00 ✓

6. When was the UK GBP End User price last increased for the Handi+?

• 10.1.2023 ✓

7. What is the replacement sensor for the Handi+?

THERE IS NO REPLACEMENT SENSOR ✓

8. What is the warranty of the MaxO2 ME monitor?

24 MONTHS ✓

9. What is the difference between the two MAX-550E oxygen sensor models and what products are they each for use with?

THE FLOW DIVERTER IS DIFFERENT

ONE IS FOR MAXO2 ME THE OTHER 'TWO IN ONE'

4 FOR USE WITH MAXO2 ME ✓

10. When does the Maxtec CE certificate expire? Date: 26.5.2024 ✓

Is this document available to be provided to customers? YES ✓

# Training Feedback Form

<b>Training Course Completed:</b> Maxtec Oxygen Monitoring Introduction			
<b>Date:</b> 12.6.23	<b>Time/Length:</b> 60 min	<b>Trainer:</b> Catrin Hollings	
<b>Content</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Was the course content presented in a logical manner?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Was the course content and material complete and comprehensive?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Will this information be useful to you in your job role?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Relevance</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Do you feel you now have a better understanding of the product/procedure/training area*?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did the course challenge your thinking and understanding of the product/procedure/training area*?	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Do you feel the training is beneficial to your team?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Trainer</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Did the trainer communicate and explain the material clearly?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did you feel the instructor was knowledgeable in the area covered?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Did the trainer encourage discussions and questions?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Comments</b>			
Do you require any further training in this area?			
If so, what would you like this training to cover?			
Further comments:			
Name: Phil Crossley			
Date: 12.6.23			

\*Please delete as applicable

Skyped

### Maxtec Oxygen Monitoring Introduction Competency Questions

**1. Why would you choose a monitor over an analyser?**

For constant monitoring, as apposed to spot checking, where high and low alarms are required ✓

**2. What is the difference between the A and AE models?**

The A has an internal sensor and the AE has an external sensor attached via a cable ✓

**3. What is a Ultramax's application?**

For testing Oxygen Concentrators ✓

**4. How do Nitrogen analysers work?**

The same way that an Oxygen analyser works, but it displays the non oxygen content of the gas ✓

**5. What is the USD distributor price for MaxO2+AE and part number?**

Part number: 0111261

Price: \$337.40 ✓

**6. When was the UK GBP End User price last increased for the Handi+?**

- 10/01/2023

*24th Nov 23*

**7. What is the replacement sensor for the Handi+?**

It does not have a replaceable sensor ✓

**8. What is the warranty of the MaxO2 ME monitor?**

24 months from date of invoice ✓

**9. What is the difference between the two MAX-550E oxygen sensor models and what products are they each for use with?**

0110452: supplied with a standard flow divertor for use with the MAXO2 ME ✓

0110454: is supplied with a special flow divertor specifically for use with the the Maxblend devices

*bleed hole*

**10. When does the Maxtec CE certificate expire? Date: 26/05/2024** ✓

Is this document available to be provided to customers? .Yes ✓

X

DDrdr  
after are to annual leave

## Training Feedback Form

<b>Training Course Completed: Maxtec Oxygen Monitoring Introduction</b>			
<b>Date:</b> _____	<b>Time/Length:</b> _____	<b>Trainer:</b> Catrin Hollings	
<b>Content</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Was the course content presented in a logical manner?			
Was the course content and material complete and comprehensive?			
Will this information be useful to you in your job role?			
<b>Relevance</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Do you feel you now have a better understanding of the product/procedure/training area*?			
Did the course challenge your thinking and understanding of the product/procedure/training area*?			
Do you feel the training is beneficial to your team?			
<b>Trainer</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Did the trainer communicate and explain the material clearly?			
Did you feel the instructor was knowledgeable in the area covered?			
Did the trainer encourage discussions and questions?			
<b>Comments</b>			
Do you require any further training in this area?			
If so, what would you like this training to cover?			
Further comments:			
Name: Ryan Swaine			
Date: 13/06/2023			

\*Please delete as applicable

Maxtec Oxygen Monitoring Introduction Competency Questions

1. Why would you choose a monitor over an analyser?

✓ A Monitor has Settable/adjustable Alarms and can be used for patient Monitoring.

2. What is the difference between the A and AE models?

✓ A Model - has a Internal Sensor  
AE Model has a external radial sensor

3. What is a Ultramax's application?

✓ Ultramax is used to measure the Concentration flow and pressure of oxygen concentrators.

4. How do Nitrogen analysers work?

✓ By Reading the oxygen in a sample and displaying the remaining number after the oxygen has been subtracted

✓ 5. What is the USD distributor price for MaxO2+AE and part number?

Part number: 0111261 Price: \$337.40

6. When was the UK GBP End User price last increased for the Handi+?

• 10th Jan 2023 ~~10th Jan 2022~~ ~~Price not~~ ~~adjusted~~ ~~generated at time~~ ~~of sale~~

✓ 7. What is the replacement sensor for the Handi+?

✓ ~~AE~~ Sensor is not replaceable.

✓ 8. What is the warranty of the MaxO2 ME monitor?

✓ 24 months

9. What is the difference between the two MAX-550E oxygen sensor models and what products are they each for use with?

✓ 0110452 - to be used with MaxO2 me  
0110454 To be used with Maxblend devices.  
The sensors have different flow diverts

✓ 10. When does the Maxtec CE certificate expire? Date: 26/5/24

✓ Is this document available to be provided to customers? yes

Did not attend.

# X Training Feedback Form

<b>Training Course Completed:</b> Maxtec Oxygen Monitoring Introduction			
<b>Date:</b>	<b>Time/Length:</b>	<b>Trainer:</b> Catrin Hollings	
<b>Content</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Was the course content presented in a logical manner?			
Was the course content and material complete and comprehensive?			
Will this information be useful to you in your job role?			
<b>Relevance</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Do you feel you now have a better understanding of the product/procedure/training area*?			
Did the course challenge your thinking and understanding of the product/procedure/training area*?			
Do you feel the training is beneficial to your team?			
<b>Trainer</b>	<b>Yes</b>	<b>No</b>	<b>Unsure</b>
Did the trainer communicate and explain the material clearly?			
Did you feel the instructor was knowledgeable in the area covered?			
Did the trainer encourage discussions and questions?			
<b>Comments</b>			
Do you require any further training in this area?			
If so, what would you like this training to cover?			
Further comments:			
Name: Cathy Green			
Date:			

\*Please delete as applicable