

# Posey Products Revision Questions (Edition 1.0)

Please review the Introduction (DOCID 116642) and Technical (DOCID 168661) training notes in the document index and answer the following questions.

Email \*

viamed.kate.griffiths@gmail.com

## Why are Posey pulse oximetry sensor wraps not used with disposable sensors? \*

- They are not adjustable.
- Disposable sensors have self-adhesive tape.
- They cause skin damage.
- They are only for paediatric use.

## Where can Posey sensor wraps be applied? \*

- Neonatal foot, infant finger, and adult wrist
- Neonatal big toe, infant foot, and paediatric or adult finger
- Adult wrist, neonatal hand, and paediatric ankle
- Paediatric toe, adult arm, and neonatal hand

**What method is used to secure Posey I.D. bracelets? \***

- Adhesive tape
- Metal clasp
- Hook and loop fasteners
- Elastic band

**Why can Posey SecureSleeves remain in place during X-ray procedures? \***

- They contain metal supports.
- They are X-ray translucent.
- They are made of plastic.
- They are disposable.

**What is the primary purpose of the Posey SecureSleeve? \***

- To secure tracheostomy tubes
- To protect against falls
- To reduce limb contractures and prevent movement during IV procedures
- To support pulse oximetry sensors

**Can Posey I.V. armboards be reused? \***

- Yes, but only for the same patient after cleaning.
- No, they are strictly single-use.
- Yes, they are designed for multiple patients.
- No, they are disposable after first use.

**Which of the following additional feature does the Posey Secure Trach Tie have? Select all correct answers. \***

- Adjustable length for comfort
- A second tie for added security
- Waterproof material
- Latex-free design

**Do any Posey products supplied by Viamed contain latex? \***

- Yes, all Posey products contain latex.
- No, all Posey products are latex-free.
- Only certain splints contain latex.
- Only Secure Trach Ties contain latex.

**What sizes are available for Posey Foam Trach Ties? \***

- Infant and paediatric only
- Small, medium, large, and extra-large
- Neonatal, paediatric, and adult
- Small and medium only

**Who is the legal manufacturer or the Posey products? \***

- Posey Ltd.
- Teledyne Analytical Industries
- Tidi Products LLC
- Viamed Ltd.

This content is neither created nor endorsed by Google.

Google Forms