

VIAMED

Suppliers of Electromedical Equipment

2014



Dowley Gap, Leeds Liverpool Canal, Bingley

January

M	6	13	20	27
T	7	14	21	28
W	1	8	15	22 29
T	2	9	16	23 30
F	3	10	17	24 31
S	4	11	18	25
S	5	12	19	26

February

M	3	10	17	24
T	4	11	18	25
W	5	12	19	26
T	6	13	20	27
F	7	14	21	28
S	1	8	15	22
S	2	9	16	23

March

M	31	3	10	17	24
T		4	11	18	25
W		5	12	19	26
T		6	13	20	27
F		7	14	21	28
S	1	8	15	22	29
S	2	9	16	23	30

April

M		7	14	21	28
T	1	8	15	22	29
W	2	9	16	23	30
T	3	10	17	24	
F	4	11	18	25	
S	5	12	19	26	
S	6	13	20	27	

May

M	5	12	19	26
T	6	13	20	27
W	7	14	21	28
T	1	8	15	22 29
F	2	9	16	23 30
S	3	10	17	24 31
S	4	11	18	25

June

M	30	2	9	16	23
T		3	10	17	24
W		4	11	18	25
T		5	12	19	26
F		6	13	20	27
S		7	14	21	28
S	1	8	15	22	29

July

M		7	14	21	28
T	1	8	15	22	29
W	2	9	16	23	30
T	3	10	17	24	31
F	4	11	18	25	
S	5	12	19	26	
S	6	13	20	27	

August

M		4	11	18	25
T		5	12	19	26
W		6	13	20	27
T		7	14	21	28
F	1	8	15	22	29
S	2	9	16	23	30
S	3	10	17	24	31

September

M	1	8	15	22	29
T	2	9	16	23	30
W	3	10	17	24	
T	4	11	18	25	
F	5	12	19	26	
S	6	13	20	27	
S	7	14	21	28	

October

M		6	13	20	27
T		7	14	21	28
W	1	8	15	22	29
T	2	9	16	23	30
F	3	10	17	24	31
S	4	11	18	25	
S	5	12	19	26	

November

M		3	10	17	24
T		4	11	18	25
W		5	12	19	26
T		6	13	20	27
F		7	14	21	28
S	1	8	15	22	29
S	2	9	16	23	30

December

M	1	8	15	22	29
T	2	9	16	23	30
W	3	10	17	24	31
T	4	11	18	25	
F	5	12	19	26	
S	6	13	20	27	
S	7	14	21	28	



VIAMED



15 Station Road · Cross Hills · Keighley · West Yorkshire · BD20 7DT · UK
Tel: +44 (0)1535 634 542 · Fax: +44 (0)1535 635 582
info@viamed.co.uk · www.viamed.co.uk



Oxygen Sensors

- Wide range of sensors available
- Sensors that utilise the latest technology
- High performance
- Reliable and accurate

0110001 C-1R Teledyne B39022 1 mm pin, x 2 200 ± 40 µA < 30 s	0110002 C-2R Teledyne B39023 1 mm pin, x 2 200 ± 40 µA < 30 s	0110013 R-13 Hudson RCI 5500 1 mm pins, x 4 0.4 ± 0.1 µA < 6 s	0110015 R-15 Datex-Ohmeda 0237-2034-700 Slip rings, 3 59 ± 11 µA < 15 s	0110016 R-17MEDV Teledyne C43690 Phone jack, 0.141" 9 – 13 mV < 6 s
0110017 R-17MED Teledyne C43690 Phone jack, 0.141" 7 – 13 mV < 6 s	0110021 R-22MEDV Dräger 83 620 30 Molex, 3 pin plug 9 – 13 mV < 6 s	Disposal service available for depleted oxygen sensors.	0110071 R-22MEDV (Matched Pair) Dräger 83 620 30 Molex, 3 pin plug 9 – 13 mV < 6 s	0110022 R-22MED Viasys C44611 Molex, 3 pin plug 7 – 13 mV < 6 s
0110023 R-23V Dräger 6850645 Slip rings, 3 24 – 32 µA < 12 s	0110024 R-24MED Siemens 64 19 332 E380E Molex, 3 pin plug 10 – 17 mV < 10 s		0110026 R-26MED Respironics 8-100412-00 Phone jack, 0.141" 7 – 13 mV < 10 s	0110034 R-34MED Maxtec R115P85 Mini power socket 12 ± 2 mV < 5 s
0110040 R-30V Datex-Ohmeda 6050-0004-110 Phone, 6P4C 9 – 12 mV < 12 s	0110041 R-41V Airshields 6735142 Phone jack, 3.5 mm 13 – 16 mV < 12 s	0110042 R-42V Stephan 130060001 AMP plug, 2 pin 9 – 11.5 mV < 12 s	0110043 R-43V Datex-Ohmeda 6600-1278-600 Slip rings, 2 23 – 35 mV < 12 s	0110044 R-44V MSA 806572 Molex, 3 pin 9 – 13 mV < 12 s
0110045 R-45V Datascope 0600-00-0002 Molex, 3 pin plug 9 – 13 mV < 12 s	0110045 R-47V Acutronic OOM204 Molex, 3 pin plug 9 – 13.5 mV (2 channels) < 12 s	For additional Oxygen Sensors see Viamed's Short Form Cross Reference www.viamed.co.uk Phone Viamed if your Oxygen Sensor or equipment is not listed		0110051 T-1 Teledyne A37016 Slip rings, 2 200 ± 40 µA < 30 s
0110054 T-4 Teledyne A37017 1 mm plug, x 3 0.4 ± 0.2 µA < 6 s	0110057 T-7 Teledyne A51327 Phone, 4 way 0.4 ± 0.1 µA < 6 s	0110080 UFO 130-1 Spacelabs 1010-1120-00 DIL, 10 pin 3.4 ± 0.4 V < 130 ms	0110081 UFO 130-2 Spacelabs 1010-1120-01 DIL, 10 pin 3.4 ± 0.4 V < 130 ms	0110403 MAX-3 City Technology MOX-3 Phone, 4 way 9 – 13 mV < 15 s
0110416 MAX-16 Nellcor PB 4-072214-00 Molex, 3 pin plug 8.5 – 15.0 mV < 15 s	0110420 MAX-14 (Fluke) Fluke R116P82-001 Molex SL, 2 pin 9 – 13 mV < 12 s	0110428 MAX-250 without 'O' Ring Maxtec R125P01-002 2 pin socket 10.0 – 15.5 mV < 15 s	0110425 MAX-250 with 'O' Ring Maxtec R125P01-007 2 pin socket 10.0 – 15.5 mV < 15 s	0110427 MAX-250+ Maxtec R125P02-011 Molex SL, 2 pin 10.0 – 15.5 mV < 15 s
0110429 MAX-250E Viasys 66045 Mini power socket 10.0 – 15.5 mV < 15 s	0110430 MAX-250ESF Maxtec R115P85 Mini power socket 12 ± 2 mV < 5 s	Viamed offer a range of Oxygen Analysers and Monitors.		0110431 MAX-250A Blease 70300001 Phone socket, 4 way 11 ± 2.5 mV < 15 s
0110433 MAX-250K Viasys 15974 Mini power plug 10.0 – 15.5 mV < 15 s	0110436 MAX-250(B)-SLE SLE N2191 2 pole socket 10.0 – 15.5 mV < 15 s	0110437 MAX-250(B)-SLE SLE N2191/10 3 pole socket 10.0 – 15.5 mV < 15 s	0110524 KE25 Viamed 100601 Unterminated 10.0 – 15.5 mV 14 ± 2 s	0110525 KE25 Siemens C79451-A3458-B55 2 pole socket 10.0 – 15.5 mV 14 ± 2 s



VIAMED

Key Main OEM
OEM Sensor Part Number
Output Connector
Sensor Output at 1013 mb, 25°C dry ambient air
Response Time to 90% of final value

Tel: 01535 634 542

Specifications Subject to Change

Part Number: 0196000

