




SpO₂ sensors and adapter cables for IntelliVue / CMS / Viridia



SpO₂ monitoring

Non-U.S. version

PHILIPS

Philips FAST SpO₂ Option A01 and Philips Sensors — 8-pin connector

With these monitors, the following sensors and adapters are validated



M1020B



M3000A/M3001A MMS



M3002A MMS X2



M8102A MP2



M8105A MP5

Optional extension cable

Sensors plug directly into instrument. Extension cable M1941A (2 m) can be used to add length to the A/B series sensors.



Adapter cable required

Adapter cable M1943A (1.1 m) or M1943AL (3 m)



Connect to these sensors



M1191A/B Adult
(2 m)



M1192A Pediatric
(1.5 m)



M1193A Neonatal
(1.5 m)



M1194A Ear Clip
(1.5 m)



M1195A Infant
(1.5 m)

Do not use extension cable



M1191AL/BL (3 m) Adult
Extra long sensor cable length




M1196A (3 m) Pediatric/Adult Clip
Extra long sensor cable length




M1191T Adult
(0.45 m)



M1192T Pediatric
(0.45 m)



M1193T Neonatal
(0.90 m)



M1196T Pediatric/
Adult Clip
(0.90 m)



M1131A Adult/
Pediatric
(0.45 m)



M1132A Infant
(0.90 m)



M1133A Neonatal/
Adult/Infant
(0.90 m)

■ Reusable

■ Single-patient-use

Philips FAST SpO₂ Option A01 and Masimo® Sensors

With these monitors, the following sensors and adapters are validated



M1020B



M3000A/M3001A MMS



M3002A MMS X2



M8102A MP2



M8105A MP5

Adapter cable required

Masimo adapter cable
M1020-61100 (3.6 m) for LNOP Sensors



Adapter cable required

Masimo adapter cable
9898031 48221 (3.0 m) for LNCS Sensors



Connect to these LNOP® Sensors



9898031 40321
Adult (DCI)



9898031 40331
Pediatric (DCIP)



9898031 40341
Adult (TC-I)



9898031 40231
Adult (Adt)



9898031 40261
Pediatric (Pdt)



9898031 40311
Infant (Inf-L)



9898031 40291
Neonatal (Neo-L)



9898031 40301
Neonatal (Neo Pt-L)

Connect to these LNCS® Sensors



9898031 48281
Adult (DC-I)



9898031 48291
Pediatric (DC-IP)



9898031 48301
Adult (TC-I)



9898031 48231
Adult (Adtx)



9898031 48241
Pediatric (Pdtx)



9898031 48251
Infant (Inf-L)



9898031 48271
Neonatal (Neo-L)



9898031 48261
Neonatal (Neo Pt-L)

Masimo SET® Technology and Masimo Sensors

With the Masimo SET-based module and monitors, the following sensors and adapters are validated:



SpO₂ Module
M1020B Opt.A03



Philips MMS X2
M3002 Opt.A03



Philips MMS
M3001A Opt.A03

Adapter cable required

Masimo adapter cable
M1020-61100 (3.6 m) for LNOP Sensors



Adapter cable required

Masimo adapter cable
9898031 48221 (3.0 m) for LNCS Sensors




Connect to these LNOP Sensors



9898031 40321
Adult (DCI)



9898031 40331
Pediatric (DCIP)




9898031 40341
Adult (TC-I)



9898031 40231
Adult (Adt)



9898031 40261
Pediatric (Pdt)



9898031 40311
Infant (Inf-L)



9898031 40291
Neonatal (Neo-L)



9898031 40271
Neonatal (Neo)



9898031 40301
Neonatal (Neo Pt-L)



9898031 40281
Neonatal (Neo Pt)

Connect to these LNCS Sensors



9898031 48281
Adult (DC-I)



9898031 48291
Pediatric (DC-IP)



9898031 48301
Adult (TC-I)



9898031 48311
Adult (TF-I)




9898031 48231
Adult (Adtx)



9898031 48241
Pediatric (Pdtx)



9898031 48251
Infant (Inf-L)



9898031 48261
Neonatal (Neo Pt-L)



9898031 48271
Neonatal (Neo-L)

Philips SpO₂ CMS / Viridia >= Rev. C (legacy HP/Agilent) 12-pin input connector

With this module, the following sensors and adapters are validated



M1020A Module

Adapter cable required

Adapter cable M1940A (2 m)
(12-pin to 8-pin connector)



Adapter cable required

Adapter cable M1900B (3 m)
(12-pin to 9-pin connector)



Adapter cable required

Adapter cable M1020-61102 (3.6 m)
for Masimo LNOP Sensors



Connect to these sensors



M1191A/B Adult
(2 m)

M1192A Pediatric
(1.5 m)



M1193A Neonatal
(1.5 m)

M1194A Ear Clip
(1.5 m)



M1195A Infant
(1.5 m)

Do not use extension cable



M1191AL/BL (3 m)
Adult
Extra long sensor
cable length

M1196A (3 m)
Pediatric/Adult Clip
Extra long sensor
cable length



M1191T Adult
(0.45 m)

M1192T Pediatric
(0.45 m)



M1193T Neonatal
(0.90 m)

M1196T Pediatric/
Adult Clip
(0.90 m)



M1131A Adult/
Pediatric
(0.45 m)

M1132A Infant
(0.90 m)




M1133A Neonatal/
Adult/Infant
(0.90 m)



9898031 40321
Adult (DCI)



9898031 40331
Pediatric (DCIP)



9898031 40341
Adult (TC-I)



9898031 40231
Adult (Adt)



9898031 40261
Pediatric (Pdt)



9898031 40311
Infant (Inf-L)



9898031 40291
Neonatal (Neo-L)



9898031 40301
Neonatal (Neo Pt-L)

Philips FAST SpO₂ Option A01 and Nellcor Sensors — 8-pin connector

With these monitors, the following sensors and adapters are validated



M1020B



M3000A/M3001A MMS



M3002A MMS X2



M8102A MP2



M8105A MP5

Adapter cable required



Adapter cable M1943A (1.1 m)
or M1943AL (3 m)

Connect to these adhesive sensors (only available from Philips in Europe and Japan)



M1901B
Neonatal foot



M1902B
Infant digit



M1903B
Pediatric digit



M1904B
Adult digit

Philips SPO₂ CMS / Viridia >= Rev. C (legacy HP/Agilent) 12-pin input connector

With these monitors, the following sensors and adapters are validated



M1020A Module

Adapter cable required



M1900B (3 m)

Connect to these adhesive sensors (only available from Philips in Europe and Japan)



M1901B
Neonatal foot



M1902B
Infant digit



M1903B
Pediatric digit



M1904B
Adult digit



© 2008 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com
fax: +31 40 27 64 887

Printed in The Netherlands
4522 962 39921 * NOV 2008

Philips Healthcare
Global Information Center
P.O. Box 1286
5602 BG Eindhoven
The Netherlands