

## **Resuscitation Cabinet Installation Instructions v9.0 Update Log**

Updated to v9.0, 14/01/14 – SJH

Large version number change from v8.1 due to Nufer specifying a completely new mounting arm for the radiant warmer.

> Headers revised to "Version number to 9.0 January 2014"

> Radiant warmer image to Ceratherm 600-3 (note: composite image with 600-3 head unit on old mounting arm as we have no image of new arm as yet - this has been requested of Daniel Urfer).

> Copyright date to 2014.

> Amended index.

> Replaced 'Overall Installation' diagram and updated 'Key Installation Measurements'.

> Reworded section: 'Mounting Heights' and 'Overhead Clearance' for clarification, including adding the statement: "If the minimum overhead clearance cannot be maintained due to the height of the ceiling, the 50 mm safety margin built into the distance between the top of the mattress and the underside of the radiant warmer head-unit can be sacrificed, providing the minimum distance of 800 mm between the top of the mattress and the underside of the radiant warmer is maintained."

> Added section: 'Mounting the System into Bespoke Furniture', including a recommendation to consider installing an additional electrical switch on the cupboard door mechanism to cut off the power to the radiant warmer when the cupboard door is closed. Also that the minimum cupboard depth required is dictated by the radiant warmer, not the cabinet.

> Edited and replaced drawing 'Cabinet Dimensions; Viewed from Rear' as the lines were hairline thickness and not printing clearly.

> Added model number and manufacturer of the mounting arm.

> Updated section 'Weights of Component Parts' to detail the weight of the new arm and the different distance from the wall when fully extended.

> Added section: 'Important Information Concerning Installation of the Radiant Warmer', to collate cautionary statements. Added a statement advising to refer to the separate Installation Guide for the arm prior to attempting installation.

> Re-drew and replaced the drawing 'Radiant Warmer Dimensions; Viewed from Side, Arm Extended' with new arm.

> Re-drew and replaced the drawing 'Radiant Warmer Dimensions; Viewed from Front, Arm Folded' with new arm.

> Drew and added the drawing 'Radiant Warmer Dimensions; Viewed from Side, Arm Folded' with new arm.

- > Amended section: 'Positioning the Warmer' to recommend the mounting arm be mounted with the bracket offset to the left so as to position the head unit centrally above the cabinet.
- > Removed the following sections as now recommending referring directly to the Installation Guide for the mounting arm:
  - Wall Bracket Mounting Instructions
  - Mounting the Arm onto the Wall Bracket
  - Mounting Arm and Pivot
  - Connecting the Mounting Arm Assembly to the Warmer
- > Drew and added 'Wall Bracket Back Plate Dimensions'
- > Added section: 'Connecting the Radiant Warmer to the Mounting Arm', including the statement: "If local regulations require it, perform an electrical safety test."
- > Replaced image of Therapy Equipment 4725 low suction controller with newer 7725 low suction controller.
- > Added statement regarding placement of the laminated Tom Thumb instruction cards in the adhesive wallet, replacing statement concerning instructions on a keychain.
- > Consumables And Parts List amended:
  - Updated part number for Therapy Equipment low suction controller to newer 7725 version.
  - Updated part number for Apgar timer bracket to rev.2.
  - Added part number for Filter for Therapy Equipment 7725 series low suction controller, single unit
  - Filters for Therapy Equipment 7725 series low suction controller, box of 10.
  - Added part number for NeoPEEP variant PEEP circuit, single patient use (without mask), box of 20.
  - Added part number for Single use silicone face mask, size 00, pack of 20.
  - Added part number for Single use silicone face mask, size 01, pack of 20.
  - Added part number for Disposable adapter for connecting NeoPEEP to Tom Thumb, single unit.
  - Added part number for Disposable adapter for connecting NeoPEEP to Tom Thumb, pack of 20.
  - Added part number for Ceratherm 600-3 operator's manual.
  - Added part number for Ceratherm 600-3 service manual.
  - Added part number for Ceratherm 600-3 supporting documents on CD.

Future revisions:

Cabinet image on front cover still has old Therapy Equipment suction controller, need to replace with new Therapy Equipment version.

Radiant warmer image on front cover has Ondal mounting arm, need to replace with new version once we have it: this has been requested from Nufer.