

[Click here to enable desktop notifications for Gmail.](#) [Learn more](#) [Hide](#)

Gmail

COMPOSE

Inbox (53)

Starred

Important

Sent Mail

Drafts (2)

Albayan Supply

AML

AML Group

Angel Healthcare...

Aramex

Armstrong Medical

Benson Medical

Bestman

Bio-Provider

Bluepoint Medic...

Bunzl (1)

Caravan Parts

[Milwaukee UK - www.uk.milwaukeetool.eu](#) - Sign Up & Win- Fuel Kit Worth £1500 Driven To Outperform!

AW: Flowsensor A(PC) autoclaving problems?

Inbox x Blueprint Medical x

09:44



Radev, N. via viamed.co.uk

to Steve, Bernd

Hi Steve,

the sterilization process at 134°C of SpiroTrueAPC was many times validated at bluepoint using steam sterilizer type WEBECO and periodically serviced by WEBECO.

For validation of the sterilizing process we used short and standard programs of the sterilizer

Short program cycle is 8min at 134°C

Standard program cycle is 26min at 134°

There were not found any problems on SpiroTrueAPC at all. (s. attached pictures)

SpiroTrueAPC was tested 200+ cycles in short program without having any problems. Based on that we decided to be officially we came out with up to 100 cycles sterilizing at 134°C on SpiroTrueAPC are possible.

Please ensure that your customer got SpiroTrueA(PC) and not SpiroTrueA(ABS). If you are supplying both types to the such mix up can be happen.

Please ensure that your customer is trained by you about sensor correct orientation into sterilizer (s. attached pictures) any other heavy stuffs over the sensors during sterilizing. As we know it is typically mistake all medical stuffs for sterilizing together with other sensitive parts in one basket and then put all together into sterilizer. Additionally it is important to ensure that the sterilizer is correctly calibrated and that there are no temperature deviations. Once the sterilizing is complete the parts are still very hot and if the ambient temperature is not correct then mechanical damages are possible. However we didn't detect any failures on SpiroTrueAPC during sterilizing at 134°C

You can check by yourself the sterilizing of SpiroTrueAPC using some of your friendly customers and keep these test results for you.