





Examples of some of the screen formats from the AX400 & MX400 design file.



Trend displays: 4 hrs, 8hrs, 24 hrs.





# VIAMED

Viamed Limited - 15 Station Road - Cross Hills  
Keighley -West Yorkshire BD20 7DT - United Kingdom  
Tel: +44 1535 634542 Fax: +44 1535 635582  
Email: info@viamed.co.uk Website: www.viamed.co.uk

## Menu structure examples

### Main Menu

Alarm settings

Data management

Setup

Exit

### Alarm settings

Auto Alarm

Upper     ↑     **50**

Lower     ↓     **18**

Exit

### Data management

Remaining     1490 min

Stored data

Delete all data

Exit

### Stored data


| Date       | Time  | Min |
|------------|-------|-----|
| 09. 18. 07 | 20:35 | 205 |
| 09. 17. 07 | 13:17 | 5   |
| 09. 16. 07 | 11:15 | 122 |
| 09. 15. 07 | 20:55 | 73  |

Exit



# VIAMED

Viamed Limited - 15 Station Road - Cross Hills  
Keighley -West Yorkshire BD20 7DT - United Kingdom  
Tel: +44 1535 634542 Fax: +44 1535 635582  
Email: info@viamed.co.uk Website: www.viamed.co.uk



### Setup

Date and Time

Brightness  
Language  
Service  
Exit

### Calibration mode

Exit

☒ Cal 21%

☐ Cal 100%