

|                              |                                   |                  |
|------------------------------|-----------------------------------|------------------|
| <b>bluepoint</b> MEDICAL     | <b>Biocompatibility</b>           | Revision<br>1    |
| <b>Internal<br/>Document</b> | <b>VersaStream Sampling Lines</b> | Seite<br>1 von 1 |

**Parts, that comes in contact with the patient during intended use**

- The skin of the patient face comes in contact with the nasal or oral/nasal tip, tubes and tube holder
- Indirect contact with patient breath circuit with airway adapter

**Available Biocompatibility Test:**

| Component                                                   | Kind                  | Report Number   | Standard    | Sample Supplier | Result of the test                                                                  |
|-------------------------------------------------------------|-----------------------|-----------------|-------------|-----------------|-------------------------------------------------------------------------------------|
| Oral/nasal tips, tubes, tube holder, splitter and connector | Biocompatibility Test | SDWH-2010-20272 | ISO 10993-5 | Hangzhou Medtec | The test result shows that of the applied sample does not induce irritation to skin |
| Airway Adapter                                              | Cytotoxicity Test     | 124734-20-A     | ISO 10993-5 | Justec          | The material is not cytotoxic                                                       |
|                                                             |                       |                 |             |                 |                                                                                     |