


## Sampling Line Kit



REORDER PN: 2013069-001



Rx ONLY U.S.



This sampling kit is intended for use with the GE Medical Systems Information Technologies **CapnoFlex LF CO2 System** with ventilator and anesthesia circuits that have an integrated airway adapter with a female Luer-Lok® connector.

### Caution:

Verify that the sampling line kit is intact and undamaged.  
Cleaning and sterilization will compromise system performance.  
Remove sample cell from monitoring system when not in use.  
Do not use on patients that cannot tolerate the removal of 50ml/min from their total minute ventilation.  
Do not position cables or tubing in any manner that may cause entanglement or strangulation.



**GE Medical Systems**  
Information Technologies  
www.gehealthcare.com

Manufactured for:  
GE Medical Systems  
Information Technologies  
8200 W. Tower Ave.  
Milwaukee, WI 53223 USA  
US Tel: (800) 558-7044  
Canada Tel: (800) 668-0732

**EC REP**

GE Medical Systems  
Information Technologies GmbH  
Munzinger Strasse 3-5  
D-79111 FREIBURG, GERMANY  
Tel: (49) 761.45.43.0

06061200

**LOT**



US PATENTS 5,693,944 D490,332 D501,802 OTHER PATENTS PENDING

2025989-001 Rev A