



VIAMED

Viamed Limited - 15 Station Road - Cross Hills
Keighley - West Yorkshire BD20 7DT - United Kingdom
Tel: +44 (0)1535 634542 Fax: +44 (0)1535 635582
Email: info@viamed.co.uk Website: www.viamed.co.uk

How to identify a Bristol Maid cot

The unit has single rib/foot at each end of the bottom of the cot. There are four holes in the bottom. There are holes at each end near the top at the centre. There is a manufactures date code moulded into the base.



How to identify an Arjo-Huntleigh cot

The unit has raised moulded handles at each end. The bottom of the cot has six strengthening ribs.

