

Software Validation of Apgar Timer

Apgar Mode Ready to Start	Buttons to press	Expected Outcome	Actual Outcome
	Start		
	Freeze		
	Stop		
	Start + Freeze		
	Start + Stop		
	Freeze + Stop		
	Start+Freeze+Stop		
Counting		Expected Outcome	Actual Outcome
	Start		
	Freeze		
	Stop		
	Start + Freeze		
	Start + Stop		
	Freeze + Stop		
	Start+Freeze+Stop		
Frozen and counting		Expected Outcome	Actual Outcome
	Start		
	Freeze		
	Stop		
	Start + Freeze		
	Start + Stop		
	Freeze + Stop		
	Start+Freeze+Stop		
Stopped		Expected Outcome	Actual Outcome
	Start		
	Freeze		
	Stop		
	Start + Freeze		
	Start + Stop		
	Freeze + Stop		
	Start+Freeze+Stop		