

<<Title>> <<Full Name>>  
<<Job Title>>  
<<Department>>  
<<Company>>  
<<Address 1>>  
<<Address 2>>  
<<City>> <<County>>  
<<Postcode>>

ACCEPTED  
DATE 6/1/16  
INITIAL SH

January 2016

Dear <<Full Name>>,

### Maxtec FLOCAP Flow Indicator & EtCO<sub>2</sub> Detector

Viamed is delighted to introduce the Maxtec FLOCAP into its Capnography range.

The FLOCAP is a single-use EtCO<sub>2</sub> & flow indication device designed for visualisation of a patient's exhaled CO<sub>2</sub> and expiratory flow. The EtCO<sub>2</sub> indication assists the user in verifying proper ET tube placement. The flow indicator helps the caregiver identify end of exhalation which may help prevent breath-stacking or overdistention and its associated health issues during resuscitation.

Please find enclosed the Maxtec FLOCAP leaflet which provides you with further technical specifications.

Should you require any further information please do not hesitate in contacting a member of our sales team on (01535) 634542.

Yours sincerely

*Catrin*

Catrin Hollings  
Marketing

