



Viamed Limited - 15 Station Road - Cross Hills
Keighley -West Yorkshire BD20 7DT - United Kingdom
Tel: +44 (0)1535 634542 Fax: +44 (0)1535 635582
Email: info@viamed.co.uk Website: www.viamed.co.uk

Maxtec Oxygen Monitors – Model Number Clarification

This document serves to clarify the correct replacement sensors required for Maxtec oxygen monitors, models **OM-25 ME** and **MAXO2 ME**.

OM-25 ME



Model **OM-25 ME** has been discontinued, it misleadingly bears the brand **maxO₂®me** on the instrument.

Oxygen sensor:

Viamed part number: 0110429
Maxtec model number: MAX-250E



MAXO2 ME



Model **MAXO2 ME** superseded the **OM-25 ME** and uses a different sensor.

Oxygen sensor:

Viamed part number: 0110452
Maxtec model number: MAX-550E

