


VIAMED

Suppliers of Electromedical Equipment

2022



Wain Wath Force on the River Swale, Keld, Yorkshire Dales

January

M	31	3	10	17	24
T		4	11	18	25
W		5	12	19	26
T		6	13	20	27
F		7	14	21	28
S	1	8	15	22	29
S	2	9	16	23	30

February

M		7	14	21	28
T	1	8	15	22	
W	2	9	16	23	
T	3	10	17	24	
F	4	11	18	25	
S	5	12	19	26	
S	6	13	20	27	

March

M		7	14	21	28
T	1	8	15	22	29
W	2	9	16	23	30
T	3	10	17	24	31
F	4	11	18	25	
S	5	12	19	26	
S	6	13	20	27	

April

M		4	11	18	25
T		5	12	19	26
W		6	13	20	27
T		7	14	21	28
F	1	8	15	22	29
S	2	9	16	23	30
S	3	10	17	24	

May

M	30	2	9	16	23
T	31	3	10	17	24
W		4	11	18	25
T		5	12	19	26
F		6	13	20	27
S		7	14	21	28
S	1	8	15	22	29

June

M		6	13	20	27
T		7	14	21	28
W	1	8	15	22	29
T	2	9	16	23	30
F	3	10	17	24	
S	4	11	18	25	
S	5	12	19	26	

July

M		4	11	18	25
T		5	12	19	26
W		6	13	20	27
T		7	14	21	28
F	1	8	15	22	29
S	2	9	16	23	30
S	3	10	17	24	31

August

M	1	8	15	22	29
T	2	9	16	23	30
W	3	10	17	24	31
T	4	11	18	25	
F	5	12	19	26	
S	6	13	20	27	
S	7	14	21	28	

September

M		5	12	19	26
T		6	13	20	27
W		7	14	21	28
T	1	8	15	22	29
F	2	9	16	23	30
S	3	10	17	24	
S	4	11	18	25	

October

M	31	3	10	17	24
T		4	11	18	25
W		5	12	19	26
T		6	13	20	27
F		7	14	21	28
S	1	8	15	22	29
S	2	9	16	23	30

November

M		7	14	21	28
T	1	8	15	22	29
W	2	9	16	23	30
T	3	10	17	24	
F	4	11	18	25	
S	5	12	19	26	
S	6	13	20	27	

December

M		5	12	19	26
T		6	13	20	27
W		7	14	21	28
T	1	8	15	22	29
F	2	9	16	23	30
S	3	10	17	24	31
S	4	11	18	25	



VIAMED

15 Station Road · Cross Hills · Keighley · West Yorkshire · BD20 7DT · UK
Tel: +44 (0)1535 634 542 · Fax: +44 (0)1535 635 582
info@viamed.co.uk · www.viamed.co.uk





TOF 3D Neuromuscular Transmission Monitor

Quantitatively monitors the level of neuromuscular transmission during surgery or in the intensive care unit, as a means of assessing the depth or degree of muscle relaxation

- Reliable
- Easy to use
- Twitch response measured by accelerometer
- Displayed values
- Downloadable trend data



TOF 3D Neuromuscular Transmission Monitor
REF: 2510091

By means of 3 dimensional acceleromyography, the TOF 3D quantitatively monitors the level of neuromuscular transmission during surgery or in the intensive care unit.

The new TOF 3D is easy to use and, in addition to the standard TOF and PTC modes, there are alternative objective and subjective modes of stimulation available.

The monitor's functions and modes of stimulation can be quickly customised to meet your specific requirements.

 2510102 Main Patient Cable REF: 5750108	 2510103 Stimulation Cable with pinch clamps REF: 5750107	 2510104 Acceleration (AMG) Transducer REF: 5750105	 2510105 Skin Temperature Sensor REF: 5750106	 2520113 Pole Mount Premium Variant A REF: 5750110	 2520114 Pole Mount Standard Variant B REF: 5750110
 2520110 Hand Adapter REF: 5750100	 2520111 Thumb Adapters (Box of 50) REF: 5750101 Can also be used on the toe	 2520112 Eye Adapters (Box of 50) REF: 5750102	 2530120 Sealing Cover (Pack of 3) REF: 5750116 For Main Patient Cable	 2520121 Sealing Cover (Pack of 3) REF: 5750109 For Data Interface Socket	

Accessories marked with a dotted border are NOT included as standard. Specifications subject to change