



15 Station Road
Cross Hills
Keighley
West Yorkshire
BD20 7DT

+44 (0)1535 634542
www.viamed.co.uk



Oxygen Analysers

OxiQuant S & Handi+

2024

OxiQuant S



Used to determine oxygen concentrations in gas mixtures, delivered by oxygen concentrators, blenders...

- Flexible manual calibration
- Auto-off (after 3 minutes)
- Response time: <13 seconds to 90% of value
- Sensor (with integrated battery) can be replaced

Handi+



This analyser is great for spot-checking oxygen percentage in hospitals, clinics and in homecare environments. Includes integrated sensor and battery.

- One-touch calibration
- Auto-off (after 80 seconds)
- Response time: <15 seconds to 90% of value
- Disposable at end of sensor life
- Overmould for extra protection