

MD300C5-3

Fingertip Pulse Oximeter

Pulse Oximeter is a very important and common device to check patient blood oxygen saturation (SpO₂) level and pulse rate. As a small, compact, simple, reliable and durable physiological monitoring device, fingertip pulse oximeter greatly enhances patient care. It is widely applied in clinic, hospital, social medical organizations, first aid, Dental operation, Oxygen bar, etc. The MD300C52 is designed specially for the children. There is a special mark to tell the children the best place to put their fingers.



TECHNICAL SPECIFICATIONS

Application Range

Adult and Pediatrics

SpO₂

| | |
|-------------------|--|
| Measurement range | 70%~99% |
| Accuracy | 70%~99%: $\pm 3\%$; 0%~69% no definition |
| Resolution | 1% |

Pulse Rate

| | |
|---------------|---|
| Measure range | 30bpm~235bpm |
| Accuracy | 30bpm~99bpm, ± 2 bpm; 100bpm~235bpm, $\pm 2\%$ |
| Resolution | 1bpm |

Display

| | |
|------------|----------------------------------|
| Type | OLED display |
| Parameters | SpO ₂ , PR, Pulse bar |

Power Requirements

Two AAA 1.5V, 600mAh alkaline batteries could be continuously operated as long as 30 hours

FEATURES

- ✧ The smallest fingertip pulse oximeter in the world
- ✧ Easy and simple operation
- ✧ Display SpO₂, PR, Pulse bar, and Plethysmogram
- ✧ Low power consumption; automatically power off
- ✧ Six display modes; adjustable brightness
- ✧ Battery-low indicator
- ✧ Suitable for use in family or hospital

Environment Requirements

| | |
|-----------------------|--|
| Operation Temperature | 5°C~40°C |
| Storage Temperature | -20°C~+55°C |
| Ambient Humidity | ≤80% no condensation in operation ≤93% no condensation in storage |
| Atmosphere pressure | 86kPa~106kPa |

Accessories

- 1pc lanyard
- 2pcs AAA-size batteries
- 1pc user manual
- 1pc silicone case
- 1pc hang bag

NOTE: The specifications are subject to changes without prior notice.

AUTHORIZED AGENT

Beijing Choice Electronic Technology Co., Ltd

Add.: Bailangyuan Building B, Rm.1127-1128, Fuxing Road A36, 100039 Beijing, PRC.
Tel: +86 10 88204188 / 88203520 FAX: +86 10 88204632
Email: market@choicemmed.com Website: www.choicemmed.com

ChoiceMMed

